



## Affordable, Fully-Manual Solution

Economical, fully manual solution carries virtually any power chair or scooter.

### Features & Benefits

- Oversized platform carries most scooters and power chairs
- No scooter or power chair modifications required
- 28" x 48" platform with 60" folding ramp
- Unit folds for storage or transit without mobility device
- Class III hitch adapter standard; Class II hitch adapter available
- Carries up to 500 lbs.

Specifications	AL001
Accommodates	Power chairs and scooters
Platform	28" x 48" platform with 60" folding ramp
Capacity	500 lbs. (227 kg)
Installed Weight	84 lbs.(38.1 kg)
Power	Manual
Securement	Manual Straps (Optional)
Hitch	Class II or III; includes Class III adaptor
Options	Class II hitch adaptor Power chair or scooter covers, strap kit
Warranty	1-year transferable