

Afiscooter-S3

The Afiscooter-S is winner of the prestigious Red Dot Design Award. It offers an impressive driving range, an ergonomic design to make the ride more comfortable, a programmable electronic controller, LCD control screen and LED headlights. The Afiscooter-S has a fully adjustable tiller and an integrated canopy for protection from sun and light rain. The Afiscooter-S has a fully adjustable orthopedic seat, full suspension and a shock-absorbing steel frame. The 3 wheel model has superb maneuverability and unrivaled ease of steering.



reddot design award
winner 2011



Afiscooter-S3 Single seat with canopy



Afiscooter-S3 Extra wide seat with wide wheels and canopy



Afiscooter-S3 Single seat with wide wheels



Afiscooter-S3 Single seat

Colors available for these models

Scan the code
to watch the video »



Afiscooters

Afiscooters develops and manufactures advanced, electrically-powered mobility scooters since 1978. Afiscooters' mobility scooters are designed with reliability and driver safety in mind. Our worldwide distribution network offers year-round local technical support. Afiscooters' scooters meet European ISO 9001-2000 and EN-12184 certification and US safety standards with FDA approval.

1-800-809-3010 (US only)

www.afiscooters.com

Features

- ▶ Fully adjustable swivel orthopedic seat with a movable headrest and armrests
- ▶ Adjustable reach configuration using a single button for forward/reverse
- ▶ Bright, easy to read LCD display control panel indicates speed, battery status, lighting and signaling
- ▶ Powerful, long-lasting, energy efficient LED headlights provide more light and a greater driving range



The sun/rain canopy is winner of the 2011 Red Dot international design award



Improved, secure lockable rear storage
Front storage for items such as glasses and mobile phone



Specifications

Specifications		AfiScooter-S3 - Single		AfiScooter-S3 - Wide seat	
Overall length		61"	155 cm	61"	155 cm
Overall width		27"	68 cm	--	--
	With wide wheels	35"	88 cm	35"	88 cm
Overall height	Without canopy	49"	125 cm	49"	125 cm
	With canopy	65"	164 cm	65"	164 cm
Ground clearance		5"	12.5 cm	5"	12.5 cm
Weight	With batteries	337 lb	153 kg	359 lb	163 kg
	Without batteries	234 lb	106 kg	256 lb	116 kg
	Add for canopy	+ 26 lb	+ 12 kg	+ 26 lb	+ 12 kg
	Add for wide wheels	+ 22 lb	+ 10 kg	+ 22 lb	+ 10 kg
Travel range		28 m	45 km	28 m	45 km
Weight capacity Inc. driver		450 lb	200 kg	450 lb	200 kg
Maximum speed		9.3 mph	15 Km/h	9.3 mph	15 Km/h
Maximun climbing angle		20%		20%	
Turning radius		56"	143 cm	56"	143 cm
Seating		18"	46 cm	33"	85 cm
Tires type		Pneumatic		Pneumatic	
Tires size	Front	16X7.5-8.0"		3X10"	
	Rear	3.5X10"		16X7.5-8.0"	
Battery type		Two 12V 73A/h		Two 12V 73A/h	
Power unit		Rear wheel drive, 24V DC motor		Rear wheel drive, 24V DC motor	
Color		Metallic Silver / Blue / Dark Grey			

Tested by manufacturer, based on actual parameters.



AFIKIM
ELECTRIC VEHICLES

1-800-809-3010
www.afiscooters.com