

# Afiscooter-S4

The Afiscooter-S is winner of the prestigious Red Dot Design Award. It offers an impressive driving range, an ergonomic design to make the ride more comfortable, a programmable electronic controller, LCD control screen and LED headlights. The Afiscooter-S has a fully adjustable tiller and an integrated canopy for protection from sun and light rain. The Afiscooter-S has a fully adjustable orthopedic seat and a shock-absorbing steel frame. The 4 wheel model suspension and steering systems provide a pleasurable car-like driving experience on long distance trips.



reddot design award  
winner 2011

## Afiscooter-S4 Single seat



**Afiscooter-S4**  
Single seat with canopy



**Afiscooter-S4**  
single seat



**Afiscooter-S4**  
Extra wide seat with canopy

Colors available for these models

Scan the code  
to watch the video >>



### Afiscooters

Afiscooters develops and manufactures advanced, electrically-powered mobility scooters since 1978. Afiscooters' mobility scooters are designed with reliability and driver safety in mind. Our worldwide distribution network offers year-round local technical support. Afiscooters' scooters meet European ISO 9001-2000 and EN-12184 certification and US safety standards with FDA approval.

1-800-809-3010 (US only)

[www.afiscooters.com](http://www.afiscooters.com)

## Features

- ▶ Fully adjustable swivel orthopedic seat with a movable headrest and armrests
- ▶ Adjustable reach configuration using a single button for forward/reverse
- ▶ Bright, easy to read LCD display control panel indicates speed, battery status, lighting and signaling
- ▶ Powerful, long-lasting, energy efficient LED headlights provide more light and a greater driving range



The sun/rain canopy is winner of the 2011 Red Dot international design award



Improved, secure lockable rear storage  
Front storage for items such as glasses and mobile phone



## Specifications

		Afiscooter-S4 - Single		Afiscooter-S4 - Wide seat	
Overall length		65"	165 cm	65"	165 cm
Overall width		27"	68 cm	33"	85 cm
	With wide wheels	35"	88 cm	35"	88 cm
Overall height	Without canopy	49"	125 cm	49"	125 cm
	With canopy	65"	164 cm	65"	164 cm
Ground clearance		5"	12.5 cm	5"	12.5 cm
Weight	With batteries	342 lb	155 kg	364 lb	165 kg
	Without batteries	238 lb	108 kg	260 lb	118 kg
	Add for canopy	+ 26 lb	+ 12 kg	+ 26 lb	+ 12 kg
	Add for wide wheels	+ 22 lb	+ 10 kg	+ 22 lb	+ 10 kg
Travel range		28 m	45 km	28 m	45 km
Weight capacity Inc. driver		450 lb	200 kg	450 lb	200 kg
Maximum speed		9.3 mph	15 Km/h	9.3 mph	15 Km/h
Maximun climbing angle		20%		20%	
Turning radius		90"	228 cm	90"	228 cm
Seating		18"	46 cm	33"	85 cm
Tires type		Pneumatic		Pneumatic	
Tires size	Front	4.1/3.5X6.0"		4.1/3.5X6.0"	
	Rear	3.5X10"		3.5X10"	
Battery type		Two 12V 73A/h		Two 12V 73A/h	
Power unit		Rear wheel drive, 24V DC motor		Rear wheel drive, 24V DC motor	
Color		Metallic Silver / Blue / Dark Grey			

Tested by manufacturer, based on actual parameters.



**AFIKIM**  
ELECTRIC VEHICLES

1-800-809-3010  
www.afiscooters.com