

# Afiscooter C3

The outstanding design of the Afiscooter-C is highlighted by these user comfort features:

A fully adjustable swiveling orthopedic seat that gives easy access and an enjoyable ride, as well as the DELTA tiller that provides easy steering and is fully adjustable to accommodate any driver.

A comfortable drive is guaranteed with the advanced full suspension system.

In addition, the Afiscooter-C features a large rear basket for convenient storage.

The three wheel model allows for perfect mobility on narrow sidewalks and easy maneuvering in shopping centres for both indoor and outdoor use.

Available Colours



## **Afikim Electric Vehicles**

Has been developing and manufacturing advanced, electrically-powered mobility scooters since 1978!

Afiscooters' mobility scooters are designed with reliability and driver safety in mind first and foremost alongside beautiful, and thoughtful design.

Afikim Electric Vehicles' worldwide distribution network offers year-round local technical support and meets European ISO 9001-2000, EN 12184 certification, as well as US safety standards with FDA approval.



**AFIKIM**  
ELECTRIC VEHICLES

## Features

- ▶ A fully adjustable orthopedic seat provides comfort while driving and swivels for ease of boarding
- ▶ Fits standard car lifts for the traveling user
- ▶ An optional canopy is available for protection from the sun
- ▶ Easy tiller adjustment for any driver



All operational controls including lights, horn and battery display are conveniently located in the center of the one-hand operation Delta handlebar



Plenty of storage space in the standard: front basket and optional large rear basket



## Specifications

		Afiscooter-C3	
<b>Overall length</b>		51"	130 cm
<b>Overall width</b>		26"	67 cm
<b>Overall height</b>	Without canopy	50"	126 cm
<b>Ground clearance</b>		4"	10 cm
<b>Weight</b>	With batteries	225 lb	102 kg
	Without batteries	161 lb	73 kg
<b>Travel range</b>		25 m	40 km
<b>Weight capacity Inc. driver</b>		330 lb*	135 kg
<b>Maximum speed</b>		9.3 mph	15 Km/h
<b>Maximun climbing angle</b>		16 %	
<b>Turning radius</b>		50"	128 cm
<b>Seating</b>		18"	46 cm
<b>Tires type</b>		Pneumatic	
<b>Tires size</b>	Front	4.1/3.5X6"	
	Rear	4.1/3.5X6"	
<b>Battery type</b>		Two 12V 45-50A/h	
<b>Power unit</b>		Rear wheel drive, 24V DC motor	
<b>Color</b>		Metallic Silver / Blue / Cherry Red	

Tested by manufacturer, based on actual parameters.

\* US Model



**1-800-809-3010** (US only)  
[www.afiscooters.com](http://www.afiscooters.com)

Scan the code to watch the video >>

