


DEALER			COLOURS INFORMATION For Office Use Only		
Dealer:			Order Date: S/N:		
Account:		P.O.#:	RTS		
Ordered By:			RTS: Phone:		
Address:			Email:		
City:		State: Zip:	SHIP TO		
Phone:		Fax:	Name: Ship Via:		
Email:			Address:		
Customer Name:			City: State: Zip:		
Customer Weight:			Phone: Fax:		

CHAIR INFORMATION	PRICE
LITTLE DIPPER K0009	Base Wheel Chair Price K0009 \$ 2,176
Adjustable FSH-RSH	Total Section A \$
Max. Capacity 250Lbs	Total Section B \$
NOTE	Total Section C \$
Main Frame : ZK55 & ZK60 High Grade A7000 Series Aluminum	Total Section D \$
Main Frame : Lifetime Warranty	Total Section E \$
*See warranty policy	Total Section F (Custom Parts Color) \$
	ORDER GRAND TOTAL \$

A FRAME COLORS & PARTS COLORS

FRAME COLORS		
Powder-Co <input type="checkbox"/> Gloss Black <input type="checkbox"/> Twilight Flake <input type="checkbox"/> Yellow <input type="checkbox"/> Magenta <input type="checkbox"/> Fire Engine Red <input type="checkbox"/> White <input type="checkbox"/> Candy Blue <input type="checkbox"/> Candy Purple <input type="checkbox"/> Candy Green <input type="checkbox"/> Candy Red <input type="checkbox"/> Candy Raspberry <input type="checkbox"/> Cobalt Blue <input type="checkbox"/> Matte Gray <input type="checkbox"/> Matte Purple <input type="checkbox"/> Matte Blue <input type="checkbox"/> Matte Green <input type="checkbox"/> Matte Black <input type="checkbox"/> Matte Red <input type="checkbox"/> Neon Red <input type="checkbox"/> Neon Green <input type="checkbox"/> Orange Sparkle <input type="checkbox"/> Mellon Yellow <input type="checkbox"/> Gun Smoke Silver <input type="checkbox"/> Pigeon Blue <input type="checkbox"/> Granny Sparkle Green <input type="checkbox"/> Sparkle Sky Blue <input type="checkbox"/> Vermilion <input type="checkbox"/> Reed Green <input type="checkbox"/> Blue Lilac		\$ STD
Special <input type="checkbox"/> Polished <input type="checkbox"/> 24 karat Gold		\$ 800 \$ CFP
PARTS COLORS (Custom Parts Color Section, Please see Section F)		
STD Color Black Color (*Suspension Springs-Silver Color)		\$ STD
Custom Color Main Bolt on Parts One Color Package	<input type="checkbox"/> Powder-Coated (.....)..	\$ 500
Including: Backrest, Foot Plate, Footrest Rings, Backrest Brackets, Saddle Stops, Forks & Knurl Nuts, Caps, Front Suspension Springs, Caster Housings Axle Plate & Axle Base ☆marked parts	<input type="checkbox"/> Polished <input type="checkbox"/> 24 karat Gold	\$ 750 \$ CFP
Total Section A		\$

B FRAME MEASUREMENTS 1/2



1. SEAT WIDTH	<input type="checkbox"/> 10" <input type="checkbox"/> 11" <input type="checkbox"/> 12" <input type="checkbox"/> 13" <input type="checkbox"/> 14" <input type="checkbox"/> 15"						\$	STD										
2. SEAT DEPTH * Only Choose One Option	Depth Adj. Frame Only	<input type="checkbox"/> Short Frame (11" Depth) Set at <input type="checkbox"/> 8" <input type="checkbox"/> 9" <input type="checkbox"/> 10" <input type="checkbox"/> 11"						\$	STD									
		<input type="checkbox"/> STD Frame (14" Depth) Set at <input type="checkbox"/> 10" <input type="checkbox"/> 11" <input type="checkbox"/> 12" <input type="checkbox"/> 13" <input type="checkbox"/> 14"						\$	STD									
<input type="checkbox"/> Positioning Frame Add (")								\$	0									
3. FOOTREST WIDTH	<input type="checkbox"/> STD (Seat Width Minus 2", No Taper Frame)						\$	STD										
	<input type="checkbox"/> Custom (") *Reverse Taper on Top Only *Measurement must be same or wider than your Seat Width.						\$	150										
		 STD *No Taper																
		 Reverse Taper on Top																
FOOT PLATE STYLE	<input type="checkbox"/> Standard single loop with adjust. height & angle. STD Foot PlateK0040						\$	STD										
	<input type="checkbox"/> Single Flip-Back-single loop (no riser avail) with adjust. height & angle. STD Foot PlateK0040						\$	150										
		<div> <div> </div> <table border="1"> <thead> <tr> <th>Your Footrest Width</th> <th>Foot Plate Size</th> </tr> </thead> <tbody> <tr> <td>8.0" - 10.5"</td> <td>W6.5"xD5"</td> </tr> <tr> <td>11.0" - 12.5"</td> <td>W8.5"xD6.5"</td> </tr> <tr> <td>13.0" - 17.0"</td> <td>W11"xD6.5"</td> </tr> </tbody> </table> </div>							Your Footrest Width	Foot Plate Size	8.0" - 10.5"	W6.5"xD5"	11.0" - 12.5"	W8.5"xD6.5"	13.0" - 17.0"	W11"xD6.5"		
Your Footrest Width	Foot Plate Size																	
8.0" - 10.5"	W6.5"xD5"																	
11.0" - 12.5"	W8.5"xD6.5"																	
13.0" - 17.0"	W11"xD6.5"																	
CUSTOM FOOT PLATE	<input type="checkbox"/> Custom Foot Plate. with adjust. height & angleK0040						\$	75										
	*Available Solid Plate Only <input type="checkbox"/> Drill 2Sets Extra Holes in Footplate																	
	 *Measurement must be 2.5" shorter than your Footrest Width																	
	<input type="checkbox"/> Custom Tubular Footrest						\$	250										
						 *Measurement must be 2.5" shorter than your Footrest Width <input type="checkbox"/> With ABS Foot Plate (\$65 up charge)												
4. FRAME ANGLE	<input type="checkbox"/> 75° <input type="checkbox"/> 80° <input type="checkbox"/> 85° <input type="checkbox"/> 90° to the ground						\$	STD										
	<input type="checkbox"/> Custom (°) to the ground						\$	NCO										

B FRAME MEASUREMENTS 2/2

5. BACKREST HEIGHT * Rigidizer bar 3.25" deep std.		<input type="checkbox"/> 8" <input type="checkbox"/> 9" <input type="checkbox"/> 10" <input type="checkbox"/> 11" <input type="checkbox"/> 12" <input type="checkbox"/> 13" <input type="checkbox"/> 14" <input type="checkbox"/> 15" <input type="checkbox"/> 16" <input type="checkbox"/> 17" <input type="checkbox"/> 18"		\$ STD																																																																																																																																																			
Type <input type="checkbox"/> Locking mechan. with adjust. height & angle \$ STD <input type="checkbox"/> Double Locking mechan. with adjust. height & angle \$ 25 <input type="checkbox"/> Folding with adjust. height & angle \$ NCO <input type="checkbox"/> Custom Backrest (ie: double rigidizer bars, custom depth, custom height, etc) \$ 150																																																																																																																																																							
6. BACKREST ANGLE		Set at (°) to the ground *STD 90°		\$ STD																																																																																																																																																			
7. FRONT SEAT HEIGHT		<input type="checkbox"/> 14.5" <input type="checkbox"/> 15.0" <input type="checkbox"/> 15.5" <input type="checkbox"/> 16.0" <input type="checkbox"/> 16.5" <input type="checkbox"/> 17.0" <input type="checkbox"/> 17.5" <input type="checkbox"/> 18.0" <input type="checkbox"/> 18.5"		\$ STD																																																																																																																																																			
Front Seat Height Ranges * = Not Available with Quick Release Caster Axle. The front seat heights are impacted by the size of your front caster.		<table border="1"> <thead> <tr> <th colspan="2"></th> <th colspan="10">Front Seat Height</th> </tr> <tr> <th colspan="2"></th> <th>14.5"</th> <th>15.0"</th> <th>15.5"</th> <th>16.0"</th> <th>16.5"</th> <th>17.0"</th> <th>17.5"</th> <th>18.0"</th> <th>18.5"</th> </tr> </thead> <tbody> <tr> <td rowspan="8">Your Casters & Forks</td> <td>3" Caster & STD Forks</td> <td><input type="checkbox"/> *</td> <td><input type="checkbox"/> *</td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td>3" Caster & Springs</td> <td></td> <td></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td>3" Caster & Frog Legs</td> <td></td> <td></td> <td></td> <td><input type="checkbox"/> *</td> <td><input type="checkbox"/> *</td> <td><input type="checkbox"/> *</td> <td><input type="checkbox"/> *</td> <td><input type="checkbox"/> *</td> <td><input type="checkbox"/></td> </tr> <tr> <td>4" Caster & STD Forks</td> <td></td> <td></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td>4" Caster & Springs</td> <td></td> <td></td> <td></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td>4" Caster & Frog Legs</td> <td></td> <td></td> <td></td> <td><input type="checkbox"/> *</td> <td><input type="checkbox"/> *</td> <td><input type="checkbox"/> *</td> <td><input type="checkbox"/> *</td> <td><input type="checkbox"/> *</td> <td><input type="checkbox"/> *</td> </tr> <tr> <td>5" Caster & STD Forks</td> <td></td> <td></td> <td></td> <td></td> <td><input type="checkbox"/> *</td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/> *</td> </tr> <tr> <td>5" Caster & Springs</td> <td></td> <td></td> <td></td> <td></td> <td><input type="checkbox"/> *</td> <td><input type="checkbox"/> *</td> <td><input type="checkbox"/> *</td> <td><input type="checkbox"/> *</td> <td><input type="checkbox"/> *</td> </tr> <tr> <td>5" Caster & Frog Legs</td> <td></td> <td></td> <td></td> <td></td> <td></td> <td><input type="checkbox"/> *</td> <td><input type="checkbox"/> *</td> <td><input type="checkbox"/> *</td> <td><input type="checkbox"/> *</td> </tr> </tbody> </table>				Front Seat Height												14.5"	15.0"	15.5"	16.0"	16.5"	17.0"	17.5"	18.0"	18.5"	Your Casters & Forks	3" Caster & STD Forks	<input type="checkbox"/> *	<input type="checkbox"/> *	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	3" Caster & Springs			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	3" Caster & Frog Legs				<input type="checkbox"/> *	<input type="checkbox"/> *	<input type="checkbox"/> *	<input type="checkbox"/> *	<input type="checkbox"/> *	<input type="checkbox"/>	4" Caster & STD Forks			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	4" Caster & Springs				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	4" Caster & Frog Legs				<input type="checkbox"/> *	<input type="checkbox"/> *	<input type="checkbox"/> *	<input type="checkbox"/> *	<input type="checkbox"/> *	<input type="checkbox"/> *	5" Caster & STD Forks					<input type="checkbox"/> *	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/> *	5" Caster & Springs					<input type="checkbox"/> *	<input type="checkbox"/> *	<input type="checkbox"/> *	<input type="checkbox"/> *	<input type="checkbox"/> *	5" Caster & Frog Legs						<input type="checkbox"/> *	<input type="checkbox"/> *	<input type="checkbox"/> *	<input type="checkbox"/> *																																		
		Front Seat Height																																																																																																																																																					
		14.5"	15.0"	15.5"	16.0"	16.5"	17.0"	17.5"	18.0"	18.5"																																																																																																																																													
Your Casters & Forks	3" Caster & STD Forks	<input type="checkbox"/> *	<input type="checkbox"/> *	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																																																																																																																																													
	3" Caster & Springs			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																																																																																																																																													
	3" Caster & Frog Legs				<input type="checkbox"/> *	<input type="checkbox"/> *	<input type="checkbox"/> *	<input type="checkbox"/> *	<input type="checkbox"/> *	<input type="checkbox"/>																																																																																																																																													
	4" Caster & STD Forks			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																																																																																																																																													
	4" Caster & Springs				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																																																																																																																																													
	4" Caster & Frog Legs				<input type="checkbox"/> *	<input type="checkbox"/> *	<input type="checkbox"/> *	<input type="checkbox"/> *	<input type="checkbox"/> *	<input type="checkbox"/> *																																																																																																																																													
	5" Caster & STD Forks					<input type="checkbox"/> *	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/> *																																																																																																																																													
	5" Caster & Springs					<input type="checkbox"/> *	<input type="checkbox"/> *	<input type="checkbox"/> *	<input type="checkbox"/> *	<input type="checkbox"/> *																																																																																																																																													
5" Caster & Frog Legs						<input type="checkbox"/> *	<input type="checkbox"/> *	<input type="checkbox"/> *	<input type="checkbox"/> *																																																																																																																																														
8. REAR SEAT HEIGHT		<input type="checkbox"/> 12.5" <input type="checkbox"/> 13.0" <input type="checkbox"/> 13.5" <input type="checkbox"/> 14.0" <input type="checkbox"/> 14.5" <input type="checkbox"/> 15.0" <input type="checkbox"/> 15.5" <input type="checkbox"/> 16.0" <input type="checkbox"/> 16.5"		\$ STD																																																																																																																																																			
Rear Seat Height Ranges The rear seat heights are impacted by the size of your Tires.		<table border="1"> <thead> <tr> <th>Rear Tire Size</th> <th>Rear Seat Height</th> </tr> </thead> <tbody> <tr> <td>20"x1"(451)</td> <td>12.5"- 15.0"</td> </tr> <tr> <td>20"x1 3/8"(451)</td> <td>13.0"- 15.5"</td> </tr> <tr> <td>22"x1"(501)</td> <td>13.5"- 16.0"</td> </tr> <tr> <td>22"x1 3/8"(501)</td> <td>14.0"- 16.5"</td> </tr> </tbody> </table>		Rear Tire Size	Rear Seat Height	20"x1"(451)	12.5"- 15.0"	20"x1 3/8"(451)	13.0"- 15.5"	22"x1"(501)	13.5"- 16.0"	22"x1 3/8"(501)	14.0"- 16.5"																																																																																																																																										
Rear Tire Size	Rear Seat Height																																																																																																																																																						
20"x1"(451)	12.5"- 15.0"																																																																																																																																																						
20"x1 3/8"(451)	13.0"- 15.5"																																																																																																																																																						
22"x1"(501)	13.5"- 16.0"																																																																																																																																																						
22"x1 3/8"(501)	14.0"- 16.5"																																																																																																																																																						
9. CENTER OF GRAVITY		<input type="checkbox"/> -1.5" <input type="checkbox"/> -1.0" <input type="checkbox"/> -0.5" <input type="checkbox"/> 0" <input type="checkbox"/> 0.5" <input type="checkbox"/> 1.0" <input type="checkbox"/> 1.5" <input type="checkbox"/> 2.0" <input type="checkbox"/> 2.5" <input type="checkbox"/> 3.0" <input type="checkbox"/> 3.5" <input type="checkbox"/> 4.0"		\$ STD																																																																																																																																																			
Center of Gravity Ranges ☆ = Not Available 5" Caster. ★ = Not Available 22" Wheel.		<table border="1"> <thead> <tr> <th colspan="2"></th> <th colspan="12">Your Center of Gravity</th> </tr> <tr> <th colspan="2"></th> <th>-1.5"</th> <th>-1.0"</th> <th>-0.5"</th> <th>0"</th> <th>0.5"</th> <th>1.0"</th> <th>1.5"</th> <th>2.0"</th> <th>2.5"</th> <th>3.0"</th> <th>3.5"</th> <th>4.0"</th> </tr> </thead> <tbody> <tr> <td rowspan="4">Short Frame</td> <td>8"</td> <td></td> <td></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/> ☆</td> <td><input type="checkbox"/> ☆</td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td>9"</td> <td></td> <td></td> <td></td> <td></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/> ☆</td> <td><input type="checkbox"/> ☆</td> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td>10"</td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/> ☆</td> <td><input type="checkbox"/> ☆</td> <td></td> <td></td> </tr> <tr> <td>11"</td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/> ☆</td> <td><input type="checkbox"/> ☆</td> </tr> <tr> <td rowspan="4">STD Frame</td> <td>10"</td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/> ☆</td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td>11"</td> <td></td> <td></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/> ☆</td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td>12"</td> <td></td> <td></td> <td></td> <td></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/> ☆</td> <td></td> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td>13"</td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/> ☆</td> <td></td> <td></td> <td></td> </tr> <tr> <td></td> <td>14"</td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/> ☆</td> </tr> </tbody> </table>				Your Center of Gravity														-1.5"	-1.0"	-0.5"	0"	0.5"	1.0"	1.5"	2.0"	2.5"	3.0"	3.5"	4.0"	Short Frame	8"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/> ☆	<input type="checkbox"/> ☆							9"					<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/> ☆	<input type="checkbox"/> ☆					10"							<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/> ☆	<input type="checkbox"/> ☆			11"									<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/> ☆	<input type="checkbox"/> ☆	STD Frame	10"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/> ☆									11"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/> ☆							12"					<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/> ☆						13"							<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/> ☆					14"								<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/> ☆	
		Your Center of Gravity																																																																																																																																																					
		-1.5"	-1.0"	-0.5"	0"	0.5"	1.0"	1.5"	2.0"	2.5"	3.0"	3.5"	4.0"																																																																																																																																										
Short Frame	8"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/> ☆	<input type="checkbox"/> ☆																																																																																																																																																
	9"					<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/> ☆	<input type="checkbox"/> ☆																																																																																																																																														
	10"							<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/> ☆	<input type="checkbox"/> ☆																																																																																																																																												
	11"									<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/> ☆	<input type="checkbox"/> ☆																																																																																																																																										
STD Frame	10"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/> ☆																																																																																																																																																		
	11"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/> ☆																																																																																																																																																
	12"					<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/> ☆																																																																																																																																															
	13"							<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/> ☆																																																																																																																																													
	14"								<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/> ☆																																																																																																																																											
<input type="checkbox"/> Amputee Axle Plate \$ 107 																																																																																																																																																							
10. CAMBER ANGLE		<input type="checkbox"/> 2° <input type="checkbox"/> 4° <input type="checkbox"/> 6° <input type="checkbox"/> 8°		\$ STD																																																																																																																																																			
11. FRONT OF SEAT SLING TO FOOTREST		<input type="checkbox"/> 12.5" <input type="checkbox"/> 13.0" <input type="checkbox"/> 13.5" <input type="checkbox"/> 14.0" <input type="checkbox"/> 14.5" <input type="checkbox"/> 15.0" <input type="checkbox"/> 15.5" <input type="checkbox"/> 16.0" <input type="checkbox"/> 16.5" <input type="checkbox"/> 17.0"		\$ STD																																																																																																																																																			
*Measurement must be 2.0" shorter than your Front Seat Height HIGH MOUNT <input type="checkbox"/> 6.5" <input type="checkbox"/> 7.0" <input type="checkbox"/> 7.5" <input type="checkbox"/> 8.0" <input type="checkbox"/> 8.5" <input type="checkbox"/> 9.0" <input type="checkbox"/> 9.5" <input type="checkbox"/> 10.0" <input type="checkbox"/> 10.5" <input type="checkbox"/> 11.0" <input type="checkbox"/> 11.5" <input type="checkbox"/> 12.0"		With Risers or Clamps \$ 65																																																																																																																																																					
12. INSIDE WHEEL CLEARANCE		<input type="checkbox"/> 3/4" * <input type="checkbox"/> 1" * <input type="checkbox"/> 1 1/4" <input type="checkbox"/> 1 1/2" <input type="checkbox"/> 1 3/4" <input type="checkbox"/> 2"		\$ STD																																																																																																																																																			
*3/4" and 1" Inside wheel clearance not available with 2 degree camber. *3/4" and 1" Inside wheel clearance not available with Armrest.																																																																																																																																																							
NOTE		Little Dipper is a custom made wheelchair. If you would like to customize your Little Dipper beyond what is available on this order form please submit your request here. Colours will determine as to whether this request can be accomplished while maintaining the safety and integrity of the wheelchair. If it can be accomplished the quote will reflect any and all custom charges, thank you.																																																																																																																																																					
		Sub Total Section B 1/2 \$ Sub Total Section B 2/2 \$ Total Section B \$																																																																																																																																																					





C WHEELS, TIRES & CASTERS, FORKS 1/2

WHEELS														
Standard Radial Spoke Wheel <input type="checkbox"/> 20" (451) <input type="checkbox"/> 22" (501) 		\$ STD												
Mag Wheel <input type="checkbox"/> 20" (451) <input type="checkbox"/> 22" (501) w/Chrome Handrims (Long Tab Only)		\$ STD												
Twister Super Strength wheel <input type="checkbox"/> 20" (451) <input type="checkbox"/> 22" (501) 		\$ 300												
Spinergy Spox 18 Spokes	<input type="checkbox"/> 22" (501) Spoke Color <input type="checkbox"/> Black <input type="checkbox"/> Red <input type="checkbox"/> Yellow <input type="checkbox"/> Blue <input type="checkbox"/> White With Silver Hub With Black Hub With Black Hub With Black Hub With Black Hub	\$ 695												
	<input type="checkbox"/> Custom Color Spoke Pattern <table><tr><th colspan="2">Left Wheel</th><th>Right Wheel</th></tr><tr><td>Spokes</td><td>Inside</td><td></td></tr><tr><td>Color</td><td>Outside</td><td></td></tr><tr><td>Hub Color</td><td colspan="2"><input type="checkbox"/> Black/<input type="checkbox"/> Silver <input type="checkbox"/> Black/<input type="checkbox"/> Silver</td></tr></table>	Left Wheel		Right Wheel	Spokes	Inside		Color	Outside		Hub Color	<input type="checkbox"/> Black/ <input type="checkbox"/> Silver <input type="checkbox"/> Black/ <input type="checkbox"/> Silver		\$ 795
	Left Wheel		Right Wheel											
	Spokes	Inside												
	Color	Outside												
Hub Color	<input type="checkbox"/> Black/ <input type="checkbox"/> Silver <input type="checkbox"/> Black/ <input type="checkbox"/> Silver													
Spinergy LX 12 Spokes	<input type="checkbox"/> 22" (501) Spoke Color <input type="checkbox"/> Black <input type="checkbox"/> Red <input type="checkbox"/> Yellow <input type="checkbox"/> Blue With Silver Hub With Black Hub With Black Hub With Black Hub With Black Hub <input type="checkbox"/> White <input type="checkbox"/> Neon Green <input type="checkbox"/> Pink <input type="checkbox"/> Orange With Black Hub With Black Hub With Black Hub With Black Hub	\$ 895												
	<input type="checkbox"/> Custom Color Spoke Pattern <table><tr><th colspan="2">Left Wheel</th><th>Right Wheel</th></tr><tr><td>Spokes</td><td>Inside</td><td></td></tr><tr><td>Color</td><td>Outside</td><td></td></tr><tr><td>Hub Color</td><td colspan="2"><input type="checkbox"/> Black/<input type="checkbox"/> Silver <input type="checkbox"/> Black/<input type="checkbox"/> Silver</td></tr></table>	Left Wheel		Right Wheel	Spokes	Inside		Color	Outside		Hub Color	<input type="checkbox"/> Black/ <input type="checkbox"/> Silver <input type="checkbox"/> Black/ <input type="checkbox"/> Silver		\$ 995
	Left Wheel		Right Wheel											
	Spokes	Inside												
	Color	Outside												
Hub Color	<input type="checkbox"/> Black/ <input type="checkbox"/> Silver <input type="checkbox"/> Black/ <input type="checkbox"/> Silver													
Pro Wheels <input type="checkbox"/> 20" (451) <input type="checkbox"/> 22" (501) w/ Aluminum Color Anodized Handrim Spoke & Handrim Color <input type="checkbox"/> Anodized Black <input type="checkbox"/> Anodized Red <input type="checkbox"/> Anodized Blue *Not available with D's Lock		\$ 825												
Wijit Drive System *Call for availability		\$ 5,488												
HANDRIMS														
Type	<input type="checkbox"/> Standard For 20", 22" Wheel Anodized Black *Not available with Mag wheel, Pro wheels 	\$ STD												
	<input type="checkbox"/> Vinyl Coated For 20", 22" Wheel *Not available with Mag wheel, Pro wheels 	\$ 142												
	<input type="checkbox"/> Natural Fit For 22" Wheel Only *not available with Mag wheel. <input type="checkbox"/> Standard Grip (Not available with short tabs.) <input type="checkbox"/> Super Grip (Not available with short tabs.) <input type="checkbox"/> Without Thumb Grip 	\$ 295 \$ 340 \$ 270												
	<input type="checkbox"/> Natural Fit (LT) For 22" Wheel Only *not available with Mag wheel. <input type="checkbox"/> Standard Grip (Not available with short tabs.) <input type="checkbox"/> Super Grip (Not available with short tabs.) <input type="checkbox"/> Without Thumb Grip 	\$ 325 \$ 370 \$ 300												
	<input type="checkbox"/> The Surge Handrims *not available with Mag wheel. 	\$ 315												
	<input type="checkbox"/> The Surge Handrims (LT) *not available with Mag wheel. 	\$ 340												
TAB LENGTH <input type="checkbox"/> Long Tabs(STD inside clearance 3/4") 		\$ STD												
<input type="checkbox"/> Short Tabs(Inside clearance 1/2") *Not available with Mag wheel 		\$ STD												
<input type="checkbox"/> Short Tab (Set closest to wheel cut off excess) *Not available with Mag Wheel, Pro Wheel.		\$ 15												
NOTE														

C WHEELS, TIRES & CASTERS, FORKS 2/2

TIRES						
<input type="checkbox"/> Low Pressure Tire <input type="checkbox"/> 20"x1-3/8" or <input type="checkbox"/> 22"x1-3/8" Gray Only						E2211/E2212 \$ STD
<input type="checkbox"/> High Pressure Tire 22"x1" Gray Tire with black wall Only						E2211/E2212 \$ NCO
<input type="checkbox"/> Solid Tire <input type="checkbox"/> 20"x1", <input type="checkbox"/> 22"x1" <input type="checkbox"/> Dark Gray <input type="checkbox"/> Black(\$30up charge)						E2218 \$ 135
<input type="checkbox"/> E2213 Flat-Free Inserts For 20"x1-3/8" or 22"x1-3/8" Tires Only						E2211/E2213 \$ 117
CASTERS						
SIZE	TIRE WIDTH	HUB COLOR	TIRE COLOR	DESCRIPTION		
<input type="checkbox"/> 3"	1"	Black Hub	Black Tire Only	Polyeurathane		\$ STD
<input type="checkbox"/> 3"	1"	Plastic Black Hub	Clear Tire Only	Lighted, Polyeurathane		\$ 65
<input type="checkbox"/> 3"	1"	Red Hub	Red Tire Only	Polyeurathane		\$ STD
<input type="checkbox"/> 3"	1"	Black Hub	Clear Tire Only	Polyeurathane		\$ STD
<input type="checkbox"/> 3"	1"	Aluminum Silver Hub	Black Tire Only	Made By Flog Legs		\$ 165
<input type="checkbox"/> 3"	1"	Aluminum Black Hub	Black Tire Only	Made By Flog Legs		\$ 165
<input type="checkbox"/> 3"	1.5"	Aluminum Silver Hub	Black Tire Only	Soft Roll Made By Flog Legs		\$ 165
<input type="checkbox"/> 3"	1.5"	Aluminum Black Hub	Black Tire Only	Soft Roll Made By Flog Legs		\$ 165
<input type="checkbox"/> 4"	1"	Plastic Black Hub	Black Tire Only			\$ STD
<input type="checkbox"/> 4"	1"	Plastic Black Hub	Clear Only	Lighted, Polyeurathane		\$ 65
<input type="checkbox"/> 4"	1"	Plastic Black Hub	Black Tire Only	Made By Flog Legs		\$ 125
<input type="checkbox"/> 4"	1"	Aluminum Silver Hub	Black Tire Only	Made By Flog Legs		\$ 165
<input type="checkbox"/> 4"	1"	Aluminum Black Hub	Black Tire Only	Made By Flog Legs		\$ 165
<input type="checkbox"/> 4"	1.5"	Composite Black Hub	Tire <input type="checkbox"/> Black <input type="checkbox"/> Blue <input type="checkbox"/> Red	Soft Roll Made By Flog Legs		\$ 125
<input type="checkbox"/> 4"	1.5"	Aluminum Silver Hub	Tire <input type="checkbox"/> Black <input type="checkbox"/> Blue <input type="checkbox"/> Red	Soft Roll Made By Flog Legs		\$ 165
<input type="checkbox"/> 4"	1.5"	Aluminum Black Hub	Black Tire Only	Soft Roll Made By Flog Legs		\$ 165
<input type="checkbox"/> 5"	1"	Composite Black Hub	Black Tire Only			\$ STD
<input type="checkbox"/> 5"	1"	Plastic Black Hub	Clear Tire Only	Lighted, Polyeurathane		\$ 65
<input type="checkbox"/> 5"	1"	Composite Black Hub	Black Tire Only	Made By Flog Legs		\$ 125
<input type="checkbox"/> 5"	1"	Aluminum Silver Hub	Black Tire Only	Made By Flog Legs		\$ 165
<input type="checkbox"/> 5"	1"	Aluminum Black Hub	Black Tire Only	Made By Flog Legs		\$ 165
<input type="checkbox"/> 5"	1.5"	Composite Black Hub	Tire <input type="checkbox"/> Black <input type="checkbox"/> Blue <input type="checkbox"/> Red	Soft Roll Made By Flog Legs		\$ 135
<input type="checkbox"/> 5"	1.5"	Aluminum Silver Hub	Tire <input type="checkbox"/> Black <input type="checkbox"/> Blue <input type="checkbox"/> Red	Soft Roll Made By Flog Legs		\$ 165
<input type="checkbox"/> 5"	1.5"	Aluminum Black Hub	Black Tire Only	Soft Roll Made By Flog Legs		\$ 165
FORKS						
<input type="checkbox"/> STD Aluminum Forks Black Only						\$ STD
Extra Options for Aluminum Forks <input type="checkbox"/> With Suspension						E1015 \$ 55
<input type="checkbox"/> With Quick Release Axles						K0073 \$ 107
<input type="checkbox"/> Frog Leg Suspension Forks Std 1" Wide (For 3", 4" x 1" width Caster and 5" lighted Caster) (Not available with Quick Release Axles)						
Color <input type="checkbox"/> Anodized Black <input type="checkbox"/> Anodized Silver						E1015 \$ 575
<input type="checkbox"/> Frog Leg Suspension Forks 1.5" Wide (For 3", 4" x 1.5" width Caster and 5" Caster) (Not available with Quick Release Axles)						
Color <input type="checkbox"/> Anodized Black <input type="checkbox"/> Anodized Silver						E1015 \$ 575
<div style="text-align: right;">Sub Total Section C 1/2</div>						\$
<div style="text-align: right;">Sub Total Section C 2/2</div>						\$
<div style="text-align: right;">Total Section C</div>						\$

D COMPONENTS

WHEEL LOCKS										
<input type="checkbox"/> Push to Lock / Anodized Black Only		\$ STD								
<input type="checkbox"/> NEW Compact Push to Lock / Anodized Black Only <input type="checkbox"/> w/ 2" Extension Handle (\$65up)		\$ 110								
<input type="checkbox"/> NEW Compact Scissor Lock / Anodized Black Only		\$ 110								
<input type="checkbox"/> D's Locks (*Not available with X-Core, Topolino WX2.5, Pro Wheels) (*Does not include installation fee)		\$ CFP								
PUSH HANDLES *If you are ordering push handles, You must select the backrest locking mechanism.										
<input type="checkbox"/> Integrated w/Fixed Backrest Height / Anodized Black Only		\$ STD								
<input type="checkbox"/> Bolt-on / Anodized Black Only Specify Height from floor -()		\$ 142								
<input type="checkbox"/> Bolt-on w/Quick Release / Anodized Black Only Specify Height from floor -()		\$ 190								
<input type="checkbox"/> NEW Height Adj. / Anodized Black Only Specify Maximum height from floor -()		\$ 295								
NEW ARM RESTS Anodized Black Only,										
<input type="checkbox"/> Fixed Tubular Armrest Adj. Angle and Height Height from seat sling Set at : <input type="checkbox"/> 8" <input type="checkbox"/> 9" <input type="checkbox"/> 10" <input type="checkbox"/> 11" <input type="checkbox"/> 12" E0973		\$ 250								
<input type="checkbox"/> Swing Away Tubular Armrest Adj. Angle and Height Height from seat sling Set at : <input type="checkbox"/> 8" <input type="checkbox"/> 9" <input type="checkbox"/> 10" <input type="checkbox"/> 11" <input type="checkbox"/> 12" E0973		\$ 250								
<input type="checkbox"/> Flip-Up Armrest <input type="checkbox"/> w/Fixed Tubular Armrest Height from seat sling Set at : <input type="checkbox"/> 8" <input type="checkbox"/> 9" <input type="checkbox"/> 10" <input type="checkbox"/> 11" <input type="checkbox"/> 12" E0973		\$ 250								
Adj. Angle <input type="checkbox"/> w/ Swing Away Tubular Armrest Height from seat sling Set at : <input type="checkbox"/> 8" <input type="checkbox"/> 9" <input type="checkbox"/> 10" <input type="checkbox"/> 11" <input type="checkbox"/> 12" E0973		\$ 250								
<input type="checkbox"/> Custom Length & Height (\$150 up charge) *Specify Height & Length → 										
ANTI-TIPPERS										
<input type="checkbox"/> Anti-Tippers (pair) w/ Quick Release / Anodized Black Only	E0971	\$ 112								
SIDE GUARDS										
<input type="checkbox"/> Rigid Side Guards - Black ABS <input type="checkbox"/> Non-Removable <input type="checkbox"/> Removable <input type="checkbox"/> Removable w/Hole (\$10 up charge)		\$ 142								
Size: <input type="checkbox"/> STD <input type="checkbox"/> Medium <input type="checkbox"/> Cut to Wheel (\$98 up charge) Measurement										
Color: Black Only										
<input type="checkbox"/> Rigid Side Guards - Aluminum <input type="checkbox"/> Non-Removable <input type="checkbox"/> Removable Front of Seat  Height <table border="1"><tr><td>7.5"</td><td>STD</td></tr><tr><td>8.25"</td><td>Medium</td></tr></table> Size <table border="1"><tr><td>7.5"</td><td>STD</td></tr><tr><td>8.25"</td><td>Medium</td></tr></table> Seat Depth 	7.5"	STD	8.25"	Medium	7.5"	STD	8.25"	Medium	Size: <input type="checkbox"/> Medium Only <input type="checkbox"/> Cut to Wheel (\$98 up charge) Measurement Color <input type="checkbox"/> Black (STD) <input type="checkbox"/> Pol. Aluminum (STD) <input type="checkbox"/> Powder Coated () (\$50 up charge) <input type="checkbox"/> Gold Plate (+CFP) 	\$ 197
7.5"	STD									
8.25"	Medium									
7.5"	STD									
8.25"	Medium									
<input type="checkbox"/> PRO-TEX Fender Guards / (Black ABS -Only) <input type="checkbox"/> Cut to Fit Specs <input type="checkbox"/> Do Not Install		\$ 175								
<input type="checkbox"/> Fabric Clothing Guards Black Nylon w/Velcro		\$ 78								
SEAT CUSHION -*Please see section E for other cushion choices.-										
<input type="checkbox"/> Foam Cushion w/Nylon Black Cover Thickness <input type="checkbox"/> 1" <input type="checkbox"/> 2" <input type="checkbox"/> 3" <input type="checkbox"/> 4"		\$ 78								
<input type="checkbox"/> Custom Cushion Cover (Std. Foam Cushion Only)-Select Fabric ()										
MISCELLANEOUS										
<input type="checkbox"/> Frame Protectors / Black Leather/Zipper Type		\$ 78								
<input type="checkbox"/> K0065 Spoke Protectors / Black Only <input type="checkbox"/> 24" <input type="checkbox"/> 25" <input type="checkbox"/> 26" K0065		\$ 107								
<input type="checkbox"/> E0978 Seat Belt w/Auto Buckle (BP) 1.5" Black Belt Only E0978		\$ 69								
<input type="checkbox"/> Tie Down Brackets 4pcs.		\$ 210								
SPECIALTY PARTS <input type="checkbox"/> Wheelchair Tray Size: <input type="checkbox"/> W14"x D12" <input type="checkbox"/> W16"x D14" <input type="checkbox"/> W18"x D14" Color: <input type="checkbox"/> Raw Alum. <input type="checkbox"/> Black Alum. <input type="checkbox"/> Polish Alum.(\$250 up charge)		\$ 250								
UPHOLSTERY CHOICE & COLOR *Positioning Back Available. Please see Section E										
<input type="checkbox"/> Standard Adjustable Ht Backrest Black Only		\$ STD								
<input type="checkbox"/> Adjustable Tension Backrest Black Only (*Not available with 8" Backrest Height)		\$ 135								
Custom Options <input type="checkbox"/> Fabrics for Fun <input type="checkbox"/> 001 <input type="checkbox"/> 002 <input type="checkbox"/> 003 <input type="checkbox"/> 004 <input type="checkbox"/> 006 <input type="checkbox"/> 010 <input type="checkbox"/> 011 . . . (\$95 up charge)										
For Std. or Adj. Style <input type="checkbox"/> 014 <input type="checkbox"/> 015 <input type="checkbox"/> 016 <input type="checkbox"/> 017 <input type="checkbox"/> 018 <input type="checkbox"/> 020 <input type="checkbox"/> 021 <input type="checkbox"/> 022 <input type="checkbox"/> 023 <input type="checkbox"/> 024 <input type="checkbox"/> 025 <input type="checkbox"/> 026 <input type="checkbox"/> 027 <input type="checkbox"/> 028										
<input type="checkbox"/> Backrest Leather Trim (outside panels) (\$77 up charge)										
<input type="checkbox"/> Calf Strap <input type="checkbox"/> Black <input type="checkbox"/> Matching Fabric (\$55 up charge) -Select Fabric() K0038		\$ 24								
<input type="checkbox"/> Matching Fabric on Padded Side Guards -Select Fabric ()		\$ 210								
Total Section D \$										

CHAIR OPTIONS 1/2

SEAT CUSHION									
SUPRACOR® WHEELCHAIR CUSHIONS									
Stimulite® Contoured Cushion		<input type="checkbox"/> W14"xD16" (CD1416) E2603						\$ 525	
*Black Upholstery only		<input type="checkbox"/> W16"xD16" (CD1616) E2603						\$ 525	
		<input type="checkbox"/> W16"xD18" (CD1618) E2603						\$ 525	
		<input type="checkbox"/> W18"xD16" (CD1816) E2603						\$ 525	
		<input type="checkbox"/> W18"xD18" (CD1818) E2603						\$ 525	
		<input type="checkbox"/> W18"xD20" (CD1820) E2603						\$ 540	
		<input type="checkbox"/> W20"xD16" (CD2016) E2603						\$ 540	
		<input type="checkbox"/> W20"xD18" (CD2018) E2603						\$ 540	
		<input type="checkbox"/> W20"xD20" (CD2020) E2603						\$ 555	
		<input type="checkbox"/> Extra Sensitive (\$30 up charge)							
		<input type="checkbox"/> Sling Style Bottom (\$45 up charge)							
		<input type="checkbox"/> Custom Size W(") x D(") (\$- CFP)							
Stimulite® Slimline Cushion		<input type="checkbox"/> W14"xD16" (SM1416) E2603						\$ 495	
*Black Upholstery only		<input type="checkbox"/> W16"xD16" (SM1616) E2603						\$ 495	
		<input type="checkbox"/> W16"xD18" (SM1618) E2603						\$ 495	
		<input type="checkbox"/> W18"xD18" (SM1818) E2603						\$ 495	
		<input type="checkbox"/> W18"xD20" (SM1820) E2603						\$ 495	
		<input type="checkbox"/> W20"xD16" (SM2016) E2603						\$ 510	
		<input type="checkbox"/> W20"xD18" (SM2018) E2603						\$ 510	
		<input type="checkbox"/> W20"xD20" (SM2020) E2603						\$ 510	
		<input type="checkbox"/> Extra Sensitive (\$30 up charge)							
		<input type="checkbox"/> Sling Style Bottom (\$45 up charge)							
		<input type="checkbox"/> Custom Size W(") x D(") (\$- CFP)							
Stimulite® Corbee Pediatric		<input type="checkbox"/> W10"xD10" (PD1010) E2603						\$ 450	
*Black Upholstery only		<input type="checkbox"/> W12"xD12" (PD1212) E2603						\$ 450	
		<input type="checkbox"/> W14"xD14" (PD1414) E2603						\$ 450	
		<input type="checkbox"/> Extra Sensitive (\$30 up charge)							
		<input type="checkbox"/> Sling Style Bottom (\$45 up charge)							
		<input type="checkbox"/> Custom Size W(") x D(") (\$- CFP)							
ROHO WHEELCHAIR CUSHIONS									
4" HIGH PROFILE® QUADTRO SELECT® CUSHION		*Black Upholstery E2622-2623						\$ 525	
Choose Size	<input type="checkbox"/> W12"xD12"	<input type="checkbox"/> W12"xD13"	<input type="checkbox"/> W13"xD13"	<input type="checkbox"/> W13"xD15"	<input type="checkbox"/> W13"xD16"	<input type="checkbox"/> W15"xD13"	<input type="checkbox"/> W15"xD15"		
	<input type="checkbox"/> W15"xD16"	<input type="checkbox"/> W15"xD18"	<input type="checkbox"/> W16"xD15"	<input type="checkbox"/> W16"xD16"	<input type="checkbox"/> W16"xD18"	<input type="checkbox"/> W16"xD20"	<input type="checkbox"/> W18"xD15"		
	<input type="checkbox"/> W18"xD16"	<input type="checkbox"/> W18"xD18"	<input type="checkbox"/> W18"xD20"	<input type="checkbox"/> W18"xD22"	<input type="checkbox"/> W20"xD16"	<input type="checkbox"/> W20"xD18"	<input type="checkbox"/> W20"xD20"		
	<input type="checkbox"/> W20"xD22"	<input type="checkbox"/> W22"xD18"	<input type="checkbox"/> W22"xD20"	<input type="checkbox"/> W22"xD22"	<input type="checkbox"/> Custom Size W(") x D(") (\$+85)				
3" MID PROFILE® QUADTRO SELECT® CUSHION		*Black Upholstery only E2622-2623						\$ 525	
Choose Size	<input type="checkbox"/> W15"xD15"	<input type="checkbox"/> W15"xD16"	<input type="checkbox"/> W16"xD16"	<input type="checkbox"/> W16"xD18"	<input type="checkbox"/> W18"xD16"	<input type="checkbox"/> W18"xD18"			
	<input type="checkbox"/> Custom Size W(") x D(") (\$+85)								
2" LOW PROFILE® QUADTRO SELECT® CUSHION		*Black Upholstery only E2622-2623						\$ 525	
Choose Size	<input type="checkbox"/> W12"xD12"	<input type="checkbox"/> W12"xD13"	<input type="checkbox"/> W13"xD13"	<input type="checkbox"/> W13"xD15"	<input type="checkbox"/> W13"xD16"	<input type="checkbox"/> W15"xD13"	<input type="checkbox"/> W15"xD15"		
	<input type="checkbox"/> W15"xD16"	<input type="checkbox"/> W15"xD18"	<input type="checkbox"/> W16"xD15"	<input type="checkbox"/> W16"xD16"	<input type="checkbox"/> W16"xD18"	<input type="checkbox"/> W16"xD20"	<input type="checkbox"/> W18"xD15"		
	<input type="checkbox"/> W18"xD16"	<input type="checkbox"/> W18"xD18"	<input type="checkbox"/> W18"xD20"	<input type="checkbox"/> W20"xD16"	<input type="checkbox"/> W20"xD18"	<input type="checkbox"/> W20"xD20"			
	<input type="checkbox"/> Custom Size W(") x D(") (\$+85)								
2" NEXUS SPIRIT® CUSHION		*Black Upholstery only E2622						\$ 525	
Choose Size	<input type="checkbox"/> W14"xD14"	<input type="checkbox"/> W14"xD16"	<input type="checkbox"/> W15"xD15"	<input type="checkbox"/> W16"xD16"	<input type="checkbox"/> W16"xD18"	<input type="checkbox"/> W17"xD17"	<input type="checkbox"/> W18"xD16"		
	<input type="checkbox"/> W18"xD18"	<input type="checkbox"/> W20"xD18"	* Custom Size Not Available						
STAR CUSHIONS									
Galaxy Cushion		*Black Upholstery only				Choose Zoning Options		E2624	\$ 498
<input type="checkbox"/> W14"xD14" (GXY 1414-2)		<input type="checkbox"/> W14"xD16" (GXY 1416-2)		⇒		<input type="checkbox"/> Standard - Dual Valve			
<input type="checkbox"/> W16"xD14" (GXY 1614-2)		<input type="checkbox"/> W16"xD16" (GXY 1616-2)				<input type="checkbox"/> Dual Valve - Front to Back			
<input type="checkbox"/> W16"xD18" (GXY 1618-2)		<input type="checkbox"/> W18"xD16" (GXY 1816-2)				<input type="checkbox"/> Four Valve			
<input type="checkbox"/> W18"xD18" (GXY 1818-2)		<input type="checkbox"/> W18"xD20" (GXY 1820-2)				<input type="checkbox"/> Simple - Single Valve			
<input type="checkbox"/> W20"xD16" (GXY 2016-2)		<input type="checkbox"/> W20"xD18" (GXY 2018-2)				<input type="checkbox"/> Four Valve - Enhanced			
<input type="checkbox"/> W20"xD20" (GXY 2020-2)						<input type="checkbox"/> Four Valve - Front to Back			
THE CURVE CUSHION		<input type="checkbox"/> W10"xD10"	<input type="checkbox"/> W10"xD11"	<input type="checkbox"/> W10"xD11"	<input type="checkbox"/> W10"xD12"	<input type="checkbox"/> W11"xD12"	<input type="checkbox"/> W12"xD12"	\$ 133	
*Black Upholstery only		<input type="checkbox"/> W12"xD13"	<input type="checkbox"/> W12"xD14"	<input type="checkbox"/> W13"xD11"	<input type="checkbox"/> W13"xD12"	<input type="checkbox"/> W13"xD13"	<input type="checkbox"/> W13"xD14"		
		<input type="checkbox"/> W14"xD14"	<input type="checkbox"/> W14"xD15"	<input type="checkbox"/> W14"xD16"	<input type="checkbox"/> W15"xD14"	<input type="checkbox"/> W15"xD15"	<input type="checkbox"/> W15"xD16"		
		<input type="checkbox"/> W16"xD14"	<input type="checkbox"/> W16"xD15"	<input type="checkbox"/> W16"xD16"	<input type="checkbox"/> W16"xD17"	<input type="checkbox"/> W16"xD18"	<input type="checkbox"/> W17"xD16"		
		<input type="checkbox"/> W17"xD17"	<input type="checkbox"/> W17"xD18"	<input type="checkbox"/> W18"xD16"	<input type="checkbox"/> W18"xD17"	<input type="checkbox"/> W18"xD18"			
NOTE									

E CHAIR OPTIONS 2/2

POSITIONING BACK			
<p>*Note : Actual height of the backrest will increase 2" or more.</p>			
<p>ADI CARBON BACKREST (Support Contour-ACTIVE 2.5") *Not Available with Welded Side Guards</p> <p>Mounting Hardware: <input type="checkbox"/> 2 Point Pro <input type="checkbox"/> Quick Release</p>			
Height	<input type="checkbox"/> LOW (10" Tall) *Backrest Height requires 12" or more *Available 13"-19" width <input type="checkbox"/> Medium (13" Tall) *Backrest Height requires 15" or more *Available 13"-19" width <input type="checkbox"/> Standard (16" Tall) *Backrest Height requires 18" or more *Available 13"-19" width <input type="checkbox"/> Custom Size (" Tall) x (" width)	Width (")	E2613 \$ 745 E2613 \$ 815 E2613 \$ 880 E2613 \$ CFP
<p>ADI ALUMINUM BACKREST (Support Contour-ACTIVE 2.5") *Not Available with Welded Side Guards</p> <p>Mounting Hardware: <input type="checkbox"/> 2 Point Pro <input type="checkbox"/> Quick Release</p>			
Height	<input type="checkbox"/> LOW (10" Tall) *Backrest Height requires 12" or more *Available 10"-20" width <input type="checkbox"/> Medium (13" Tall) *Backrest Height requires 15" or more *Available 12"-20" width <input type="checkbox"/> Standard (16" Tall) *Backrest Height requires 18" or more *Available 14"-20" width <input type="checkbox"/> Tall (18" Tall) *Backrest Height requires 20" or more *Available 16"-20" width <input type="checkbox"/> X-Tall (20" Tall) *Backrest Height requires 22" or more *Available 16"-20" width <input type="checkbox"/> Custom Size (" Tall) x (" width)	Width (")	E2613 \$ 610 E2613 \$ 625 E2613 \$ 640 E2613 \$ 655 E2613 \$ 670 E2613 \$ CFP
<p>ADI ALUMINUM BACKREST (Support Contour-DEEP 4") *Not Available with Welded Side Guards *Not Available with Lateral Supports</p> <p>Mounting Hardware: <input type="checkbox"/> 2 Point Pro <input type="checkbox"/> Quick Release</p>			
Height	<input type="checkbox"/> LOW (10" Tall) *Backrest Height requires 12" or more *Available 10"-20" width <input type="checkbox"/> Medium (13" Tall) *Backrest Height requires 15" or more *Available 12"-20" width <input type="checkbox"/> Standard (16" Tall) *Backrest Height requires 18" or more *Available 14"-20" width <input type="checkbox"/> Tall (18" Tall) *Backrest Height requires 20" or more *Available 16"-20" width <input type="checkbox"/> X-Tall (20" Tall) *Backrest Height requires 22" or more *Available 16"-20" width <input type="checkbox"/> Custom Size (" Tall) x (" width)	Width (")	E2615 \$ 710 E2615 \$ 725 E2615 \$ 740 E2615 \$ 755 E2615 \$ 770 E2615 \$ CFP
<p>ADI LATERAL SUPPORTS (Lateral Supports are compatible with only active 2.5" ADI positioning back.)</p> <p><input type="checkbox"/> Stationary Pad Pad Size <input type="checkbox"/> Small (3" H x 3" D) <input type="checkbox"/> Medium(3" H x 4" D) <input type="checkbox"/> Large(4.5" H x 4.5" D) <input type="checkbox"/> Large(4" H x 6" D)</p> <p><input type="checkbox"/> Pair \$ 320</p> <p><input type="checkbox"/> Left side only or <input type="checkbox"/> Right side only \$ 160</p> <p><input type="checkbox"/> Swing-Out Pad Size <input type="checkbox"/> Small (3" H x 3" D) <input type="checkbox"/> Medium(3" H x 4" D) <input type="checkbox"/> Large(4.5" H x 4.5" D) <input type="checkbox"/> Large(4" H x 6" D)</p> <p><input type="checkbox"/> Pair \$ 520</p> <p><input type="checkbox"/> Left side only or <input type="checkbox"/> Right side only \$ 260</p>			
OTHER			
<input type="checkbox"/> ADI Transfer Board <input type="checkbox"/> 21" (ADV21) <input type="checkbox"/> 24" (ADV24) <input type="checkbox"/> 29"			\$ 62
<input type="checkbox"/> ADI Push and Transfer Gloves Leather (one size fits)			\$ 60
<input type="checkbox"/> ADI Carbon Solid Seat (Wide: 18" Max, Depth: 20") ※Colours Side Guards Not Available with ADI Carbon Solid Seat			\$ 625
With / <input type="checkbox"/> ADI Folding Carbon Side Guard			\$ 400
With / <input type="checkbox"/> ADI Fixed Carbon Side Guard			\$ 375
<p>VARILITE® Seating Systems Available</p> <p style="text-align: right;">Please visit for all options. www.varilite.com</p>			CFP
NOTE			
Sub Total Section E 1/2			\$
Sub Total Section E 2/2			\$
Total Section E			\$

F SPECIAL PARTS COLORS (Please allow an additional 1-2 weeks delivery)
PARTS COLORS

Powder-Coated	Gloss Black	Twilight Flake	Yellow	Magenta	Fire Engine Red	White
	Candy Blue	Candy Purple	Candy Green	Candy Red	Candy Raspberry	Cobalt Blue
	Matte Gray	Matte Purple	Matte Blue	Matte Green	Matte Black	Matte Red
	Neon Red	Neon Green	Orange Sparkle	Mellon Yellow	Gun Smoke Silver	Pigeon Blue
		Granny Sparkle Green	Sparkle Sky Blue	Vermilion	Reed Green	Blue Lilac


Special Colors Polished Aluminum
24 karat Gold Plated

Please Choose the Color for each Part

*Parts not selected will be black color. (Suspension Spring: Silver Color)

*Minimum 3 parts selected per order.

*PC=Powder-Coated



Saddle Stops

<input type="checkbox"/> PC ()	\$40
<input type="checkbox"/> Polished.....	\$115
<input type="checkbox"/> Gold Plated.....	\$CFP

Backrest

<input type="checkbox"/> PC ()	\$125
<input type="checkbox"/> Polished.....	\$225
<input type="checkbox"/> Gold Plated.....	CFP

Backrest Brackets

<input type="checkbox"/> PC ()	\$90
<input type="checkbox"/> Polished.....	\$115
<input type="checkbox"/> Gold Plated.....	\$CFP

Axle Plates

<input type="checkbox"/> PC ()	\$65
<input type="checkbox"/> Polished.....	\$85
<input type="checkbox"/> Gold Plated.....	\$CFP

Caps

<input type="checkbox"/> PC ()	\$20
<input type="checkbox"/> Polished.....	N/A
<input type="checkbox"/> Gold Plated.....	N/A

Caster Housing

<input type="checkbox"/> PC ()	\$150
<input type="checkbox"/> Polished.....	\$280
<input type="checkbox"/> Gold Plated.....	N/A

Footrest Rings

<input type="checkbox"/> PC ()	\$20
<input type="checkbox"/> Polished.....	\$40
<input type="checkbox"/> Gold Plated.....	\$CFP

Foot Plate Clamps

<input type="checkbox"/> PC ()	\$50
<input type="checkbox"/> Polished.....	\$75
<input type="checkbox"/> Gold Plated.....	\$CFP

Footrest Tubes

<input type="checkbox"/> PC ()	\$65
<input type="checkbox"/> Polished.....	\$85
<input type="checkbox"/> Gold Plated.....	\$CFP

Foot Plate

<input type="checkbox"/> PC ()	\$45
<input type="checkbox"/> Polished.....	\$90
<input type="checkbox"/> Gold Plated.....	\$CFP

Front Springs

<input type="checkbox"/> PC ()	\$20
<input type="checkbox"/> Polished.....	N/A
<input type="checkbox"/> Gold Plated.....	N/A


Forks & Knurl Nuts

<input type="checkbox"/> PC ()	\$65
<input type="checkbox"/> Polished.....	\$90
<input type="checkbox"/> Gold Plated.....	\$CFP


Total Section F \$

1. SEAT WIDTH


Measure from outside of seat tube to the outside of the opposite seat tube.

**2. SEAT DEPTH**


Measure from at the front bottom of backrest tube to the front edge of seat sling.

**3. FOOTREST WIDTH**


Measure from inside of front frame tube to inside of opposite front frame tube.

**4. FRAME ANGLE**


Measure angle frame front to floor.

**5. BACKREST HEIGHT**


Measure from top of back post to seat sling at front bottom of backrest tube.

**6. BACKREST ANGLE**


Measure angle from front of back post to floor.

**7. FRONT SEAT HEIGHT**

Measure from floor to top of seat sling at front edge.


**8. REAR SEAT HEIGHT**

Measure from floor to top of seat sling at front bottom of backrest tube.


**9. CENTER OF GRAVITY**

Measure from center of back tube at the folding pivot to the center of the rear axle.


(If welding backrest is selected, Measure from center of back tube at the seat tube to the center of the rear axle.)

**10. CAMBER ANGLE**

Measure angle from wheel to vertical at floor.

**11. FRONT OF SEAT SLING TO FOOTREST****ADULT CHAIR**

Measure from top of seat sling at front edge to top of foot plate at front edge.

**PEDIATRIC CHAIR**

Measure vertical from top of seat sling at front edge to top of foot plate.

**12. INSIDE WHEEL CLEARANCE**

Measure from outside of seat back tube to inside of rear tire.

(If side guards are selected, measure from outside of side guards to inside of rear tire.)

