

SABER


ORDER FORM
Ver. 9.0 Rev-072513Tel: 951-808-9131 FAX: 951-808-9949
Toll Free: 800-892-8998 www.colourswheelchair.com

DEALER			COLOURS INFORMATION For Office Use Only		
Dealer:			Order Date: S/N:		
Account:		P.O.#:	RTS		
Ordered By:			RTS: Phone:		
Address:			Email:		
City:	State:	Zip:	SHIP TO		
Phone:	Fax:		Name: Ship Via:		
Email:			Address:		
Customer Name:			City: State: Zip:		
Customer Weight:			Phone: Fax:		
CHAIR INFORMATION			PRICE		
SABER K0005			Base Wheel Chair Price K0005 \$ 2,145		
Adjustable FSH-RSH			Total Section A \$		
Max. Capacity 250Lbs			Total Section B \$		
NOTE			Total Section C \$		
Main Frame : ZK55 & ZK60 High Grade A7000 Series Aluminum			Total Section D \$		
Main Frame : Lifetime Warranty			Total Section E \$		
*See warranty policy			Total Section F(Custom Parts Color) \$		
			ORDER GRAND TOTAL \$		

A FRAME COLORS & PARTS COLORS


FRAME COLORS						
Powder-Co	☐ Gloss Black	☐ Twilight Flake	☐ Yellow	☐ Magenta	☐ Fire Engine Red	☐ White
	☐ Candy Blue	☐ Candy Purple	☐ Candy Green	☐ Candy Red	☐ Candy Raspberry	☐ Cobalt Blue
	☐ Matte Gray	☐ Matte Purple	☐ Matte Blue	☐ Matte Green	☐ Matte Black	☐ Matte Red
	☐ Neon Red	☐ Neon Green	☐ Orange Sparkle	☐ Mellon Yellow	☐ Gun Smoke Silver	☐ Pigeon Blue
	☐ Copper	☐ Granny Sparkle Green	☐ Sparkle Sky Blue	☐ Vermilion	☐ Reed Green	☐ Blue Lilac
						\$ STD
Special	☐ Polished					\$ 800
	☐ 24 karat Gold					\$ CFP
PARTS COLORS						
STD Color	Black Color (*Suspension Springs-Silver Color)					\$ STD
Custom Color Main Bolt on Parts One Color Package				☐ Powder-Coated ()		\$ 500
Includes Backrest, Foot Plate, Footrest Rings, Backrest Brackets, Saddle Stops				☐ Polished		\$ 750
Forks & Knurl Nuts, Caps, Front Suspension Springs, Axle Tube, Axle Clamps, Axle Plates, Axle Bases, Camber Discs ☆marked parts				☐ 24 karat Gold		\$ CFP
Total Section A						\$

B FRAME MEASUREMENTS 1/2



1. SEAT WIDTH	☐ 14" ☐ 15" ☐ 16" ☐ 17" ☐ 18" ☐ 19"						\$	STD																																																									
2. SEAT DEPTH	Depth Adj. Frame		☐ Short Frame (16" Depth) Set at ☐ 14" ☐ 15" ☐ 16" ☐ Long Frame (19" Depth) Set at ☐ 17" ☐ 18" ☐ 19" 				\$	STD																																																									
3. FOOTREST WIDTH	☐ STD (Seat Width Minus 5", Taper on Bottom Frame) *No Custom Modifications						\$	STD																																																									
Seat Width Footrest Width (Seat Width Minus 5")																																																																	
FOOT PLATE STYLE	☐ Standard single loop with adjust. height & angle. STD Foot Plate K0040 ☐ Single Flip-Back-single loop (no riser avail) with adjust. height & angle. STD Foot Plate . . K0040						\$	STD																																																									
STD Foot Plate Size Chart *No Custom Modifications <table border="1" style="flex: 2;"> <thead> <tr> <th>Your Seat Width</th> <th>Foot Plate Size</th> </tr> </thead> <tbody> <tr> <td>14"-15"</td> <td>W6.5"xD5"</td> </tr> <tr> <td>16"-18"</td> <td>W8.5"xD6.5"</td> </tr> <tr> <td>19"</td> <td>W11"xD6.5"</td> </tr> </tbody> </table>									Your Seat Width	Foot Plate Size	14"-15"	W6.5"xD5"	16"-18"	W8.5"xD6.5"	19"	W11"xD6.5"																																																	
Your Seat Width	Foot Plate Size																																																																
14"-15"	W6.5"xD5"																																																																
16"-18"	W8.5"xD6.5"																																																																
19"	W11"xD6.5"																																																																
4. FRAME BEND ANGLE	☐ 95° only *No Custom Modifications						\$	STD																																																									
Frame Bend Angle Frame Angle to the Ground FRAME ANGLE to the Ground <table border="1"> <thead> <tr> <th colspan="2" rowspan="2"></th> <th colspan="6">Your Seat Depth</th> </tr> <tr> <th>14"</th> <th>15"</th> <th>16"</th> <th>17"</th> <th>18"</th> <th>19"</th> </tr> </thead> <tbody> <tr> <th rowspan="6" style="writing-mode: vertical-rl; transform: rotate(180deg);">Seat Bucket</th> <th>0"</th> <td>95°</td> <td>95°</td> <td>95°</td> <td>95°</td> <td>95°</td> <td>95°</td> </tr> <tr> <th>1"</th> <td>91°</td> <td>91°</td> <td>91°</td> <td>92°</td> <td>92°</td> <td>92°</td> </tr> <tr> <th>2"</th> <td>87°</td> <td>87°</td> <td>88°</td> <td>88°</td> <td>89°</td> <td>89°</td> </tr> <tr> <th>3"</th> <td>83°</td> <td>83°</td> <td>84°</td> <td>85°</td> <td>85°</td> <td>86°</td> </tr> <tr> <th>4"</th> <td>78°</td> <td>80°</td> <td>81°</td> <td>81°</td> <td>82°</td> <td>83°</td> </tr> <tr> <th>5"</th> <td>74°</td> <td>76°</td> <td>77°</td> <td>78°</td> <td>79°</td> <td>80°</td> </tr> </tbody> </table>											Your Seat Depth						14"	15"	16"	17"	18"	19"	Seat Bucket	0"	95°	95°	95°	95°	95°	95°	1"	91°	91°	91°	92°	92°	92°	2"	87°	87°	88°	88°	89°	89°	3"	83°	83°	84°	85°	85°	86°	4"	78°	80°	81°	81°	82°	83°	5"	74°	76°	77°	78°	79°	80°
		Your Seat Depth																																																															
		14"	15"	16"	17"	18"	19"																																																										
Seat Bucket	0"	95°	95°	95°	95°	95°	95°																																																										
	1"	91°	91°	91°	92°	92°	92°																																																										
	2"	87°	87°	88°	88°	89°	89°																																																										
	3"	83°	83°	84°	85°	85°	86°																																																										
	4"	78°	80°	81°	81°	82°	83°																																																										
	5"	74°	76°	77°	78°	79°	80°																																																										

B FRAME MEASUREMENTS 2/2

5. BACKREST HEIGHT * Rigidizer bar 3.25" deep std.	☐ 8" ☐ 9" ☐ 10" ☐ 11" ☐ 12" ☐ 13" ☐ 14" ☐ 15" ☐ 16" ☐ 17" ☐ 18" Adj. Range Adj. Range Adj. Range  Type ☐ Locking mechan. with adjust. height & angle \$ STD ☐ Double Locking mechan. with adjust. height & angle \$ 25 ☐ Folding with adjust. height & angle \$ NCO *No Custom Modifications	\$ STD																																																																																																														
6. BACKREST ANGLE	Set at (°) to the ground *STD 90°	\$ STD																																																																																																														
7. FRONT SEAT HEIGHT STD Front Seat Height Rang No Custom Modifications *=Not Available with Quick	The front seat heights are impacted by the size of your front caster. <table border="1" style="width: 80%; border-collapse: collapse; text-align: center;"> <thead> <tr> <th colspan="2" rowspan="2"></th><th colspan="8">Front Seat Height</th></tr> <tr> <th>16.5"</th><th>17.0"</th><th>17.5"</th><th>18.0"</th><th>18.5"</th><th>19.0"</th><th>19.5"</th><th>20.0"</th><th>20.5"</th></tr> </thead> <tbody> <tr><td rowspan="8" style="writing-mode: vertical-rl; transform: rotate(180deg);">Your Casters & Forks</td><td>3"Caster & STD Forks</td><td>☐*</td><td>☐</td><td>☐</td><td>☐</td><td>☐</td><td>☐</td><td>☐</td><td></td><td></td></tr> <tr><td>3"Caster & Springs</td><td>☐*</td><td>☐</td><td>☐</td><td>☐</td><td>☐</td><td>☐</td><td>☐</td><td></td><td></td></tr> <tr><td>3"Caster & Frog Legs</td><td></td><td></td><td>☐*</td><td>☐*</td><td>☐*</td><td>☐*</td><td>☐*</td><td></td><td></td></tr> <tr><td>4"Caster & STD Forks</td><td>☐*</td><td>☐</td><td>☐</td><td>☐</td><td>☐</td><td>☐</td><td>☐</td><td>☐</td><td></td></tr> <tr><td>4"Caster & Springs</td><td></td><td></td><td>☐</td><td>☐</td><td>☐</td><td>☐</td><td>☐</td><td>☐</td><td></td></tr> <tr><td>4"Caster & Frog Legs</td><td></td><td></td><td>☐*</td><td>☐*</td><td>☐*</td><td>☐*</td><td>☐*</td><td>☐*</td><td></td></tr> <tr><td>5"Caster & STD Forks</td><td></td><td>☐*</td><td>☐*</td><td>☐*</td><td>☐</td><td>☐</td><td>☐</td><td>☐</td><td>☐*</td></tr> <tr><td>5"Caster & Springs</td><td></td><td></td><td></td><td>☐*</td><td>☐*</td><td>☐*</td><td>☐*</td><td>☐*</td><td>☐*</td></tr> <tr><td>5"Caster & Frog Legs</td><td></td><td></td><td></td><td></td><td>☐*</td><td>☐*</td><td>☐*</td><td>☐*</td><td>☐*</td></tr> </tbody> </table>			Front Seat Height								16.5"	17.0"	17.5"	18.0"	18.5"	19.0"	19.5"	20.0"	20.5"	Your Casters & Forks	3"Caster & STD Forks	☐ *	☐	☐	☐	☐	☐	☐			3"Caster & Springs	☐ *	☐	☐	☐	☐	☐	☐			3"Caster & Frog Legs			☐ *	☐ *	☐ *	☐ *	☐ *			4"Caster & STD Forks	☐ *	☐	☐	☐	☐	☐	☐	☐		4"Caster & Springs			☐	☐	☐	☐	☐	☐		4"Caster & Frog Legs			☐ *	☐ *	☐ *	☐ *	☐ *	☐ *		5"Caster & STD Forks		☐ *	☐ *	☐ *	☐	☐	☐	☐	☐ *	5"Caster & Springs				☐ *	☐ *	☐ *	☐ *	☐ *	☐ *	5"Caster & Frog Legs					☐ *	☐ *	☐ *	☐ *	☐ *	\$ STD
				Front Seat Height																																																																																																												
		16.5"	17.0"	17.5"	18.0"	18.5"	19.0"	19.5"	20.0"	20.5"																																																																																																						
Your Casters & Forks	3"Caster & STD Forks	☐ *	☐	☐	☐	☐	☐	☐																																																																																																								
	3"Caster & Springs	☐ *	☐	☐	☐	☐	☐	☐																																																																																																								
	3"Caster & Frog Legs			☐ *	☐ *	☐ *	☐ *	☐ *																																																																																																								
	4"Caster & STD Forks	☐ *	☐	☐	☐	☐	☐	☐	☐																																																																																																							
	4"Caster & Springs			☐	☐	☐	☐	☐	☐																																																																																																							
	4"Caster & Frog Legs			☐ *	☐ *	☐ *	☐ *	☐ *	☐ *																																																																																																							
	5"Caster & STD Forks		☐ *	☐ *	☐ *	☐	☐	☐	☐	☐ *																																																																																																						
	5"Caster & Springs				☐ *	☐ *	☐ *	☐ *	☐ *	☐ *																																																																																																						
5"Caster & Frog Legs					☐ *	☐ *	☐ *	☐ *	☐ *																																																																																																							
8. REAR SEAT HEIGHT	☐ 14.5" ☐ 15.0" ☐ 15.5" ☐ 16.0" ☐ 16.5" ☐ 17.0" ☐ 17.5" ☐ 18.0" ☐ 18.5" ☐ 19.0" . . . \$ STD Rear Seat Height Option ☐ w/Colours Air Ride System Black Color \$ 1,580 *Including ;Fox Racing Shox 300 HP Pump ☐ Powder-Coated () (\$380 up charge) ☐ PolishedAluminum (\$800 up charge) *No Custom Modifications <table border="1" style="width: 100%; border-collapse: collapse; text-align: center;"> <thead> <tr><th colspan="2">STD Rear Seat Height Ranges</th></tr> <tr><th>Rear Tire Size</th><th>Rear Seat Height</th></tr> </thead> <tbody> <tr><td>24"x1"(540)</td><td>14.5"- 18.0"</td></tr> <tr><td>24"x1 3/8"(540)</td><td>15.0"- 18.0"</td></tr> <tr><td>25"x1"(559)</td><td>15.5"- 18.0"</td></tr> <tr><td>24"x2"(540)</td><td>15.5"- 18.5"</td></tr> <tr><td>26"x1"(590)</td><td>15.5"- 19.0"</td></tr> <tr><td>25"x1.95"(559)</td><td>16.0"- 19.0"</td></tr> </tbody> </table> *No Custom Modifications <table border="1" style="width: 100%; border-collapse: collapse; text-align: center;"> <thead> <tr><th colspan="2">Rear Seat Height Ranges w/ Air Ride System</th></tr> <tr><th>Rear Tire Size</th><th>Rear Seat Height</th></tr> </thead> <tbody> <tr><td>24"x1"(540)</td><td>16.0"- 18.0"</td></tr> <tr><td>24"x1 3/8"(540)</td><td>16.5"- 18.0"</td></tr> <tr><td>25"x1"(559)</td><td>16.5"- 18.0"</td></tr> <tr><td>24"x2"(540)</td><td>17.0"- 18.5"</td></tr> <tr><td>26"x1"(590)</td><td>17.0"- 19.0"</td></tr> <tr><td>25"x1.95"(559)</td><td>17.5"- 19.0"</td></tr> </tbody> </table> REAR SUSPENSION RATING ☐ Set up 50psi (For 80 lbs-110 lbs weight) ☐ Set up 60psi (For 110 lbs-140 lbs weight) ☐ Set up 80psi (For 140 lbs-170 lbs weight) ☐ Set up 100psi (For 170 lbs-200 lbs weight) ☐ Set up 120psi (For 200 lbs-230 lbs weight) ☐ Set up 150psi (For 230 lbs-250 lbs weight)	STD Rear Seat Height Ranges		Rear Tire Size	Rear Seat Height	24"x1"(540)	14.5"- 18.0"	24"x1 3/8"(540)	15.0"- 18.0"	25"x1"(559)	15.5"- 18.0"	24"x2"(540)	15.5"- 18.5"	26"x1"(590)	15.5"- 19.0"	25"x1.95"(559)	16.0"- 19.0"	Rear Seat Height Ranges w/ Air Ride System		Rear Tire Size	Rear Seat Height	24"x1"(540)	16.0"- 18.0"	24"x1 3/8"(540)	16.5"- 18.0"	25"x1"(559)	16.5"- 18.0"	24"x2"(540)	17.0"- 18.5"	26"x1"(590)	17.0"- 19.0"	25"x1.95"(559)	17.5"- 19.0"	\$ STD \$ 1,580																																																																														
STD Rear Seat Height Ranges																																																																																																																
Rear Tire Size	Rear Seat Height																																																																																																															
24"x1"(540)	14.5"- 18.0"																																																																																																															
24"x1 3/8"(540)	15.0"- 18.0"																																																																																																															
25"x1"(559)	15.5"- 18.0"																																																																																																															
24"x2"(540)	15.5"- 18.5"																																																																																																															
26"x1"(590)	15.5"- 19.0"																																																																																																															
25"x1.95"(559)	16.0"- 19.0"																																																																																																															
Rear Seat Height Ranges w/ Air Ride System																																																																																																																
Rear Tire Size	Rear Seat Height																																																																																																															
24"x1"(540)	16.0"- 18.0"																																																																																																															
24"x1 3/8"(540)	16.5"- 18.0"																																																																																																															
25"x1"(559)	16.5"- 18.0"																																																																																																															
24"x2"(540)	17.0"- 18.5"																																																																																																															
26"x1"(590)	17.0"- 19.0"																																																																																																															
25"x1.95"(559)	17.5"- 19.0"																																																																																																															
9. CENTER OF GRAVITY	☐ 1.0" ☐ 1.5" ☐ 2.0" ☐ 2.5" ☐ 3.0" ☐ 3.5" ☐ 4.0" ☐ 4.5" ☐ 5.0" ☐ 5.5" . . .	\$ STD																																																																																																														
Center of Gravity Ranges ★=Excludes 24" x 2", 26" x 1", 25" x 1.95" tires are selected. ☆=Not available with Air Ride System	<table border="1" style="width: 100%; border-collapse: collapse; text-align: center;"> <thead> <tr> <th colspan="2" rowspan="2"></th><th colspan="8">Your Center of Gravity</th></tr> <tr> <th>1.0"</th><th>1.5"</th><th>2.0"</th><th>2.5"</th><th>3.0"</th><th>3.5"</th><th>4.0"</th><th>4.5"</th><th>5.0"</th><th>5.5"</th></tr> </thead> <tbody> <tr><td rowspan="6" style="writing-mode: vertical-rl; transform: rotate(180deg);">Your Seat Depth</td><td>14"</td><td>☐</td><td>☐★</td><td>☐★</td><td></td><td></td><td></td><td></td><td></td><td></td><td></td></tr> <tr><td>15"</td><td>☐☆</td><td>☐☆</td><td>☐</td><td>☐★</td><td>☐★</td><td></td><td></td><td></td><td></td><td></td></tr> <tr><td>16"</td><td>☐☆</td><td>☐☆</td><td>☐☆</td><td>☐</td><td>☐</td><td>☐★</td><td>☐★</td><td></td><td></td><td></td></tr> <tr><td>17"</td><td>☐☆</td><td>☐☆</td><td>☐</td><td>☐</td><td>☐</td><td>☐</td><td>☐</td><td>☐★</td><td>☐★</td><td></td></tr> <tr><td>18"</td><td>☐☆</td><td>☐☆</td><td>☐</td><td>☐</td><td>☐</td><td>☐</td><td>☐</td><td>☐</td><td>☐</td><td>☐★</td></tr> <tr><td>19"</td><td>☐☆</td><td>☐☆</td><td>☐☆</td><td>☐</td><td>☐</td><td>☐</td><td>☐</td><td>☐</td><td>☐</td><td>☐</td></tr> </tbody> </table>			Your Center of Gravity								1.0"	1.5"	2.0"	2.5"	3.0"	3.5"	4.0"	4.5"	5.0"	5.5"	Your Seat Depth	14"	☐	☐ ★	☐ ★								15"	☐ ☆	☐ ☆	☐	☐ ★	☐ ★						16"	☐ ☆	☐ ☆	☐ ☆	☐	☐	☐ ★	☐ ★				17"	☐ ☆	☐ ☆	☐	☐	☐	☐	☐	☐ ★	☐ ★		18"	☐ ☆	☐ ☆	☐	☐	☐	☐	☐	☐	☐	☐ ★	19"	☐ ☆	☐ ☆	☐ ☆	☐	☐	☐	☐	☐	☐	☐	\$ STD																							
				Your Center of Gravity																																																																																																												
		1.0"	1.5"	2.0"	2.5"	3.0"	3.5"	4.0"	4.5"	5.0"	5.5"																																																																																																					
Your Seat Depth	14"	☐	☐ ★	☐ ★																																																																																																												
	15"	☐ ☆	☐ ☆	☐	☐ ★	☐ ★																																																																																																										
	16"	☐ ☆	☐ ☆	☐ ☆	☐	☐	☐ ★	☐ ★																																																																																																								
	17"	☐ ☆	☐ ☆	☐	☐	☐	☐	☐	☐ ★	☐ ★																																																																																																						
	18"	☐ ☆	☐ ☆	☐	☐	☐	☐	☐	☐	☐	☐ ★																																																																																																					
	19"	☐ ☆	☐ ☆	☐ ☆	☐	☐	☐	☐	☐	☐	☐																																																																																																					
10. CAMBER ANGLE	☐ 2° ☐ 4° 	\$ STD																																																																																																														
11. FRONT OF SEAT SLING TO FRONT OF FOOTREST	☐ 14.5" ☐ 15.0" ☐ 15.5" ☐ 16.0" ☐ 16.5" ☐ 17.0" ☐ 17.5" ☐ 18.0" ☐ 18.5" ☐ 19.0" . .	\$ STD																																																																																																														
HIGH MOUNT With Risers or Clamps	☐ 9.5" ☐ 10.0" ☐ 10.5" ☐ 11.0" ☐ 11.5" ☐ 12.0" ☐ 12.5" ☐ 13.0" ☐ 13.5" ☐ 14.0" 	\$ 65																																																																																																														
12. INSIDE WHEEL CLEARANCE	☐ 3/4" * ☐ 1" ☐ 1 1/4" ☐ 1 1/2" ☐ 1 3/4" ☐ 2" *3/4" Inside wheel clearance not available with 2 degree camber. (with the exception of Y-360 wheels) *3/4" Inside wheel clearance not available with Armrest *3/4" Inside wheel clearance not available with D's Locks	\$ STD																																																																																																														
NOTE																																																																																																																
Sub Total Section B 1/2 \$																																																																																																																
Sub Total Section B 2/2 \$																																																																																																																
Total Section B \$																																																																																																																




C WHEELS, TIRES & CASTERS, FORKS 1/2

WHEELS																
Standard Radial Spoke Wheel ☐ 24" (540) Only				\$ STD												
Heat-Treated Wheel ☐ 24" (540) ☐ 25" (559) ☐ 26" (590)				\$ 250												
Cross Knot Laced Spokes ☐ 24" (540) ☐ 25" (559) ☐ 26" (590)				\$ 250												
Twister Super Strength Wheel ☐ 24" (540) ☐ 25" (559) ☐ 26" (590)				\$ 300												
SUN Fusion 16 Wheel ☐ 24" (540) ☐ 25" (559)				\$ 410												
16 Black Spokes Hub Color ☐ Black ☐ Silver																
Y-360 Billet Wheels ☐ 24" (540) ☐ 25" (559) *Max. Capacity 220Lbs *Camber 0- 2° max *Not for Sports or Rugged use Style ☐ Saw Blade ☐ Hearts(25"only) ☐ T-Style(25"only)				\$ 999												
Y-360 Billet Wheels MOD Style ☐ 25" (559) Only *Max. Capacity 220Lbs *Camber 0- 2° max *Not for Sports or Rugged use Style ☐ Tribal Sun ☐ Girls ☐ Eagle ☐ Holy Cross ☐ Hono Specify Powder Coat Color for Insert Only () *Choose from Frame Colors				\$ 1,499												
Extra Options for Y-360 ☐ Spinners Clover (\$650 up charge) ☐ Powder-Coated Spokes Only () *Choose from Frame Colors (\$150 up charge)																
X-Core (3 Spokes comp) ☐ 25" (559) Only *Camber 0- 2° max *Not available with D's Lock Wheel Color ☐ Black				\$ 743												
☐ Yellow ☐ Blue ☐ Red ☐ White				\$ 943												
Spinergy Spox 18 Spokes	☐ 24" (540) ☐ 25" (559) ☐ 26" (590) Spoke Color ☐ Black ☐ Red ☐ Yellow ☐ Blue ☐ White With Silver Hub With Black Hub With Black Hub With Black Hub With Black Hub			\$ 695												
	☐ Custom Color Spoke Pattern															
	<table><tr><td>Spokes Color</td><td>Inside</td><td></td><td></td></tr><tr><td></td><td>Outside</td><td></td><td></td></tr><tr><td>Hub Color</td><td></td><td>☐ Black/☐ Silver</td><td>☐ Black/☐ Silver</td></tr></table>			Spokes Color	Inside				Outside			Hub Color		☐ Black/ ☐ Silver	☐ Black/ ☐ Silver	\$ 795
	Spokes Color	Inside														
	Outside															
Hub Color		☐ Black/ ☐ Silver	☐ Black/ ☐ Silver													
Spinergy LX 12 Spokes	☐ 24" (540) ☐ 25" (559) ☐ 26" (590) Spoke Color ☐ Black ☐ Red ☐ Yellow ☐ Blue With Silver Hub With Black Hub With Black Hub With Black Hub ☐ White ☐ Neon Green ☐ Pink ☐ Orange With Black Hub With Black Hub With Black Hub With Black Hub			\$ 895												
	☐ Custom Color Spoke Pattern															
	<table><tr><td>Spokes Color</td><td>Inside</td><td></td><td></td></tr><tr><td></td><td>Outside</td><td></td><td></td></tr><tr><td>Hub Color</td><td></td><td>☐ Black/☐ Silver</td><td>☐ Black/☐ Silver</td></tr></table>			Spokes Color	Inside				Outside			Hub Color		☐ Black/ ☐ Silver	☐ Black/ ☐ Silver	\$ 995
	Spokes Color	Inside														
	Outside															
Hub Color		☐ Black/ ☐ Silver	☐ Black/ ☐ Silver													
Topolino Cabon Core WX2.5 Wheels Black Spokes Only ☐ 24" (540) ☐ 25" (559) *Max. Capacity 220Lbs				\$ 1,095												
*Not available with D's Lock Nipples Color ☐ Red ☐ Yellow																
Pro Wheels ☐ 24" (540) ☐ 25" (559)				\$ 825												
w/ Aluminum Color Anodized Handrim Spoke & Handrim Color ☐ Anodized Black ☐ Anodized Red ☐ Anodized Blue *Not available with D's Lock																
HANDRIMS																
Type	☐ Standard For 24", 25" or 26" Wheel Anodized Black (*Not available with Pro Wheels)			\$ STD												
	☐ Vinyl Coated For 24", 25", 26" Wheel (*Not available with Pro Wheels)			\$ 142												
	☐ Natural Fit For 24", 25" Wheel															
	☐ Standard Grip (Not available with short tabs. Not available with Y-360 Billet Wheels, X-Core, Topolino Wheels, Pro Wheels, Mag Wheels) .			\$ 295												
	☐ Super Grip (Not available with short tabs. Not available with Y-360 Billet Wheels, X-Core, Topolino Wheels, Pro Wheels, Mag Wheels) . . .			\$ 340												
	☐ Without Thumb Grip (Not available with Y-360 Billet Wheels, Not available with Pro Wheels, Mag Wheels)			\$ 270												
	☐ Natural Fit (LT) For 24", 25" Wheel															
	☐ Standard Grip (Not available with short tabs. Not available with Y-360 Billet Wheels, X-Core, Topolino Wheels, Pro Wheels, Mag Wheels) .			\$ 325												
	☐ Super Grip (Not available with short tabs. Not available with Y-360 Billet Wheels, X-Core, Topolino Wheels, Pro Wheels, Mag Wheels) . . .			\$ 370												
	☐ Without Thumb Grip (Not available with Y-360 Billet Wheels, Not available with Pro Wheels, Mag Wheels)			\$ 300												
	☐ The Surge Handrims * Not available with Y-360 Billet Wheels * w / X-Core Wheels required Installation fee (\$30) flat fee			\$ 315												
	☐ The Surge Handrims (LT) * Not available with Y-360 Billet Wheels * w / X-Core Wheels required Installation fee (\$30) flat fee			\$ 340												
☐ Spinergy Flex Rim Available only on 24" and 25" Spinergy LX.			\$ 495													
☐ Projections Handrims *Not available with Y-360 Billet Wheels, X-Core, Topolino Wheels, Pro Wheels, Mag Wheels) . . E0967 *Call For Availability			\$ CFP													
TAB LENGTH ☐ Long Tab (STD inside clearance 3/4") * Not available with X-Core Wheels				\$ STD												
☐ Short Tab (Inside clearance 1/2") *Not available with Ultra Grip and Natural Fit Handrim with thumb grip				\$ STD												
☐ Short Tab (Set closest to wheel cut off excess) * Not available with Y-360 wheel, Pro Wheels, X-Core, Pro Wheels				\$ 15												

C WHEELS, TIRES & CASTERS, FORKS 2/2

TIRES					
☐ High Pressure Tire	☐ 24"x1", ☐ 25"x1" or ☐ 26"x1"	Gray Tire with black wall Only		E2211/E2212	\$ STD
☐ High Pressure Tire	☐ 24"x1", ☐ 25"x1" or ☐ 26"x1"	☐ Blue ☐ Red		E2211/E2212	\$ 65
☐ Low Pressure Tire	☐ 24"x1-3/8" or ☐ 26"x1-3/8"	Gray Only		E2211/E2212	\$ STD
☐ E2211Schwalbe Marathon Plus High Pressure Tire	☐ 24"x1", ☐ 25"x1" or ☐ 26"x1"	Black Only		E2211/E2212	\$ 99
☐ Solid Tire	☐ 24"x1", ☐ 25"x1" or ☐ 26"x1"	☐ Dark Gray ☐ Black(\$30up charge)		E2218	\$ 135
☐ E2211 Moderate V-Track Knobby Tire	24"x1-3/8" (37-540)	Gray Only		E2211/E2212	\$ 117
☐ E2211 Aggressive Knobby Tire	24"x 2" (50-540)	Gray Only		E2211/E2212	\$ 117
☐ E2211 Aggressive Knobby Tire	26"x 1.95" (50-559)	Black Only		E2211/E2212	\$ 117
☐ E2213 Flat-Free Inserts For 24"x1-3/8" Tires Only				E2211/E2213	\$ 117
CASTERS					
SIZE	TIRE WIDTH	HUB COLOR	TIRE COLOR	DESCRIPTION	
☐ 3"	1"	Black Hub	Black Tire Only	Polyeurethane	\$ STD
☐ 3"	1"	Plastic Black Hub	Clear Tire Only	Lighted, Polyurethane	\$ 65
☐ 3"	1"	Red Hub	Red Tire Only	Polyurethane	\$ STD
☐ 3"	1"	Black Hub	Clear Tire Only	Polyurethane	\$ STD
☐ 3"	1"	Aluminum Silver Hub	Black Tire Only	Made By Flog Legs	\$ 165
☐ 3"	1"	Aluminum Black Hub	Black Tire Only	Made By Flog Legs	\$ 165
☐ 3"	1.5"	Aluminum Silver Hub	Black Tire Only	Soft Roll Made By Flog Legs	\$ 165
☐ 3"	1.5"	Aluminum Black Hub	Black Tire Only	Soft Roll Made By Flog Legs	\$ 165
☐ 4"	1"	Plastic Black Hub	Black Tire Only		\$ STD
☐ 4"	1"	Plastic Black Hub	Clear Only	Lighted, Polyurethane	\$ 65
☐ 4"	1"	Plastic Black Hub	Black Tire Only	Made By Flog Legs	\$ 125
☐ 4"	1"	Aluminum Silver Hub	Black Tire Only	Made By Flog Legs	\$ 165
☐ 4"	1"	Aluminum Black Hub	Black Tire Only	Made By Flog Legs	\$ 165
☐ 4"	1.5"	Composite Black Hub	Tire ☐ Black ☐ Blue ☐ Red	Soft Roll Made By Flog Legs	\$ 135
☐ 4"	1.5"	Aluminum Silver Hub	Tire ☐ Black ☐ Blue ☐ Red	Soft Roll Made By Flog Legs	\$ 165
☐ 4"	1.5"	Aluminum Black Hub	Black Tire Only	Soft Roll Made By Flog Legs	\$ 165
☐ 5"	1"	Composite Black Hub	Black Tire Only		\$ STD
☐ 5"	1"	Plastic Black Hub	Clear Tire Only	Lighted, Polyurethane	\$ 65
☐ 5"	1"	Composite Black Hub	Black Tire Only	Made By Flog Legs	\$ 125
☐ 5"	1"	Aluminum Silver Hub	Black Tire Only	Made By Flog Legs	\$ 165
☐ 5"	1"	Aluminum Black Hub	Black Tire Only	Made By Flog Legs	\$ 165
☐ 5"	1.5"	Composite Black Hub	Tire ☐ Black ☐ Blue ☐ Red	Soft Roll Made By Flog Legs	\$ 135
☐ 5"	1.5"	Aluminum Silver Hub	Tire ☐ Black ☐ Blue ☐ Red	Soft Roll Made By Flog Legs	\$ 165
☐ 5"	1.5"	Aluminum Black Hub	Black Tire Only	Soft Roll Made By Flog Legs	\$ 165
FORKS					
☐ STD Aluminum Forks	Black Only				\$ STD
Extra Options for Aluminum Forks	☐ With Suspension Springs		E1015	\$ 55	
	☐ With Quick Release Axles		K0073	\$ 107	
☐ Frog Leg Suspension Forks Std 1" Wide (For 3", 4" x 1" width Caster and 5" STD Caster, 5" lighted Caster)	Color ☐ Anodized Black ☐ Anodized Silver		E1015	\$ 575	(Not available with Quick Release Axles)
☐ Frog Leg Suspension Forks 1.5" Wide (For 3", 4" x 1.5" width Caster and 5" x 1" & 5" x 1.5" Casters)	Color ☐ Anodized Black ☐ Anodized Silver		E1015	\$ 575	(Not available with Quick Release Axles)
Sub Total Section C 1/2					\$
Sub Total Section C 2/2					\$
Total Section C					\$


D COMPONENTS

WHEEL LOCKS		
☐ Push to Lock / Anodized Black Only		\$ STD
☐ NEW Compact Push to Lock / Anodized Black Only ☐ w/ 2" Extension Handle (\$65up)		\$ 110
☐ NEW Compact Scissor Lock / Anodized Black Only		\$ 110
☐ D's Locks (*Not available with Air Ride System and X-Core, Topolino WX2.5, Pro Wheels, Magic wheels) (*Does not include installation fee)		\$ CFP
PUSH HANDLES *If you are ordering push handles, You must select the backrest locking mechanism.		
☐ Integrated w/Fixed Backrest Height / Anodized Black Only		\$ STD
☐ Bolt-on / Anodized Black Only Specify Height from floor -()		\$ 142
☐ Bolt-on w/Quick Release / Anodized Black Only Specify Height from floor -()		\$ 190
☐ NEW Height Adj. / Anodized Black Only Specify Maximum height from floor -()		\$ 295
NEW ARM RESTS Anodized Black Only,		
☐ Fixed Tubular Armrest	Adj. Angle and Height Height from seat sling Set at : ☐ 8" ☐ 9" ☐ 10" ☐ 11" ☐ 12" E0973	\$ 250
☐ Swing Away Tubular Armrest	Adj. Angle and Height Height from seat sling Set at : ☐ 8" ☐ 9" ☐ 10" ☐ 11" ☐ 12" E0973	\$ 250
☐ Flip-Up Armrest	☐ w/Fixed Tubular Armrest Height from seat sling Set at : ☐ 8" ☐ 9" ☐ 10" ☐ 11" ☐ 12" E0973	\$ 250
Adj. Angle ☐ w/ Swing Away Tubular Armrest	Height from seat sling Set at : ☐ 8" ☐ 9" ☐ 10" ☐ 11" ☐ 12" E0973	\$ 250
☐ Custom Length & Height (\$150 up charge) *Specify Height & Length →		
		
ANTI-TIPPERS		
☐ Anti-Tippers (pair) w/ Quick Release / Anodized Black Only		E0971 \$ 112
SIDE GUARDS		
☐ Rigid Side Guards - Black ABS	☐ Non-Removable ☐ Removable ☐ Removable w/Hole (\$10 up charge) Size: ☐ STD ☐ Medium ☐ Cut to Wheel (\$98 up charge) Measurement	\$ 142
☐ Rigid Side Guards - Aluminum	☐ Non-Removable ☐ Removable Size: ☐ Medium Only ☐ Cut to Wheel (\$98 up charge) Measurement Color ☐ Black (STD) ☐ Pol. Aluminum (STD) ☐ Powder Coated () (\$50 up charge) ☐ Gold Plate (+CFP)	\$ 197
		
☐ PRO-TEX Fender Guards / (Black ABS -Only) ☐ Cut to Fit Specs ☐ Do Not Install		\$ 175
☐ Fabric Clothing Guards Black Nylon w/Velcro		\$ 78
SEAT CUSHION *Please see section E for other cushion choices.-		
☐ Foam Cushion w/Nylon Black Cover Thickness ☐ 1" ☐ 2" ☐ 3" ☐ 4"		\$ 78
☐ Custom Cushion Cover (Std. Foam Cushion Only) -Select Fabric () (\$90 up charge)		
MISCELLANEOUS		
☐ Frame Protectors / Black Leather/Zipper Type		\$ 78
☐ K0065 Spoke Protectors / Black Only ☐ 24" ☐ 25" ☐ 26"		K0065 \$ 107
☐ E0978 Seat Belt w/Auto Buckle (BP) 1.5" Black Belt Only		E0978 \$ 69
☐ Tie Down Brackets 4pcs.		\$ 210
SPECIALTY PARTS ☐ Wheelchair Tray Size: ☐ W14"x D12" ☐ W16"x D14" ☐ W18"x D14" Color: ☐ Raw Alum. ☐ Black Alum. ☐ Polish Alum. (\$250 up charge)		
UPHOLSTERY CHOICE & COLOR *Positioning Back Available. Please see Section E		
☐ Standard Adjustable Ht Backrest Black Only		\$ STD
☐ Adjustable Tension Backrest Black Only (*Not available with 8" Backrest Height)		\$ 135
Custom Options	☐ Fabrics for Fun ☐ 001 ☐ 002 ☐ 003 ☐ 004 ☐ 006 ☐ 010 ☐ 011 . . . (\$95 up charge) For Std. or Adj. Style ☐ 014 ☐ 015 ☐ 016 ☐ 017 ☐ 018 ☐ 020 ☐ 021 ☐ 022 ☐ 023 ☐ 024 ☐ 025 ☐ 026 ☐ 027 ☐ 028 ☐ Backrest Leather Trim (outside panels) (\$77 up charge)	
☐ Calf Strap ☐ Black ☐ Matching Fabric (\$55 up charge) -Select Fabric()		K0038 \$ 24
☐ Matching Fabric on Padded Side Guards -Select Fabric ()		\$ 210
Total Section D \$		

CHAIR OPTIONS 1/2

SEAT CUSHION													
SUPRACOR® WHEELCHAIR CUSHIONS													
Stimulite® Contoured Cushion *Black Upholstery only		☐ W14"xD16" (CD1416)						E2603	\$ 525				
		☐ W16"xD16" (CD1616)						E2603	\$ 525				
		☐ W16"xD18" (CD1618)						E2603	\$ 525				
		☐ W18"xD16" (CD1816)						E2603	\$ 525				
		☐ W18"xD18" (CD1818)						E2603	\$ 525				
		☐ W18"xD20" (CD1820)						E2603	\$ 540				
		☐ W20"xD16" (CD2016)						E2603	\$ 540				
		☐ W20"xD18" (CD2018)						E2603	\$ 540				
		☐ W20"xD20" (CD2020)						E2603	\$ 555				
				☐ Extra Sensitive (\$30 up charge) ☐ Sling Style Bottom (\$45 up charge) ☐ Custom Size W(") x D(") (\$- CFP)									
Stimulite® Slimline Cushion *Black Upholstery only		☐ W14"xD16" (SM1416)						E2603	\$ 495				
		☐ W16"xD16" (SM1616)						E2603	\$ 495				
		☐ W16"xD18" (SM1618)						E2603	\$ 495				
		☐ W18"xD18" (SM1818)						E2603	\$ 495				
		☐ W18"xD20" (SM1820)						E2603	\$ 495				
		☐ W20"xD16" (SM2016)						E2603	\$ 510				
		☐ W20"xD18" (SM2018)						E2603	\$ 510				
		☐ W20"xD20" (SM2020)						E2603	\$ 510				
				☐ Extra Sensitive (\$30 up charge) ☐ Sling Style Bottom (\$45 up charge) ☐ Custom Size W(") x D(") (\$- CFP)									
Stimulite® Corbee Pediatric *Black Upholstery only		☐ W10"xD10" (PD1010)						E2603	\$ 450				
		☐ W12"xD12" (PD1212)						E2603	\$ 450				
		☐ W14"xD14" (PD1414)						E2603	\$ 450				
		☐ Extra Sensitive (\$30 up charge) ☐ Sling Style Bottom (\$45 up charge) ☐ Custom Size W(") x D(") (\$- CFP)											
ROHO WHEELCHAIR CUSHIONS													
Choose Size	☐ 4" HIGH PROFILE® QUADTRO SELECT® CUSHION		*Black Upholstery					E2622-2623	\$ 525				
	☐ W12"xD12"	☐ W12"xD13"	☐ W13"xD13"	☐ W13"xD15"	☐ W13"xD16"	☐ W15"xD13"	☐ W15"xD15"						
	☐ W15"xD16"	☐ W15"xD18"	☐ W16"xD15"	☐ W16"xD16"	☐ W16"xD18"	☐ W16"xD20"	☐ W18"xD15"						
	☐ W18"xD16"	☐ W18"xD18"	☐ W18"xD20"	☐ W18"xD22"	☐ W20"xD16"	☐ W20"xD18"	☐ W20"xD20"						
	☐ W20"xD22"	☐ W22"xD18"	☐ W22"xD20"	☐ W22"xD22"	☐ Custom Size W(") x D(") (\$+85)								
Choose Size	☐ 3" MID PROFILE® QUADTRO SELECT® CUSHION		*Black Upholstery only					E2622-2623	\$ 525				
	☐ W15"xD15"	☐ W15"xD16"	☐ W16"xD16"	☐ W16"xD18"	☐ W18"xD16"	☐ W18"xD18"	☐ Custom Size W(") x D(") (\$+85)						
Choose Size	☐ 2" LOW PROFILE® QUADTRO SELECT® CUSHION		*Black Upholstery only					E2622-2623	\$ 525				
	☐ W12"xD12"	☐ W12"xD13"	☐ W13"xD13"	☐ W13"xD15"	☐ W13"xD16"	☐ W15"xD13"	☐ W15"xD15"						
	☐ W15"xD16"	☐ W15"xD18"	☐ W16"xD15"	☐ W16"xD16"	☐ W16"xD18"	☐ W16"xD20"	☐ W18"xD15"						
	☐ W18"xD16"	☐ W18"xD18"	☐ W18"xD20"	☐ W20"xD16"	☐ W20"xD18"	☐ W20"xD20"	☐ Custom Size W(") x D(") (\$+85)						
Choose Size	☐ 2" NEXUS SPIRIT® CUSHION		*Black Upholstery only					E2622	\$ 525				
	☐ W14"xD14"	☐ W14"xD16"	☐ W15"xD15"	☐ W16"xD16"	☐ W16"xD18"	☐ W17"xD17"	☐ W18"xD16"						
☐ W18"xD18"										☐ W20"xD18"	* Custom Size Not Available		
STAR CUSHIONS													
Galaxy Cushion *Black Upholstery only		Choose Zoning Options					E2624	\$ 498					
☐ W14"xD14" (GXY 1414-2)		☐ W14"xD16" (GXY 1416-2)		⇒		☐ Standard - Dual Valve							
☐ W16"xD14" (GXY 1614-2)		☐ W16"xD16" (GXY 1616-2)				☐ Dual Valve - Front to Back							
☐ W16"xD18" (GXY 1618-2)		☐ W18"xD16" (GXY 1816-2)				☐ Four Valve							
☐ W18"xD18" (GXY 1818-2)		☐ W18"xD20" (GXY 1820-2)				☐ Simple - Single Valve							
☐ W20"xD16" (GXY 2016-2)		☐ W20"xD18" (GXY 2018-2)				☐ Four Valve - Enhanced							
☐ W20"xD20" (GXY 2020-2)						☐ Four Valve - Front to Back							
THE CURVE CUSHION		☐ W10"xD10"	☐ W10"xD11"	☐ W10"xD11"	☐ W10"xD12"	☐ W11"xD12"	☐ W12"xD12"	\$ 133					
*Black Upholstery only		☐ W12"xD13"	☐ W12"xD14"	☐ W13"xD11"	☐ W13"xD12"	☐ W13"xD13"	☐ W13"xD14"						
		☐ W14"xD14"	☐ W14"xD15"	☐ W14"xD16"	☐ W15"xD14"	☐ W15"xD15"	☐ W15"xD16"						
		☐ W16"xD14"	☐ W16"xD15"	☐ W16"xD16"	☐ W16"xD17"	☐ W16"xD18"	☐ W17"xD16"						
		☐ W17"xD17"	☐ W17"xD18"	☐ W18"xD16"	☐ W18"xD17"	☐ W18"xD18"							
NOTE													

E CHAIR OPTIONS 2/2

POSITIONING BACK		
*Note : Actual height of the backrest will increase 2" or more. 		
ADI CARBON BACKREST (Support Contour-ACTIVE 2.5") *Not Available with Welded Side Guards		
Mounting Hardware: ☐ 2 Point Pro ☐ Quick Release		
Height	☐ LOW (10" Tall) *Backrest Height requires 12" or more *Available 13"-19" width Width (") E2613 \$ 745 ☐ Medium (13" Tall) *Backrest Height requires 15" or more *Available 13"-19" width Width (") E2613 \$ 815 ☐ Standard (16" Tall) *Backrest Height requires 18" or more *Available 13"-19" width Width (") E2613 \$ 880 ☐ Custom Size (" Tall) x (" width) E2613 \$ CFP	
ADI ALUMINUM BACKREST (Support Contour-ACTIVE 2.5") *Not Available with Welded Side Guards		
Mounting Hardware: ☐ 2 Point Pro ☐ Quick Release		
Height	☐ LOW (10" Tall) *Backrest Height requires 12" or more *Available 10"-20" width Width (") E2613 \$ 610 ☐ Medium (13" Tall) *Backrest Height requires 15" or more *Available 12"-20" width Width (") E2613 \$ 625 ☐ Standard (16" Tall) *Backrest Height requires 18" or more *Available 14"-20" width Width (") E2613 \$ 640 ☐ Tall (18" Tall) *Backrest Height requires 20" or more *Available 16"-20" width Width (") E2613 \$ 655 ☐ X-Tall (20" Tall) *Backrest Height requires 22" or more *Available 16"-20" width Width (") E2613 \$ 670 ☐ Custom Size (" Tall) x (" width) E2613 \$ CFP	
ADI ALUMINUM BACKREST (Support Contour-DEEP 4") *Not Available with Welded Side Guards *Not Available with Lateral Supports		
Mounting Hardware: ☐ 2 Point Pro ☐ Quick Release		
Height	☐ LOW (10" Tall) *Backrest Height requires 12" or more *Available 10"-20" width Width (") E2615 \$ 710 ☐ Medium (13" Tall) *Backrest Height requires 15" or more *Available 12"-20" width Width (") E2615 \$ 725 ☐ Standard (16" Tall) *Backrest Height requires 18" or more *Available 14"-20" width Width (") E2615 \$ 740 ☐ Tall (18" Tall) *Backrest Height requires 20" or more *Available 16"-20" width Width (") E2615 \$ 755 ☐ X-Tall (20" Tall) *Backrest Height requires 22" or more *Available 16"-20" width Width (") E2615 \$ 770 ☐ Custom Size (" Tall) x (" width) E2615 \$ CFP	
ADI LATERAL SUPPORTS (Lateral Supports are compatible with only active 2.5" ADI positioning back.)		
☐ Stationary Pad Pad Size ☐ Small (3" H x 3" D) ☐ Medium(3" H x 4" D) ☐ Large(4.5" H x 4.5" D) ☐ Large(4" H x 6" D)		
	☐ Pair \$ 320 ☐ Left side only or ☐ Right side only \$ 160	
☐ Swing-Out Pad Size ☐ Small (3" H x 3" D) ☐ Medium(3" H x 4" D) ☐ Large(4.5" H x 4.5" D) ☐ Large(4" H x 6" D)		
	☐ Pair \$ 520 ☐ Left side only or ☐ Right side only \$ 260	
OTHER		
☐ ADI Transfer Board ☐ 21" (ADV21) ☐ 24" (ADV24) ☐ 29"		\$ 62
☐ ADI Push and Transfer Gloves Leather (one size fits)		\$ 60
☐ ADI Carbon Solid Seat (Wide: 18" Max, Depth: 20") ※Colours Side Guards Not Available with ADI Carbon Solid Seat		\$ 625
With / ☐ ADI Folding Carbon Side Guard		\$ 400
With / ☐ ADI Fixed Carbon Side Guard		\$ 375
VARILITE® Seating Systems Available		
Please visit for all options. www.varilite.com		CFP
NOTE		
Sub Total Section E 1/2 \$		
Sub Total Section E 2/2 \$		
Total Section E \$		

F SPECIAL PARTS COLORS (Please allow an additional 1-2 weeks delivery)

PARTS COLORS						
Powder-Coated	Gloss Black	Twilight Flake	Yellow	Magenta	Fire Engine Red	White
	Candy Blue	Candy Purple	Candy Green	Candy Red	Candy Raspberry	Cobalt Blue
	Matte Gray	Matte Purple	Matte Blue	Matte Green	Matte Black	Matte Red
	Neon Red	Neon Green	Orange Sparkle	Mellon Yellow	Gun Smoke Silver	Pigeon Blue
	Granny Sparkle Green	Sparkle Sky Blue	Vermilion	Reed Green	Blue Lilac	
Special Colors	Polished Aluminum 24 karat Gold Plated *Please allow an additional 4 weeks delivery					
Please Choose the Color for each Part *Parts not selected will be black color. (Suspension Spring: Silver Color) *Minimum 3 parts selected per order.						

*PC=Powder-Coated

Saddle Stops

☐ PC () \$40

☐ Polished . . . \$115

☐ Gold Plated . . . \$CFP

Backrest

☐ PC () \$125

☐ Polished . . . \$225

☐ Gold Plated . . . \$CFP

Backrest Brackets

☐ PC () \$90

☐ Polished . . . \$115

☐ Gold Plated . . . \$CFP

Foot Plate

☐ PC () \$45

☐ Polished . . . \$90

☐ Gold Plated . . . \$CFP

Footrest Rings

☐ PC () \$20

☐ Polished . . . \$40

☐ Gold Plated . . . \$CFP

Footrest Tubes

☐ PC () \$65

☐ Polished . . . \$85

☐ Gold Plated . . . \$CFP

Foot Plate Clamps

☐ PC () \$50

☐ Polished . . . \$75

☐ Gold Plated . . . \$CFP

Caps

☐ PC () \$20

☐ Polished . . . N/A

☐ Gold Plated . . . N/A

Caster Housings

☐ PC () \$125

☐ Polished . . . \$160

☐ Gold Plated . . . \$CFP

Front Springs

☐ PC () \$20

☐ Polished . . . N/A

☐ Gold Plated . . . N/A

Forks & Knurl Nuts

☐ PC () \$65

☐ Polished . . . \$90

☐ Gold Plated . . . \$CFP

Axle Plates

☐ PC () \$65

☐ Polished . . . \$85

☐ Gold Plated . . . \$CFP

Camber Plates

☐ PC () \$50

☐ Polished . . . \$75

☐ Gold Plated . . . \$CFP

Axle Tube

☐ PC () \$50

☐ Polished . . . \$75

☐ Gold Plated . . . \$CFP

Axle Bases

☐ PC () \$50

☐ Polished . . . \$75

☐ Gold Plated . . . \$CFP

Axle Clamps

☐ PC () \$50

☐ Polished . . . \$75

☐ Gold Plated . . . \$CFP

Camber Plate

☐ PC () \$50

☐ Polished . . . \$75

☐ Gold Plated . . . \$CFP

Axle Plate

☐ PC () \$60

☐ Polished . . . \$180

☐ Gold Plated . . . \$CFP

Suspension Main Frame

☐ PC () \$150

☐ Polished . . . \$400

☐ Gold Plated . . . \$CFP

Upper A-Arms

☐ PC () \$100

☐ Polished . . . \$300

☐ Gold Plated . . . \$CFP

Lower A-Arms

☐ PC () \$100


☐ Polished . . . \$280

☐ Gold Plated . . . \$CFP


Total Section F \$

1. SEAT WIDTH


Measure from outside of seat tube to the outside of the opposite seat tube.

**2. SEAT DEPTH**


Measure from at the front bottom of backrest tube to the front edge of seat sling.

**3. FOOTREST WIDTH**


Measure from inside of front frame tube to inside of opposite front frame tube.

**4. FRAME ANGLE**


Measure angle frame front to floor.

**5. BACKREST HEIGHT**


Measure from top of back post to seat sling at front bottom of backrest tube.

**6. BACKREST ANGLE**


Measure angle from front of back post to floor.

**7. FRONT SEAT HEIGHT**

Measure from floor to top of seat sling at front edge.


**8. REAR SEAT HEIGHT**

Measure from floor to top of seat sling at front bottom of backrest tube.


**9. CENTER OF GRAVITY**

Measure from center of back tube at the folding pivot to the center of the rear axle.


(If welding backrest is selected, Measure from center of back tube at the seat tube to the center of the rear axle.)

**10. CAMBER ANGLE**

Measure angle from wheel to vertical at floor.

**11. FRONT OF SEAT SLING TO FOOTREST****ADULT CHAIR**

Measure from top of seat sling at front edge to top of foot plate at front edge.

**PEDIATRIC CHAIR**

Measure vertical from top of seat sling at front edge to top of foot plate.

**12. INSIDE WHEEL CLEARANCE**

Measure from outside of seat back tube to inside of rear tire.

(If side guards are selected, measure from outside of side guards to inside of rear tire.)

