


Dealer Information		Ship To Information	
Account #		Ship Via	
Dealer Name		Name	
Order Date		Address	
Ordered By		City	
Phone		Phone	
Customer Name		Ship Date	Order# :



FRAME MEASUREMENTS

VERVE FOLDING K0009, E1229 \$3,200
(MAX CAPACITY 220Lbs)

*No custom measurement available on this chair.

A SEAT WIDTH

13.4" (34cm) 14.2" (36cm) 14.5" (37cm) STD
 15" (38cm) 15.4" (39cm) 15.75" (40cm)
 16.2" (41cm) 16.5" (42cm) 17" (43cm) 17.3" (44cm)

B SEAT DEPTH

() Choose From 13.5" (34cm)-17" (43cm) STD

C FRONT SEAT HEIGHT

() Choose From 15.75" (40cm)-18" (46cm) STD

D REAR SEAT HEIGHT

() Choose From 14.25" (36cm)-17" (43cm) STD

E SEAT BACK HEIGHT

() Choose From 12" (30cm) -15.75" (40cm) STD

F CENTER OF GRAVITY

1" (2.5cm) 1.5" (3.5cm) 1.8" (4.5cm) STD
 2.2" (5.5cm) 2.6" (6.5cm) 3" (7.7cm)

G BACK ANGLE () °To Floor STD

H FOOTRESTS HEIGHT-Center Joint Folding
Angle & Height Adj. Plastic Cover- Choose Setting
() Choose From 13.5" (34cm) -17" (43cm) STD
*Distance from Seat to Footrest

I REAR WHEEL SPACING (1/4 inch increments)

1" (2.5cm) 1.25" (3cm) 1.5" (3.5cm) 1.75" (4.5cm) STD
 0.5"* (1.3cm) 0.75"* (2cm) *Camber 3°Wheels Only STD

J REAR WHEEL CAMBER

0° 3° STD

K FOOTREST WIDTH

±0" - 2.75" (-7cm) - 3.75" (-9.5cm) STD

L FRAME ANGLE To Floor

90° (STD) 95° 100° STD

M BASE LENGTH

() ≥ **F** + 17" STD

N FRAME LENGTH

() ≥ **B** + 3.5" STD

†*Please fill it in when you have not decided the angle of L.

SIDE GUARDS SHAPE (Welded Side Guard)

SHAPE-I **O**=() **P**=() STD
 Q=() or **S**=()
 SHAPE-II **Q**=() STD
 R=()
 No Side Guards STD

CHAIR OPTIONS		REAR WHEELS - with Quick Release Axles	
COVER OF SIDE GUARDS		<input type="checkbox"/> 22" (520) Standard Wheels (Tangent Cross Spokes)	STD
<input type="checkbox"/> Aluminum Cover (The same color as frame.)	STD	<input type="checkbox"/> 24" (540) Standard Wheels (Tangent Cross Spokes)	STD
<input type="checkbox"/> Nylon Cover Black	STD	24"(540) SPINERGY SPOX LX -Choose Color	\$895
ARMRESTS		<input type="checkbox"/> Black <input type="checkbox"/> White	
<input type="checkbox"/> Omit Armrests	STD	24" (540) Carbon Spinning Wheels -Choose Color	\$900
Tubular Armrests Fixed Frame- Choose Setting		<input type="checkbox"/> Black <input type="checkbox"/> White	
<input type="checkbox"/> 11.5"(28cm) <input type="checkbox"/> 12"(30cm) <input type="checkbox"/> 12.5"(32cm)	\$150	<input type="checkbox"/> Omit Rear Wheels	-\$50
(Distance from Seat to Armrest)		NOTES	
Tubular Armrests Flip Back - Choose Setting		HANDRIMS	
<input type="checkbox"/> 11.5"(28cm) <input type="checkbox"/> 12"(30cm) <input type="checkbox"/> 12.5"(32cm)	\$150	Handrim Distance from Wheel - Choose Setting	
(Distance from Seat to Armrest)		<input type="checkbox"/> 3/8" (1.0cm) <input type="checkbox"/> 5/8" (1.5cm) <input type="checkbox"/> 6/8" (2.0cm) <input type="checkbox"/> 1" (2.5cm) *	STD
BACK UPHOLSTERY		*1" : Not Available with Carbon Spinning Wheels	
<input type="checkbox"/> Tension Adjustable w/Fabric Cushioning ,Black	STD	Handrim Material	
CUSHION		<input type="checkbox"/> Aluminum Hard Anodized Coated	STD
High Density Latex w/Nylon Cover <input type="checkbox"/> 2" (5cm) <input type="checkbox"/> 3" (8cm)	\$195	<input type="checkbox"/> Vinyl Coated *	\$200
WHEEL LOCKS		*Vinyl Coated : Not Available with Carbon Spinning Wheels	
<input type="checkbox"/> Scissor Lock	STD	REAR TIRES & INNER TUBES	
<input type="checkbox"/> Scissor Lock w/Return Spring	STD	22"(25-520) Panaracer (S) Tires w/Woods Valve	STD
<input type="checkbox"/> Toggle - Push/Pull to Lock	\$50	<input type="checkbox"/> Gray <input type="checkbox"/> Blue	
ACCESSORIES		24"(20-540) Panaracer (S) Tires w/Schrader Valve	STD
<input type="checkbox"/> Push Handles	\$75	<input type="checkbox"/> Gray <input type="checkbox"/> Blue	
<input type="checkbox"/> Anti-Tips(Folding)	\$350	<input type="checkbox"/> Omit Rear Tires & Tubes	-\$10
<input type="checkbox"/> Frame Protectors	\$30	NOTES	
<input type="checkbox"/> Calf Strap	\$30	COLOR OF BOLT ON PARTS (Excluding the frame)	
<input type="checkbox"/> Universal Net	\$50	<input type="checkbox"/> Aluminum Silver Anodized Coated	STD
<input type="checkbox"/> Seat Belt	\$50	<input type="checkbox"/> Aluminum Black Anodized Coated	\$50
NOTES		(*1 Including Cross Frame, Wheel Locks, Caster Folks, Footrest, Armrest)	
FRONT WHEELS (3" or 4" caster only. Not Available with 5" and 6")		(*2 The color of each part cannot be specified.)	
<input type="checkbox"/> 3" x 0.8"width Soft Roll Caster Wheels	STD	FRAME COLOR Choose From Paint -Standard Colors	
<input type="checkbox"/> 4" x 0.8"width Soft Roll Caster Wheels	STD	<input type="checkbox"/> Paint Color Name = ()	STD
<input type="checkbox"/> 4" x 1.25"width Soft Roll Caster Wheels	STD	<input type="checkbox"/> Custom Color Paint	Call For Quote
<input type="checkbox"/> 4" x 0.8"width Soft Roll Caster Wheels with Alum Slv. Hub	\$125		
<input type="checkbox"/> 4" x 1.25"width Soft Roll Caster Wheels with Alum Slv. Hub	\$125	Page 1 Grand Total	
<input type="checkbox"/> Omit Front Wheels	-\$20	Page 2 Grand Total	
NOTES		Order Grand Total	
		NOTES	