



Model: X8 Extreme 4x4



201 Growth Parkway Angola, IN 46703

U.S Suggested Retail Price List
Apr-15

Phone: 800-327-0681

Fax: 260-665-3047

Email: iim@mobility-usa.com

www.mobility-usa.com

Quote:

PO #:

Bill To:

Dealer: _____ Contact: _____

Address: _____

City: _____ State: _____ Zip: _____

Phone: _____ Fax: _____

Email: _____ Mark For: _____

Ship To:

Same as Bill To:

Dealer: _____ Contact: _____

Address: _____

City: _____ State: _____ Zip: _____

Phone: _____ Fax: _____

Email: _____ Mark For: _____

Client Information and Specifications:

Demo Performed: Yes No

Client Name: _____ A) Top of Shoulders: _____

B) Chest Depth: _____

Diagnosis: _____ C) Chest Width: _____

D) Seat Depth: _____


Funding Source: _____ E) Top of Head: _____

F) Elbow to Hand: _____

Client Age: _____ G) Knee to Foot: _____

Weight _____ H) Hip Width: _____

Height: _____ I) Seat Pan to Elbow: _____



X8 Extreme Features and Specifications

Joystick/Power Module: Penny and Giles R-net

Length: 45.25"

Width: 28"

Turning Radius: 53"

Weight limit: 400 lbs - 340lbs with lift

Motors: (4) 4 Pole, 24 Volt Motors, 700W each, Gear in Line

Batteries: (2) Group 24, 73 Amp hr

Drive Tires: Black Pneumatic Low Pressure 14" (145/70-6)

Top Speed: 6.2 mph

*performance will vary depending on weight and power options

Warranty: 3 Year Limited on Frames

1 Year Limited on Electronics



Other Products by:

V6 Frontier AT

V6 Frontier Hybrid

V6 Front Wheel Drive

V4 Frontier Rear Wheel Drive

Product Weight- No power seat frame, MPS seat, standard legrests and batteries: 305 lbs

Product Weight- including power tilt, seat elevator, power recline back, power elevating legrests and batteries: 410 lbs

**weights vary depending upon options

X8 Extreme Base Options

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> X8	X8 Base with shroud set up for rehab seat or MPS seat w/ solid seat pan, charger,Pneumatic tires, and standard positioning belt <u>Must Choose:</u> <input type="checkbox"/> Rehab Seat or <input type="checkbox"/> MPS Seat w/ Solid Seat Pan	\$15,200.00	
<input type="checkbox"/> X8C	X8 Base set up for MPS Seat, charger, pneumatic tires, and standard postioning belt	\$15,200.00	

COWL COLOR Choose one

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 3200BK	Midnight Black	N/C	
<input type="checkbox"/> 3200R	Fire Engine Red	N/C	
<input type="checkbox"/> 3200BL	Ocean Blue	N/C	
<input type="checkbox"/> 3200G	Clover Green	N/C	
<input type="checkbox"/> 3200T	Tangerine	N/C	
<input type="checkbox"/> X8CAMO	Realtree AP HD	\$450.00	A9207
<input type="checkbox"/> X8CC	One Custom Color	\$500.00	A9208

BATTERIES

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> BAT	Group 24 Sealed LA/AGM Batteries - installed, two required	\$285.00/ea.	E2363
<input type="checkbox"/> BAT-G	Group 24 Sealed Gel Batteries - installed, two required	\$345.00/ea.	E2363
<input type="checkbox"/> BCSTD	Battery Charger	N/C	E2367

REVERSE DIRECTIONAL STEERING LOCK

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 1325PG	Directional Lock - used to reverse in a straight line	\$1,375.00	

NOTES:

SEATING SYSTEM

Part Number	Description	Price	HCPSC Code
Seating Systems below include: Either MPS Seat Back or Rehab Back Canes with a MPS seat cushion (MPS seat only available with MPS back) or a solid rehab seat pan, height adjustable armrest. (N/C if ordered with a power option)			
<input type="checkbox"/> 2799	MPS Seating System or Rehab Seating and No Power Seating Options	\$2,900.00	
<input type="checkbox"/> 2851	50° Power Tilt in Space *	\$5,050.00	E1002
<input type="checkbox"/> ASKIT	135° Integrated Shear Reduction (ISR) Power Recline Rehab Seating System*	\$4,975.00	E1008
<input type="checkbox"/> 1337	135° Power Recline MPS Seating System*	\$3,500.00	E1003
<input type="checkbox"/> LS	Lift - 12" Seat Elevator - 340 lb. Capacity*	\$4,750.00	E2300
<input type="checkbox"/> 2850	Lift & Tilt Combination - 12" Seat Elevator* - 340 lb. Capacity 50° Power TIS*	\$9,500.00	

Reduced Drive Capability (Option allows for driving capability at 20% speed while elevated and/or tilted)

<input type="checkbox"/> RDC	50° Power TIS RDC	\$700.00	K0108
<input type="checkbox"/> RDC	12" Seat Elevator RDC	\$700.00	K0108

*Must be ordered with upgraded electronics SFCK or MFCK depending on number of options selected

POWER SEAT OPTIONS

Part Number	Description	Price	HCPSC Code
<input type="checkbox"/> SFCK	<u>Single Seat Function Control Kit</u> For operation of a single seat function through the wheelchair electronics via the drive controls and/or membrane switch.	\$1,400.00	E2310
<input type="checkbox"/>	Function Through Joystick	N/C	
<input type="checkbox"/> MSC5-P-S	External 5 Button Box - placed between the armpad and the joystick	\$600.00	
<input type="checkbox"/> MFCK	<u>Multiple Seat Function Control Kit</u> MUST be selected if ordering two or more seat functions. This device is used to operate two or more seat functions	\$2,850.00	E2311
<input type="checkbox"/>	Function Through Joystick	N/C	
<input type="checkbox"/> MSC5-P-S	External 5 Button Box - placed between the armpad and the joystick	\$600.00	

SEAT SIZES (width x depth) Select one of the following:

Part Number	Description	Price	HCPSC Code
MPS SEAT SIZE:			
<input type="checkbox"/> 16 x 16	MPS	N/C	
<input type="checkbox"/> 16 x 18	MPS	N/C	
<input type="checkbox"/> 18 x 16	MPS	N/C	
<input type="checkbox"/> 18 x 18	MPS	N/C	
<input type="checkbox"/> 18 x 20	MPS	N/C	
<input type="checkbox"/> 20 x 18	MPS	N/C	
<input type="checkbox"/> 20 x 20	MPS	N/C	
<input type="checkbox"/> 22 x 18	MPS	\$350.00	K0108
<input type="checkbox"/> 22 x 20	MPS	\$350.00	K0108
<input type="checkbox"/> 23 x 18	MPS	\$350.00	K0108
<input type="checkbox"/> 23 x 20	MPS	\$350.00	K0108

REHAB OR FLAT PAN SEAT SIZE:

<input type="checkbox"/> 2792/2791	Seat Width x Depth	N/C	
	W: _____ L: _____		
	(14"-23") (14"-22")		

SEATBACK OPTIONS

Part Number	Description	Price	HCPSC Code
MPS Options			
<input type="checkbox"/> 1360	Standard MPS Seat Back - 16"-20" Wide - 22" Tall	N/C	
<input type="checkbox"/> 1520	Extra Tall MPS Seat Back - 16"-20" Wide - 26" Tall	\$415.00	K0108
<input type="checkbox"/> 1519	Vinyl Standard MPS Seat Back - 16"-20" Wide - 22" Tall	\$225.00	
<input type="checkbox"/> 1520V	Vinyl Extra Tall MPS Seat Back - 16"-20" Wide - 26" Tall	\$600.00	K0108
<input type="checkbox"/> 1361	Standard MPS Seat Back Wide - 22"-23" Wide - 22" Tall	\$250.00	K0108
<input type="checkbox"/> 1361T	Extra Tall and Wide MPS Seat Back - 22"-23" Wide - 26" Tall	\$500.00	K0108
<input type="checkbox"/> 1361V	Vinyl Standard MPS Seat Back - 22"-23" Wide - 22" Tall	\$300.00	K0108
<input type="checkbox"/> 1361VT	Vinyl Extra Tall and Wide MPS Seat Back - 22"-23" Wide - 26" Tall	\$700.00	K0108

Rehab Options

<input type="checkbox"/> 1367	Standard Rehab Back	N/C	
<input type="checkbox"/> 3086	Flip Forward Rehab Back - Quick release for folding back canes forward	\$900.00	
<input type="checkbox"/> 1560	Manual Recline Rehab Back	\$825.00	K0108

REHAB SEATING OPTIONS

Part Number	Description	Price	HCPSC Code
<input type="checkbox"/> CU-FN	Curve Cushion - Molded foam w/ medial & lateral thigh supports - size: ___ W x ___ L	\$86.00	E2601
<input type="checkbox"/> 46G	Elements Cushion - Single layer foam with quadra gel pack - size: ___ W x ___ L	\$86.00	E2601
<input type="checkbox"/> M2-F	M2 Cushion-Quadra gel with pelvic, lateral & medial thigh supports - size: ___ W x ___ L	\$305.00	E2607
<input type="checkbox"/> BK-M	ACTA Back - with moldable center section and high resiliency wings	\$640.00	E2613
<input type="checkbox"/> BK-R	M2 Cushion-Quadra gel with pelvic, lateral & medial thigh supports - size: ___ W x ___ L Radius Back - with high resiliency foam	\$636.00	E2611
<input type="checkbox"/> BK-V	circle height: 12" 14" 16" 18" 20" 22" Visco Back - with high resiliency foam	\$578.00	E2611
	circle height: 12" 14" 16" 18" 20" 22"		

HEADREST OPTIONS

Part Number	Description	Price	HCPSC Code
<input type="checkbox"/> 1370	MPS Headrest - for MPS seat back only	\$295.00	
<input type="checkbox"/> HD-MF	Comfort Company Moldable Headrest Pad	\$250.00	E0955
<input type="checkbox"/> HD-BH	Removable Multi-Axis Headrest Hardware	\$250.00	E1028

TIRES

Part Number	Description	Price	HCPSC Code
Drive Tires			
<input type="checkbox"/> TS4	Tire Sealant - for all four	\$60.00	
Spare Tires			
<input type="checkbox"/> 1575C	All Terrain Tire Complete w/ Rim - each <input type="checkbox"/> Single <input type="checkbox"/> Pair	\$210.00	
<input type="checkbox"/> TS	Tire Sealant - per tire	\$15.00	
<input type="checkbox"/> 1595	Jack - power wheelchair tire jack	\$150.00	A9207

ARMREST STYLE OPTIONS

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 1386/87	Standard height adjustable, removable armrest	N/C	K0017
<input type="checkbox"/> FUA	Flip-Up Armrest	\$650.00	
<input type="checkbox"/> PT0031	Parallelogram Reclining Flip Up	\$750.00	
<input type="checkbox"/> 3072/3073	Heavy Duty Flip Back Articulating - height and angle adjustable	\$875.00/ea.	

ARMREST PAD OPTIONS

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> AP3	3" x 12" Desk length leatherette	N/C	K0020
<input type="checkbox"/> AP4	4" x 16" Tapered leatherette	\$40.00/ea.	K0020
<input type="checkbox"/> GAP	Trough style armrest- full length w/hand pad opposite joystick <input type="checkbox"/> Right <input type="checkbox"/> Left	\$185.00/ea	K0106
<input type="checkbox"/> GAPS	Trough style armrest- short without hand pad joystick side <input type="checkbox"/> Right <input type="checkbox"/> Left	\$185.00/ea.	K0106
<input type="checkbox"/> AP75PR	Magic Mobility Arm Pad - 3" x 12" rounded edges	\$105.00/ea.	K0020
<input type="checkbox"/> AP7550PR	Magic Mobility Arm Pad - 3" x 14" rounded edges	\$115.00/ea.	K0020
<input type="checkbox"/> Custom	Custom Arm Pads: _____ _____	RQ	

LEGREST OPTIONS

Part Number	Description	Price	HCPCS Code
<u>NonPower Legrests</u>			
<input type="checkbox"/> 1347	Center mount, fixed, one piece, height adjustable	N/C	K0040
<input type="checkbox"/> 1548	Center mount, flip-up, angle adjustable, one piece, height adjustable	\$250.00	K0040
<input type="checkbox"/> 1355	Center mount, flip-up, angle adjustable, two piece, height adjustable	\$265.00	K0040
<input type="checkbox"/> 1371	Side mount, flip-up, two piece, height adjustable	\$280.00	K0040
<input type="checkbox"/> AMP	Amputee support <input type="checkbox"/> Right <input type="checkbox"/> Left	\$355.00/ea	K0108
<u>Power Legrests</u>			
<input type="checkbox"/> 1886	Side Mount Power Elevating Leg Rest- Swingaway <input type="checkbox"/> Independent <input type="checkbox"/> Combined	\$3,850.00	E1010
<input type="checkbox"/> 2232PG	Center Mount Power Elevating Articulating Legrest	\$3,800.00	K0108
<u>Power Footplate Options</u> - For 2232PG			
<input type="checkbox"/> 22578	Standard Two Piece footplate - 5"x 10.5"	N/C	K0041
<input type="checkbox"/> PT-0109	Extra Large Two Piece Footplate - may not fit with all seat sizes - 5.5" x 11.5"	\$500.00	K0041
<input type="checkbox"/> PT-0117	Extra Large Two Piece Footplate with 2" Heel Guard - may not fit with all seat sizes	\$650.00	K0041
<u>Legrest Accessories</u>			
<input type="checkbox"/> 1550	Fixed Center Mount Calf Pads - NonPowered Legrests	\$430.00	E0995
<input type="checkbox"/> SMCP	Swing Away Side Mount Calf Pads - for 1371 - w/breathable fabric	\$460.00	E0995
<input type="checkbox"/> CMCP	Center Mount Calf Pads - for 2232PG - Fixed	\$450.00	E0995
<input type="checkbox"/> 827A	Foot Straps - pair	\$68.00	E0952
<input type="checkbox"/> HL	Heel Loops	\$95.00	E0951
<input type="checkbox"/> CS	Calf Straps (side mount only)	\$90.00	K0038
<input type="checkbox"/> SP102M	Calf Straps (side mount only) BODYPOINT Aeromesh M-36"	\$90.00	K0038
<input type="checkbox"/> SP102L	Calf Straps (side mount only) BODYPOINT Aeromesh L-40"	\$90.00	K0039

LATERAL SUPPORT OPTIONS

Part Number	Description	Price	HCPCS Code
<u>MPS Seating</u>			
<input type="checkbox"/> LT1-SA	Comfort Company Single Swing away Lateral Hardware <input type="checkbox"/> Right <input type="checkbox"/> Left	\$400.00/ea	E0955
<input type="checkbox"/> LT1-L/RS4W4T	Comfort Company Single Lateral Pad <input type="checkbox"/> Right <input type="checkbox"/> Left	\$53.00/ea	E1028
<u>Rehab Seating</u>			
<input type="checkbox"/> LT1-SA	Comfort Company Single Swing away Lateral Hardware <input type="checkbox"/> Right <input type="checkbox"/> Left	\$263.00/ea	E1028
<input type="checkbox"/> LT1-RGD	Comfort Company Single Rigid Lateral Hardware <input type="checkbox"/> Right <input type="checkbox"/> Left	\$250.00	E1028
<input type="checkbox"/> LT1-RS4W4T	Comfort Company Single Lateral Pad Rigid, 5" W x 4" T, Stretch-Air Cover <input type="checkbox"/> Right <input type="checkbox"/> Left	\$53.00/ea	

POSITIONING OPTIONS

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 1509	Hip Guides - set	\$300.00	E0956
<input type="checkbox"/> I2010	Knee adductors - set, multi-axis, removable with 4.5"x6" pad	\$380.00	K0108

SEAT RESTRAINT OPTIONS: Select Only One

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 1608	Push Button Style Seat Pos. Strap - 64"	N/C	E0978
<input type="checkbox"/> 1608XL	Push Button Style Seat Pos. Strap - 80"	\$50.00	E0978
<input type="checkbox"/> ABSB	Airline Buckle Seat Pos. Strap	\$125.00	E0978
<input type="checkbox"/> EB205L-B1/HW312-	BODYPOINT Evoflex Pelvic Stabilizer L And Mounting Bracket Kit	\$145.00	E0978
<input type="checkbox"/> WC19X86	WC 19 Transit Kit - is approved according to ANSI/RESNA WC19 - includes lap belt 600 mm	\$1,250.00	
<input type="checkbox"/> WC19X87	WC 19 Transit Kit - is approved according to ANSI/RESNA WC19 - includes lap belt 700 mm	\$1,250.00	

CHEST HARNESS OPTIONS

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 1401	Loop Style Positioning Chest Strap	\$105.00	E0906
<input type="checkbox"/> SH103M	Chest Strap BODYPOINT Monoflex, Rear-Pull M	\$110.00	E0906
<input type="checkbox"/> PBH	Padded Butterfly Positioning Chest Strap	\$145.00	E0906
<input type="checkbox"/> SH331L-B1	BODYPOINT Stayflex, Butterfly, Narrow Front-Pull L without zipper	\$168.00	E0906

JOYSTICK OPTIONS

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> PG-JSM	Standard Rnet Joystick Module - No LCD display or mono jacks	N/C	E2323
<input type="checkbox"/> PG-JSM509	Color LCD Joystick Display - with mono jack ports for external switches	\$450.00	E2323

JOYSTICK HANDLE TYPE:

<input type="checkbox"/> RJK	Cone-Shaped - Standard	N/C	E2323
<input type="checkbox"/> 1502	Goal Post	\$140.00	E2323
<input type="checkbox"/> PC102A	BodyPoint U Shaped w/Flex Shaft 3" 4"	\$143.00	E2323
<input type="checkbox"/> PC108A	BodyPoint Mushroom Head Rubber Grip Handle 2" 2.5"	\$134.00	E2323

JOYSTICK MOUNT:

<input type="checkbox"/> 1490R	Standard Fixed - Depth Adjustable Right Mount	N/C	
<input type="checkbox"/> 1490L	Standard Fixed - Depth Adjustable Left Mount	N/C	
<input type="checkbox"/> 1500P	Swing Away Joystick - Right	\$480.00	E1028
<input type="checkbox"/> 1501P	Swing Away Joystick - Left	\$480.00	E1028
<input type="checkbox"/> MM	Stealth Midline Mount	\$568.00	E1028
	<input type="checkbox"/> Right <input type="checkbox"/> Left		

JOYSTICK ACCESSORIES:

<input type="checkbox"/> PG-AM	Attendant Control Joystick	\$1,100.00	E2331
<input type="checkbox"/> TBB	Tash Buddy Button	\$145.00	K0108
<input type="checkbox"/> ESBLK	Egg Switch	\$110.00	
<input type="checkbox"/> CPS9	Swing Away Egg Switch Hardware - Stealth	\$400.00	
<input type="checkbox"/> JSLK	Joystick Lock- Electronically locks chair from driving	\$175.00	
<input type="checkbox"/> RVAL	Reverse Alarm - Beeps when backing up	\$175.00	

ALTERNATIVE DRIVE SYSTEMS:

<input type="checkbox"/> IDH300-1	Stealth i-Drive Head Array	\$5,175.00	E0955
<input type="checkbox"/> F8675	Stealth Mounting Bracket and Egg Switch	\$305.00	E2330
<input type="checkbox"/> TWB480-17S	Stealth Headrest Hardware	\$270.00	E1028
<input type="checkbox"/> HDR-4M-S-C	Switch-It Head Array - Standard three switch head array	\$4,650.00	
<input type="checkbox"/> PASL 105	ASL Proximity Head Array - ASL 300 Egg Switch	\$5,775.00	E2330
<input type="checkbox"/> PASL 606	ASL Adjustable Removable Mounting Hardware	\$215.00	E1028
<input type="checkbox"/> PASL 611	ASL Adjustable Removable Mounting Bracket	\$215.00	E1028
<input type="checkbox"/> PG-CK	Chin Control - w/bib mount and adapter	\$2,232.00	
<input type="checkbox"/> SNP	Sip and Puff	\$1,300.00	E2325
<input type="checkbox"/> PG-OMNI	Omni Display - 4" Color Display For Alternative Drive Systems	\$1,500.00	
	<input type="checkbox"/> Right <input type="checkbox"/> Left		
<input type="checkbox"/> ARMS260	Omni adjustable swing away hardware - Stealth	\$295.00	E1028

ENVIRONMENTAL CONTROLS:

<input type="checkbox"/> PG-BTM2	R-Net Bluetooth Module 2	\$944.00	
----------------------------------	--------------------------	----------	--

NOTES:

ACCESSORIES

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 1605	Waterproof Wheelchair Cover	\$125.00	A9270
<input type="checkbox"/> 1606	Waterproof Joystick Cover	\$40.00	A9270
<input type="checkbox"/> 1603	Medical Necessity Bag - Velcro's to armrest	\$30.00	A9270
<input type="checkbox"/> 1604	Deluxe Backpack	\$75.00	A9270
<input type="checkbox"/> PT0052	Cup Holder - self leveling	\$100.00	A9270
<input type="checkbox"/> A02H	Articulating Oxygen Holder	\$800.00	E2208
<input type="checkbox"/> PT0035	Cane Holder - May not work w/some options	\$350.00	E2207
<input type="checkbox"/> 1396	Removable Clear Polycarb Table - May not work w/ some options	\$525.00	E0980
<input type="checkbox"/> 1137	Luggage Rack	\$270.00	A9270
<input type="checkbox"/> 2855	Lighting Kit Mounted on Seat- 2 headlights, 2 Red Marker Lights w/Turn Signals	\$995.00	A9270
<input type="checkbox"/> PT0036	Lift Brackets- Used with 4-Point Boom Lift	\$375.00	A9270
<input type="checkbox"/> PWB	Power Buddy- Includes USB Charging Port and Reading Light	\$120.00	A9270
<input type="checkbox"/> PT0133	Cell Phone Holder	\$225.00	A9270
<input type="checkbox"/> WCR630	6' Multifolding Ramp - Traction Tape	\$535.00	
<input type="checkbox"/> ADD	Hand Held Programmer	\$1,100.00	
<input type="checkbox"/> LGSVC	Lift Gate Service - For facilities without a loading dock	\$105.00	

*Discounts do not apply

NOTES:

Note: All specification and prices are subject to change without notice. Prices listed are only valid if ordering with a wheelchair. The HCPCS codes provided are not intended to be billing or legal advice, rather Innovation In Motion's interpretation of the code definitions. Use of the codes provided does not ensure coverage or payment for the item. For coverage information, verify the policy of the appropriate payer. It is recommended to contact Innovation In Motion with any questions regarding proper ordering of the chair at 1-800-327-0681.