



FEATURES



Manual Recline



Fixed Arms



3" Dual Casters with Rear-Central Lock and Central Front Brake




500 lb.. Weight Capacity




Heat + Massage

SPACE REQUIREMENTS


TOP VIEW




SIDE VIEW



FRONT VIEW



INFINITE RECLINE POSITIONS



SPECIFICATIONS

Manual recline	Yes
Weight capacity	500 lb.
Back height (from floor/from seat)	45.75"/28.5"
Overall width (no tables)	34.5"
Overall width (2 optional tables down/up)	38.5"/56.5"
Seat width	26"
Seat height/depth	19.75"/20.75"
Arm height (from floor, highest point)	26.75"
Overall upright depth (with handle)	40"
Maximum recline depth	76"
Dual thermoplastic casters (diameter)	3"

ACCESSORIES



Arm Caps




Side Tables



Headrest Cover & Pillow Set



Head Pillow & Lumbar Support



IV Pole