

# The Drop Arm Care Cliner

champion.



## FEATURES



Manual Recline



Fixed Padded Drop-Down Armrests



5" Nylon or Steel Casters



Standard & XL Sizes



Heat + Massage

## SPECIFICATIONS

	Standard	XL
Seat Height	21" (53.3 cm)	21" (53.3 cm)
Seat Width	21" (53.3 cm)	25" (63.5 cm)
Seat Depth	20" (51 cm)	20" (51 cm)
Overall Height	49" (124.5 cm)	49" (124.5 cm)
Overall Width*	33.75" (85.73 cm)	37.75" (95.9 cm)
Back Height above Seat	31.5" (80 cm)	31.5" (80 cm)
Overall Length (Upright)	40" (101.6 cm)	43" (109.2 cm)
Overall Length (Reclined)	76.5" (194.3 cm)	76.5" (194.3 cm)
Armrest Height	27" (68.6 cm)	27" (68.6 cm)
Weight	105 lb. (47 kg)	111 lb. (50.4 kg)
Weight Capacity	275 lb. (124.7 kg)	450 lb. (204 kg)

\*Measured without tables and with arms closed.

## ACCESSORIES



Side Tables




Push Handle



Leg Extensions



Headrest Protector



Head Pillow & Lumbar Support