


		COLOURS INFORMATION For Office Use Only	
Dealer:		Order Date:	S/N:
Account:	P.O.#:	RTS	
Ordered By:		RTS:	Phone:
Address:		Email:	
City:	State: Zip:	SHIP TO	
Phone:	Fax:	Name:	Ship Via:
Email:		Address:	
Customer Name:		City:	State: Zip:
Customer Weight:		Phone:	Fax:

CHAIR INFORMATION	PRICE	
GENESIS	SABER FOLDING	BASE PRICE
Folding, Adjustable FSH-RSH		\$ 2,780
Max. Capacity 265Lbs(120kg)	Total Section A	
Main Frame : ZK60 High Grade A7000 Series Aluminum	Total Section B	
Main Frame : Lifetime Warranty *See warranty policy	Total Section C	
	Total Section D	
	Total Section E	
	ORDER GRAND TOTAL	

A FRAME COLORS & PARTS COLOR)


FRAME COLORS

Powder-Coated	☐ Gloss Black	☐ Twilight Flake	☐ Yellow	☐ Magenta	☐ Fire Engine Red	☐ White	\$ STD
	☐ Candy Blue	☐ Candy	☐ Candy Green	☐ Candy Red	☐ Candy Raspberry	☐ Cobalt Blue	
	☐ Matte Gray	☐ Matte Purple	☐ Matte Blue	☐ Matte Green	☐ Matte Black	☐ Matte Red	
	☐ Neon Red	☐ Neon Green	☐ Orange Sparkle	☐ Mellon Yellow	☐ Gun Smoke	☐ Pigeon Blue	
	☐ Granny Sparkle Green	☐ Sparkle Sky Blue	☐ Vermilion	☐ Reed Green	☐ Blue Lilac		
Special Color	☐ Polished Aluminum Color					\$ 800	



Total Section A \$

B FRAME MEASUREMENTS 1/2



1.Seat Width ☐ 12.5"(32.0cm) ☐ 13 3/8"(34.0cm) ☐ 14 1/8"(36.0cm) ☐ 15.0"(38.0cm) ☐ 15 3/4"(40.0cm)
☐ 16.5"(42.0cm) ☐ 17 1/4"(44.0cm) ☐ 18 1/8"(46.0cm)

\$ STD

2.Frame Depth &
3.Seat Depth

		Seat Depth				
Frame Depth	STD Frame 20"(51cm)	☐ 13 3/8"(34.0cm)	☐ 14 1/8"(36.0cm)	☐ 15.0"(38.0cm)	☐ 15 3/4"(40.0cm)	
	Long Frame 22 3/4"(58cm)	☐ 15 3/4"(40.0cm)	☐ 16.5"(42.0cm)	☐ 17 1/4"(44.0cm)	☐ 18 1/8"(46.0cm)	☐ 18 7/8"(48cm)

\$ STD

\$ STD

4.Footrest Style &
Footrest Width

Footrest Style : Center Joint Folding Foot Plate only Adj. Angle & Height
Footrest Width : Your seat width minus 3.5" only

\$ STD

		Your Seat Width						
Footrest Width		12.5"(32.0cm)	13 3/8"(34.0cm)	14 1/8"(36.0cm)	15.0"(38.0cm)	15 3/4"(40.0cm)	16.5"(42.0cm)	17 1/4"(44.0cm)
		18 1/8"(46.0cm)						
	9.0"(23.0cm)	9 7/8"(25.0cm)	10 5/8"(27.0cm)	11 3/8"(29.0cm)	12 1/8"(31.0cm)	13"(33.0cm)	13 3/4"(35.0cm)	14.5"(37.0cm)

5.Frame Bend Angle

95° Only

\$ STD

Frame Angle to the Ground

		Your Seat Depth							
Seat Bucket		13 3/8"(34.0cm)	14 1/8"(36.0cm)	15.0"(38.0cm)	15 3/4"(40.0cm)	16.5"(42.0cm)	17 1/4"(44.0cm)	18 1/8"(46.0cm)	18 7/8"(48.0cm)
	0"(0cm)	95°	95°	95°	95°	95°	95°	95°	95°
	0.5"(1.2cm)	93°	93°	93°	93°	93°	93°	93°	93°
	1.0"(2.5cm)	91°	91°	91°	91°	92°	92°	92°	92°
	1.5"(3.8cm)	89°	89°	89°	90°	90°	90°	90°	90°
	2.0"(5.1cm)	86°	87°	87°	88°	88°	88°	89°	89°
	2.5"(6.4cm)	84°	85°	85°	86°	86°	87°	87°	87°
	3.0"(7.6cm)	82°	83°	83°	84°	85°	85°	85°	86°
	3.5"(8.9cm)	80°	81°	81°	82°	83°	83°	84°	84°
	4.0"(10.2cm)	78°	79°	79°	80°	81°	82°	82°	83°

6.Backrest Height

☐ 12.5"(32.0cm) ☐ 13 3/8"(34.0cm) ☐ 14 1/8"(36.0cm)

Adj. range

\$ STD

☐ 15.0"(38.0cm) ☐ 15 3/4"(40.0cm) ☐ 16.5"(42.0cm)

Adj. range

7.Backrest Angle

☐ Fixed Backrest 87° from seat sling

\$ STD

NOTE

Sub Total Section B 1/2 \$

B FRAME MEASUREMENTS 2/2

8. Front Seat Height

- ☐ 16.5"(42.0cm) *1 ☐ 17.0"(43.0cm) *1 ☐ 17 1/4"(44.0cm) ☐ 17 3/4"(45.0cm)
☐ 18 1/8"(46.0cm) ☐ 18 1/2"(47.0cm) ☐ 18 7/8"(48.0cm) *2 ☐ 19 1/4"(49.0cm) *2
☐ 19 5/8"(50.0cm) *2
***1: STD Frame Only**
***2: Long Frame Only**

\$ STD

Front Seat height Range

		Front Seat Height								
		16.5"(42.0cm)	17.0"(43.0cm)	17 1/4"(44.0cm)	17 3/4"(45.0cm)	18 1/8"(46.0cm)	18 1/2"(47.0cm)	18 7/8"(48.0cm)	19 1/4"(49.0cm)	19 5/8"(50.0cm)
Caster & Forks	3" Caster & STD Forks	☐ STD FRAME ONLY	☐ STD FRAME ONLY	☐	☐	☐	☐ LONG FRAME ONLY	☐ LONG FRAME ONLY		
	4" Caster & STD Forks		☐ STD FRAME ONLY	☐ STD FRAME ONLY	☐	☐	☐ LONG FRAME ONLY	☐ LONG FRAME ONLY	☐ LONG FRAME ONLY	
	5" Caster & STD Forks			☐ STD FRAME ONLY	☐ STD FRAME ONLY	☐	☐ LONG FRAME ONLY	☐ LONG FRAME ONLY	☐ LONG FRAME ONLY	☐ LONG FRAME ONLY

9. Rear Seat Height

- ☐ 14.5"(37.0cm) ☐ 15.0"(38.0cm) ☐ 15 3/8"(39.0cm) ☐ 15 3/4"(40.0cm) ☐ 16 1/8"(41.0cm)
☐ 16.5"(42.0cm) ☐ 17.0"(43.0cm) ☐ 17 1/4"(44.0cm) ☐ 17 3/4"(45.0cm)

\$ STD

Rear Seat height Range

		Rear Seat Height								
		14.5"(37.0cm)	15.0"(38.0cm)	15 3/8"(39.0cm)	15 3/4"(40.0cm)	16 1/8"(41.0cm)	16.5"(42.0cm)	17.0"(43.0cm)	17 1/4"(44.0cm)	17 3/4"(45.0cm)
Rear Tire Size	24"x1"(540)	☐	☐	☐	☐	☐	☐			
	24"x1 3/8"(540)		☐	☐	☐	☐	☐			
	25"x1"(559)		☐	☐	☐	☐	☐	☐		
	24"x2"(540)			☐	☐	☐	☐	☐		
	26"x1"(590)			☐	☐	☐	☐	☐	☐	
	25"x1.95"(559)				☐	☐	☐	☐	☐	☐

10. Center of Gravity

- ☐ 1 1/8"(3.0cm) ☐ 1.5"(4.0cm) ☐ 2.0 (5.0cm) ☐ 2 3/8"(6.0cm)
☐ 2 3/4"(7.0cm) ☐ 3 1/8"(8.0cm)

\$ STD

11. Camber Angle


- ☐ 0° ☐ 2° ☐ 4°

\$ STD

12. Seat to Footrest

()
 *Adj. range : 13.0" - 17.5" (33cm-45cm)

*Measure vertical from top of seat sling at front edge to top of foot plate.



\$ STD

13. Inside Wheel Clearance

- ☐ 3/4"(2.0cm) ☐ 1.0"(2.5cm) ☐ 1 1/4"(3.0cm)

\$ STD

NOTE

Sub Total Section B 1/2 \$

Sub Total Section B 2/2 \$

Total Section B \$

WHEELS

GENESIS: Page 4 of 7

C WHEELS, TIRES & CASTERS, FORKS 2/2
TIRES

☐ High Pressure Tire ☐ 24"x1", ☐ 25"x1" or ☐ 26"x1" Gray Tire with black wall Only	E2211/E2212	\$ STD
☐ High Pressure Tire ☐ 24"x1", ☐ 25"x1" or ☐ 26"x1" ☐ Blue ☐ Red ☐ Yellow	E2211/E2212	\$ 85
☐ Low Pressure Tire ☐ 24"x1-3/8" or ☐ 26"x1-3/8" Gray Only	E2211/E2212	\$ STD
☐ E2211Schwalbe Marathon Plus High Pressure Tire ☐ 24"x1", ☐ 25"x1" or ☐ 26"x1" Black Only	E2211/E2212	\$ 99
☐ Solid Tire ☐ 24"x1", ☐ 25"x1" or ☐ 26"x1" ☐ Dark Gray ☐ Black(\$30up charge)	E2218	\$ 115
☐ E2211 Moderate V-Track Knobby Tire 24"x1-3/8" (37-540) Gray Only	E2211/E2212	\$ 117
☐ E2211 Aggressive Knobby Tire 24"x 2" (50-540) Black Only	E2211/E2212	\$ 117
☐ E2211 Aggressive Knobby Tire 26"x 1.95" (50-559) Black Only	E2211/E2212	\$ 117
☐ E2213 Flat-Free Inserts For 24"x1-3/8" Tires Only	E2211/E2213	\$ 117

CASTERS

SIZE	TIRE WIDTH	HUB COLOR	TIRE COLOR	DESCRIPTION	
☐ 3"	1"	Black Hub	Black Tire Only	Polyeurathane	\$ STD
☐ 3"	1"	Plastic Black Hub	Clear Tire Only	Lighted, Polyeurathane	\$ 65
☐ 3"	1"	Red Hub	Red Tire Only	Polyeurathane	\$ STD
☐ 3"	1"	Black Hub	Clear Tire Only	Polyeurathane	\$ STD
☐ 3"	1"	Aluminum Silver Hub	Black Tire Only	Made By Frog Legs	\$ 165
☐ 3"	1"	Aluminum Black Hub	Black Tire Only	Made By Frog Legs	\$ 165
☐ 3"	1.5"	Aluminum Silver Hub	Black Tire Only	Soft Roll Made By Frog Legs	\$ 165
☐ 3"	1.5"	Aluminum Black Hub	Black Tire Only	Soft Roll Made By Frog Legs	\$ 165
☐ 4"	1.5"	Plastic Black Hub	Black Tire Only		\$ STD
☐ 4"	1"	Plastic Black Hub	Clear Only	Lighted, Polyeurathane	\$ 65
☐ 4"	1"	Aluminum Silver Hub	Black Tire Only	Made By Frog Legs	\$ 165
☐ 4"	1"	Aluminum Black Hub	Black Tire Only	Made By Frog Legs	\$ 165
☐ 4"	1.5"	Aluminum Silver Hub	Tire ☐ Black ☐ Blue ☐ Red	Soft Roll Made By Frog Legs	\$ 165
☐ 4"	1.5"	Aluminum Black Hub	Black Tire Only	Soft Roll Made By Frog Legs	\$ 165
☐ 5"	1"	Composite Black Hub	Black Tire Only		\$ STD
☐ 5"	1"	Plastic Black Hub	Clear Tire Only	Lighted, Polyeurathane	\$ 65
☐ 5"	1"	Aluminum Silver Hub	Black Tire Only	Made By Frog Legs	\$ 165
☐ 5"	1"	Aluminum Black Hub	Black Tire Only	Made By Frog Legs	\$ 165
☐ 5"	1.5"	Aluminum Silver Hub	Tire ☐ Black ☐ Blue ☐ Red	Soft Roll Made By Frog Legs	\$ 165
☐ 5"	1.5"	Aluminum Black Hub	Black Tire Only	Soft Roll Made By Frog Legs	\$ 165

FORKS


☐ STD Aluminum Forks Black Only	\$ STD
---	--------

Sub Total Section C 1/2 \$

Sub Total Section C 2/2 \$

Total Section C \$

D COMPONENTS

WHEEL LOCKS										
☐ Push to Lock / Anodized Black Only										\$ STD
☐ NEW Compact Push to Lock / Anodized Black Only ☐ w/ 2" Extension Handle (\$65up)										\$ 110
☐ NEW Compact Scissor Lock / Anodized Black Only										\$ 110
PUSH HANDLES										
☐ Integrated w/Fixed Backrest Height / Anodized Black Only										\$ STD
☐ Omit Push Handle										\$ NCO
NEW ARM RESTS Anodized Black Only										
☐ Fixed Tubular Armrest Adj. Angle and Height Height from seat sling Set at : ☐ 8" ☐ 9" ☐ 10" ☐ 11" ☐ 12" E0973										\$ 250
☐ Swing Away Tubular Armrest Adj. Angle and Height Height from seat sling Set at : ☐ 8" ☐ 9" ☐ 10" ☐ 11" ☐ 12" E0973										\$ 250
☐ Flip-Up Armrest ☐ w/Fixed Tubular Armrest Height from seat sling Set at : ☐ 8" ☐ 9" ☐ 10" ☐ 11" ☐ 12" E0973										\$ 250
Adj. Angle ☐ w/ Swing Away Tubular Armrest Height from seat sling Set at : ☐ 8" ☐ 9" ☐ 10" ☐ 11" ☐ 12" E0973										\$ 250
☐ Custom Length (\$150 up charge) *Specify Length 										
ANTI-TIPPERS										
☐ Anti-Tippers (pair) w/ Quick Release / Anodized Black Only E0971										\$ 112
☐ Folding Anti-Tippers (pair) w/ Quick Release / Black Only E0971										\$ 358
SIDE GUARDS										
☐ Adj. Height Rigid Side Guards - Black ABS Non-Removable Style Only										\$ 258
SEAT CUSHION -*Please see section E for other cushion choices.-										
☐ Foam Cushion w/Nylon Black Cover Thickness ☐ 1" ☐ 2" ☐ 3" ☐ 4"										\$ 78
☐ Custom Cushion Cover (Std. Foam Cushion Only) -Select Fabric () (\$90 up charge)										
MISCELLANEOUS										
☐ Frame Protectors / Black Leather/Zipper Type										\$ 78
☐ K0065 Spoke Protectors / Black Only ☐ 24" ☐ 25" ☐ 26" K0065										\$ 107
☐ E0978 Seat Belt w/Auto Buckle (BP) 1.5" Black Belt Only E0978										\$ 69
UPHOLSTERY CHOICE & COLOR *Positioning Back Available. Please see Section E										
☐ Adjustable Tension Backrest Dark Grey Color Only										\$ STD
Custom Options	☐ Fabrics for Fun ☐ 001 ☐ 002 ☐ 003 ☐ 004 ☐ 006 ☐ 010 ☐ 011 . . . (\$95 up charge) ☐ For Std. or Adj. Style ☐ 014 ☐ 015 ☐ 016 ☐ 017 ☐ 018 ☐ 020 ☐ 021 ☐ 022 ☐ 023 ☐ 024 ☐ 025 ☐ 026 ☐ 027 ☐ 028									
☐ Calf Strap ☐ Black ☐ Matching Fabric (\$55 up charge) -Select Fabric() K0038										\$ 24
☐ Matching Fabric on Padded Side Guards -Select Fabric ()										\$ 210
NOTE										
Total Section D \$										

E CHAIR OPTIONS

SEAT CUSHION

SUPRACOR® WHEELCHAIR CUSHIONS

Stimulite® Contoured Cushion *Black Upholstery only	☐ W14"xD16" (CD1416)		E2603	\$ 525
	☐ W16"xD16" (CD1616)		E2603	\$ 525
	☐ W16"xD18" (CD1618)		E2603	\$ 525
	☐ W18"xD16" (CD1816)		E2603	\$ 525
	☐ W18"xD18" (CD1818)		E2603	\$ 525
	☐ W18"xD20" (CD1820)		E2603	\$ 540
	☐ W20"xD16" (CD2016)		E2603	\$ 540
	☐ W20"xD18" (CD2018)		E2603	\$ 540
	☐ W20"xD20" (CD2020)		E2603	\$ 555
	☐ Extra Sensitive (\$30 up charge) ☐ Sling Style Bottom (\$45 up charge) ☐ Custom Size W(") x D(") (\$- CFP)			

Stimulite® Slimline Cushion *Black Upholstery only	☐ W14"xD16" (SM1416)		E2603	\$ 495
	☐ W16"xD16" (SM1616)		E2603	\$ 495
	☐ W16"xD18" (SM1618)		E2603	\$ 495
	☐ W18"xD18" (SM1818)		E2603	\$ 495
	☐ W18"xD20" (SM1820)		E2603	\$ 495
	☐ W20"xD16" (SM2016)		E2603	\$ 510
	☐ W20"xD18" (SM2018)		E2603	\$ 510
	☐ W20"xD20" (SM2020)		E2603	\$ 510
	☐ Extra Sensitive (\$30 up charge) ☐ Sling Style Bottom (\$45 up charge) ☐ Custom Size W(") x D(") (\$- CFP)			

Stimulite® Corbee Pediatric *Black Upholstery only	☐ W10"xD10" (PD1010)		E2603	\$ 450
	☐ W12"xD12" (PD1212)		E2603	\$ 450
	☐ W14"xD14" (PD1414)		E2603	\$ 450
	☐ Extra Sensitive (\$30 up charge) ☐ Sling Style Bottom (\$45 up charge) ☐ Custom Size W(") x D(") (\$- CFP)			

ROHO WHEELCHAIR CUSHIONS

☐ 4" HIGH PROFILE® QUADTRO SELECT® CUSHION				*Black Upholstery				E2622-2623	\$ 525
Choose Size	☐ W12"xD12"	☐ W12"xD13"	☐ W13"xD13"	☐ W13"xD15"	☐ W13"xD16"	☐ W15"xD13"	☐ W15"xD15"		
	☐ W15"xD16"	☐ W15"xD18"	☐ W16"xD15"	☐ W16"xD16"	☐ W16"xD18"	☐ W16"xD20"	☐ W18"xD15"		
	☐ W18"xD16"	☐ W18"xD18"	☐ W18"xD20"	☐ W18"xD22"	☐ W20"xD16"	☐ W20"xD18"	☐ W20"xD20"		
	☐ W20"xD22"	☐ W22"xD18"	☐ W22"xD20"	☐ W22"xD22"	☐ Custom Size W(") x D(") (\$+85)				

☐ 3" MID PROFILE® QUADTRO SELECT® CUSHION				*Black Upholstery only		E2622-2623	\$ 525
Choose Size	☐ W15"xD15"	☐ W15"xD16"	☐ W16"xD16"	☐ W16"xD18"	☐ W18"xD16"	☐ W18"xD18"	
	☐ Custom Size W(") x D(") (+\$85)						

☐ 2" LOW PROFILE® QUADTRO SELECT® CUSHION				*Black Upholstery only		E2622-2623	\$ 525
Choose Size	☐ W12"xD12"	☐ W12"xD13"	☐ W13"xD13"	☐ W13"xD15"	☐ W13"xD16"	☐ W15"xD13"	☐ W15"xD15"
	☐ W15"xD16"	☐ W15"xD18"	☐ W16"xD15"	☐ W16"xD16"	☐ W16"xD18"	☐ W16"xD20"	☐ W18"xD15"
	☐ W18"xD16"	☐ W18"xD18"	☐ W18"xD20"	☐ W20"xD16"	☐ W20"xD18"	☐ W20"xD20"	
	☐ Custom Size W(") x D(") (+\$85)						

☐ 2" NEXUS SPIRIT® CUSHION				*Black Upholstery only		E2622	\$ 525
Choose Size	☐ W14"xD14"	☐ W14"xD16"	☐ W15"xD15"	☐ W16"xD16"	☐ W16"xD18"	☐ W17"xD17"	☐ W18"xD16"
	☐ W18"xD18"	☐ W20"xD18"	* Custom Size Not Available				

STAR CUSHIONS

Galaxy Cushion *Black Upholstery only			Choose Zoning Options		E2624	\$ 498
☐ W14"xD14" (GXY 1414-2)	☐ W14"xD16" (GXY 1416-2)		☐ Standard - Dual Valve			
☐ W16"xD14" (GXY 1614-2)	☐ W16"xD16" (GXY 1616-2)		☐ Dual Valve - Front to Back			
☐ W16"xD18" (GXY 1618-2)	☐ W18"xD16" (GXY 1816-2)		☐ Four Valve			
☐ W18"xD18" (GXY 1818-2)	☐ W18"xD20" (GXY 1820-2)		☐ Simple - Single Valve			
☐ W20"xD16" (GXY 2016-2)	☐ W20"xD18" (GXY 2018-2)		☐ Four Valve - Enhanced			
☐ W20"xD20" (GXY 2020-2)			☐ Four Valve - Front to Back			

THE CURVE CUSHION *Black Upholstery only				☐ W10"xD10"	☐ W10"xD11"	☐ W10"xD11"	☐ W10"xD12"	☐ W11"xD12"	☐ W12"xD12"	\$ 133
				☐ W12"xD13"	☐ W12"xD14"	☐ W13"xD11"	☐ W13"xD12"	☐ W13"xD13"	☐ W13"xD14"	
				☐ W14"xD14"	☐ W14"xD15"	☐ W14"xD16"	☐ W15"xD14"	☐ W15"xD15"	☐ W15"xD16"	
				☐ W16"xD14"	☐ W16"xD15"	☐ W16"xD16"	☐ W16"xD17"	☐ W16"xD18"	☐ W17"xD16"	
				☐ W17"xD17"	☐ W17"xD18"	☐ W18"xD16"	☐ W18"xD17"	☐ W18"xD18"		

NOTE