## AQUATIC LIFTS



# Lolo<sup>™</sup> & Gallatin<sup>™</sup> WP 400





The Lolo and Gallatin Pool Lifts are ADA compliant Water Powered pool lifts that are self-operable from the deck and water. The Lolo is designed for in-ground pools while the Gallatin can accommodate more shallow applications and spas. Both will provide years of low maintenance pool access.





Gallatin



## Accessories

- Booster pump
- Protective cover
- Spineboard attachment
- Pre-construction and retro-fit anchors
- Waterbox

#### **Features**

- 400 pound lifting capacity with
  55 psi water pressure
- Quiet and stable operation
- Spacious seat with excellent lumbar support and ease of transfer
- Stainless steel construction (electropolished 304L grade)
- Accommodates setback of 18" to 30" and water draft up to 12"
- Easy to install, even into existing concrete decks and anchors

### Includes

- 5 year structural warranty
- Spectra Shield®coating for maximum corrosion protection
- Built in wheels for portability
- Adjustable swing out footrest
- Flip-up armrests with non-slip grip
- Padded headrest
- Seat belt