

188 EZEE LIFE™ MODEL 188 – ALUMINUM BARIATRIC SHOWER & COMMODE CHAIR, 600 LB CAPACITY

The EZee Life™ Model 188 is a lightweight, rust-resistant 2-in-1 bariatric shower and commode chair with a standard seat with a 24" internal width between the arms. Its durable aluminum frame is ideal for wet environments and provides reliable support for larger users. These aluminum shower/commode chairs facilitate smooth client movement, enabling care providers to handle heavier individuals with ease.

Key Safety Features:

- Flip-back padded armrests with red safety lockouts to help prevent accidental client falls and allow for safe and easy side transfers.
- Anti-tippers for enhanced rear stability, reducing the risk of tipping and improving overall safety.
- Patient safety belt to ensure proper positioning, support, and added security during use.
- Calf strap to support the legs and help maintain a stable 90-degree seated position, promoting comfort and posture.

Features:

- Soft to touch molded seat, back & arm rest pads
- Commode hole insert
- Flip-back, half-length padded armrest for comfort with red button lockout for client safety
- Locking swivel casters
- Adjustable height swing away footrests
- Flip-up composite footplates are soft to the touch and stay warm, enhancing client comfort
- Push button seat belt
- Strong tubular aluminum frame: Heavy-duty reinforced lightweight aluminum chair designed for bariatric clients
- Adjustable removable anti-tippers for client safety



Specifications:

Width between the Arms:	25" (63.5 cm)
Casters:	5" twin wheel casters
Seat Width:	22" (55.88 cm)
Seat Depth:	20" (50.8 cm)
Seat Height:	19", 20", 21" (adjustable)
Seat to Floor Clearance:	16.3", 17.3", 18.3" (41.4 x 43.9 x 46.5 cm)
Seat Height Range:	19" x 20" x 21" (48.2 x 50.7 x 53.2 cm)
Back Height:	8" (20 cm)
Arm Height:	11" (30 cm)
Arm Pad Length:	9½" (24cm)
Footrests:	Yes – Swing-away/Removable
Weight Capacity:	600 lbs
Frame Material:	Aluminum
Weight:	43.54 lbs (19.70 kg)
Box Dimensions:	30" x 26" x 17" (75 x 65.5 x 42 cm)

