

**198**

## EZEE LIFE™ MODEL 198 BARIATRIC ALUMINUM TILT SHOWER & COMMODE CHAIR – 440 LB CAPACITY

The EZee Life™ Model 198 is a lightweight, rust-resistant 2-in-1 bariatric tilting shower and commode chair, built to support users up to 200 kg (440 lbs). With a standard seat and a spacious 22" width between the arms, it provides comfort and stability for larger individuals. Its durable aluminum frame is built to withstand wet environments, ensuring reliable long-term use. Designed for ease of transfers, the Model 198 helps caregivers handle heavier users safely and efficiently


### Features:

- Flip-back, half-length padded armrests with safety lockouts
- Adjustable/removable anti-tippers for rear stability
- Push-button safety belt for secure positioning
- Calf strap for leg support
- Adjustable, flexible headrest
- 90–120° tilt range for comfort and posture
- Soft-touch molded seat, back, arms, and footplates
- Warm-retaining, flip-up footplates
- Locking swivel casters for safe, smooth movement
- Height-adjustable casters with locking pin
- Adjustable-height, swing-away footrests
- Commode hole insert and removable pail included
- Lightweight, heavy-duty tubular aluminum frame



**440LB  
CAPACITY**  
22" INTERNAL WIDTH CHAIR

|                          |                                                |
|--------------------------|------------------------------------------------|
| Width between the Arms:  | 22.63" (57.5 cm)                               |
| Casters:                 | 5" swivel locking                              |
| Seat Width:              | 22" (55.88 cm)                                 |
| Seat Depth:              | 20" (50.8 cm)                                  |
| Seat Height:             | 21.85" x 22.83" x 23.82" (adjustable)          |
| Seat to Floor Clearance: | 19.29" x 20.28" x 21.26" (49 x 51.5 x 54 cm)   |
| Seat Height Range:       | 21.85" x 22.83" x 23.82" (55.5 x 58 x 60.5 cm) |
| Back Height:             | 18.7" (47.5 cm)                                |
| Arm Height               | 8" (20 cm)                                     |
| Arm Pad Length:          | 11.81" (30 cm)                                 |
| Footrests:               | Yes – Swing-away/Removable                     |
| Weight Capacity:         | 440 lbs                                        |
| Frame Material:          | Aluminum                                       |
| Weight:                  | 45.63 lbs (20.7 kg)                            |
| Box Dimensions:          | 29" x 26" x 24" (72 x 66 x 61 cm)              |



Release red safety lock handle



Half-Length Padded Armrest in Locked Position



Safety anti-tippers with height adjustment and rubber tip



Half-length padded armrest flipped up