



| SETBACK - W.D. | SETBACK - WATER DRAFT GUIDE | | | | | | | | | | | | |
|------------------|-----------------------------|---------|---------|---------|---------|---------|---------|---------|---------|---------|----------|----------|----------|
| | 0" W.D. | 1" W.D. | 2" W.D. | 3" W.D. | 4" W.D. | 5" W.D. | 6" W.D. | 7" W.D. | 8" W.D. | 9" W.D. | 10" W.D. | 11" W.D. | 12" W.D. |
| 6" - 22" SETBACK | X | X | X | X | X | X | X | X | X | X | X | X | X |
| 23" SETBACK | X | X | X | X | X | X | X | X | X | X | X | X | |
| 24" SETBACK | X | X | X | X | X | X | X | X | X | X | X | | |
| 25" SETBACK | X | X | X | X | X | X | X | X | X | X | | | |
| 26" SETBACK | X | X | X | X | X | X | X | X | X | | | | |
| 27" SETBACK | X | X | X | X | X | X | X | | | | | | |
| 28" SETBACK | X | X | X | X | X | | | | | | | | |
| 29" SETBACK | X | X | X | X | | | | | | | | | |
| 30" SETBACK | X | X | X | | | | | | | | | | |
| 31" SETBACK | X | | | | | | | | | | | | |
| 32" SETBACK | X | | | | | | | | | | | | |
| 33" SETBACK | X | | | | | | | | | | | | |

SPECTRUM
Aquatics

A PLAYCORE Company

7100 Spectrum Lane Missoula, Montana 59808
800.791.8056 www.spectrumaquatics.com

ENGINEERING DATA

| | |
|-----------|-----------------------------------|
| PRODUCT: | HORIZON BP 450 DELUXE |
| MODEL #: | 2010102-DLX |
| MATERIAL: | SEE PARTS LIST |
| REVISION: | A Monday, July 6, 2020 9:11:53 AM |

NOTE:

△ DIMENSIONS VARY WITH SPECIFIC APPLICATIONS. SPECIFICATIONS ARE NOMINAL AND ARE SUBJECT TO CHANGE. PLEASE CONTACT SPECTRUM AQUATICS FOR CUSTOM APPLICATIONS

◆ 450 LB LIFTING CAPACITY ◆ TWO WATERTIGHT REMOTE CONTROLS ◆ ADJUSTABLE FOOTREST FOR EASE OF TRANSFER ◆ FLIP UP ARM REST WITH NON-SLIP GRIP FOR SAFETY AND COMFORT ◆ ACCOMMODATES SETBACK OF 6" TO 33" AND WATER DRAFT UP TO 12"