



	SETBACK - WATER DRAFT GUIDE												
SETBACK - W.D.	0" W.D.	1" W.D.	2" W.D.	3" W.D.	4" W.D.	5" W.D.	6" W.D.	7" W.D.	8" W.D.	9" W.D.	10" W.D.	11" W.D.	12" W.D.
6" - 22" SETBACK	X	X	X	X	X	X	X	X	X	X	X	X	X
23" SETBACK	X	X	X	X	X	X	X	X	X	X	X	X	
24" SETBACK	X	X	X	X	X	X	X	X	X	X	X		
25" SETBACK	X	X	X	X	X	X	X	X	X	X			
26" SETBACK	X	X	X	X	X	X	X	X	X				
27" SETBACK	X	X	X	X	X	X	X						
28" SETBACK	X	X	X	X	X								
29" SETBACK	X	X	X	X									
30" SETBACK	X	X	X										
31" SETBACK	X												
32" SETBACK	X												
33" SETBACK	X												

SPECTRUM
Aquatics

A PLAYCORE Company

7100 Spectrum Lane
800.791.8056

Missoula, Montana 59808
www.spectrumaquatics.com

ENGINEERING DATA

PRODUCT:	HORIZON BP 450
MODEL #:	2010102
MATERIAL:	SEE PARTS LIST
REVISION:	A Monday, July 6, 2020 9:10:25 AM

NOTE:

△ DIMENSIONS VARY WITH SPECIFIC APPLICATIONS. SPECIFICATIONS ARE NOMINAL AND ARE SUBJECT TO CHANGE. PLEASE CONTACT SPECTRUM AQUATICS FOR CUSTOM APPLICATIONS

◆ 450 LB LIFTING CAPACITY ◆ TWO WATERTIGHT REMOTE CONTROLS ◆ ADJUSTABLE FOOTREST FOR EASE OF TRANSFER ◆ FLIP UP ARM REST WITH NON-SLIP GRIP FOR SAFETY AND COMFORT ◆ ACCOMMODATES SETBACK OF 6" TO 33" AND WATER DRAFT UP TO 12"