

Horizon Long Reach[™] BP 300 & Deluxe



The Horizon Long Reach 300 and 300 Deluxe are wireless battery powered lifts, self operable from the deck and water. With a robust operating system, these lifts can accommodate a wide variety of decks and gutter setbacks for both in-ground and above ground pools.



*Anchor sold separately

Includes

- 5 year structural warranty
- 2 year electronics warranty
- 24 volt battery and charger
- Spectra Shield[®] coating for corrosion protection
- 2 wireless remote controls
- Flip-up armrests with non-slip grip
- Padded headrest (Deluxe models)
- Adjustable footrest
- Seat belt

Features

- 300 pound lifting capacity
- Quiet and stable operation
- 90 lifts per battery charge
- Easy to install, even into existing concrete decks and anchors
- Rotationally and vertically powered with 360° continuous power rotation
- Accommodates setback of 6" to 40" and water draft up to 14"
- Stainless steel construction (electropolished 304L grade)
- Spacious seat offers excellent lumbar support and ease of transfer

Accessories

- Protective cover
- Anchor with threaded cap and key
- Extra 24 volt battery and charger
- Extra wireless remote
- Spine board attachment and stability strap
- Powder coating available
- Transport cart