



FEATURES




Manual Recline



Swing-Away Arms



Fixed Front Casters
with Accessible Lock
tabs




500 lb. Weight Capacity




Heat + Massage

SPACE REQUIREMENTS


TOP VIEW




TOP VIEW WITH SWING-AWAY ARMS



FRONT VIEW




SIDE VIEW



SPECIFICATIONS

Manual recline	Yes
Weight capacity	500 lb.
Back height (from floor/from seat)	51.75"/28.5"
Overall width (no tables)	34.5"
Overall width (2 optional tables down/up)	38.5"/56.5"
Seat width	26"
Seat height/depth	25.5"/20.75"
Arm height (from floor)	32.5"
Overall upright depth	41.5"
Trendelenburg	74"
Casters (diameter)	5"

ACCESSORIES



Arm Caps



Side Tables



Foot Tray



Push Handle



Headrest Cover &
Pillow Set