LiftOperator® Intelligent Control



A standard feature with all S.R.Smith ADA compliant pool lifts, the LiftOperator is the most technologically advanced pool lift control on the market.

- Emergency stop button sounds an audible alarm and can stop the lift
- · Bright LED color-coded indicator shows battery life
- Service notification light meets ADA guidelines for scheduled maintenance, promoting regular service checks (based on 3–4 months average usage)
- USB port transfers usage data to a USB memory stick for review on PC


		Shipping			
Model No.	Description	Weight	Length	Width	Height
400-7001	2-Button Control Box	5 lbs 2.5kg	16" 40.6cm	5" 12.7cm	5" 12.7cm
400-7001-BC*	CEC LiftOperator, 2-Button (2018 California compliant)				
400-7000-BC*	CEC LiftOperator, 4-Button (2018 California compliant)				
400-7000	4-Button Control Box				
400-7000AX	aXs2 4-Button Control Box				
400-7000AX-BC	aXs2 4-Button Control Box (2018 California compliant)				

• 2-Button control for multiLift and ML300

Notes:

• 4-Button control for PAL, PAL2, Splash! and aXs2 lifts

Components, Accessories & Upgrade Kits



^{*}As of 1.2.18, all new pool lift orders into California feature LiftOperators with California Energy Commission (CEC) compliant battery charging systems. The 2018 CEC compliant LiftOperator control boxes, batteries and charging systems are not compatible with the standard LiftOperator units. Customers in California will be required to provide the serial number of their pool lift when ordering replacement LiftOperator components. For more information, please visit srsmith.com/CEC2018.

Upgrade Kits — from LINAK to LiftOperator;
Upgrade Kit includes: Mounting Plate, Lock,
Battery, Charger and LiftOperator Control Unit







^{**}Tan covers to be launched in March 2018. Call for availability.