

A PAYCORE Company





The Motion Trek 400 and 400 Deluxe are ADA compliant battery powered lifts that are self operable from the deck and water with a robust operating system. These lifts are designed for inground pools and spas.



*Anchor sold separately

Includes

- 5 year structural warranty
- 2 year electronics warranty
- 24 volt battery and charger
- Spectra Shield®coating for maximum corrosion protection
- Wired handset
- Seatbelt
- Flip-down armrests
- Adjustable padded headrest (Deluxe models)
- Swing-out adjustable footrest (Deluxe models)

Features

- 400 pound lifting capacity
- 90 lifts per battery charge
- Spacious seat with excellent lumbar support and ease of transfer
- Rotationally and vertically powered with 360° continuous power rotation
- Stainless steel construction (electro-polished 304L grade)
- Accommodates setback of 6" to 28" and water draft up to 13"
- Easy to install, even into existing concrete decks and anchors
- Powder coating availability to match any pool setting

Accessories

- Wheelchair attachment
- Spine board
- Spine board attachment and stability strap
- Protective cover
- Bronze 1.90" ID anchor with threaded cap and key
- Extra 24 volt battery and charger
- Lift maintenance kit
- Transport cart

800-791-8056 www.spectrumproducts.com