

BESTLIFT

Bestcare Foldable

NEW! PL400EF | FOLDABLE | PORTABLE

The Bestlift PL400EF represents a new standard in full body lifts for individuals who travel or are seeking storage efficincy in their home. PL400EF assembles (or folds) with 3 easy steps allowing for full operation or preparing for storage/transport.



- Light weight and easy to move
- 400lb capacity
- Optional scale available
- Optional spare battery and wall charger available
- Easy to service and disassemble



| **PL400EF**




DIMENSIONS

Bestcare Foldable

Key dimensions of PL400EF	
maximum capacity	400Lb / 182Kg
external base width (closed)	25.25 in
internal leg width (opened)	32 in
upright mast height	47.5 in
overall base length	40.75 in
base height	4.5 in
minimum lifting height with spreader bar	27.5
maximum lifting height with spreader bar	67.75
net weight	90Lb / 41Kg
shipping weight	104Lb / 47Kg

* SPECIFICATIONS PROVIDED ARE A GENERAL GUIDELINE AND ARE SUBJECT TO CHANGE WITHOUT NOTICE

TOP VIEW OF PL400EF



SIDE VIEW OF PL400EF







2920 PACIFIC DR., NORCROSS, GA 30071 WWW.BESTCARELLC.COM 678-679-6690