

## **FEATURES**

- Full platform provides excellent stability
- Durable welded construction
- Anti-slip, high traction surface
- Accommodates wheelchairs and scooters with various wheel configurations
- Safety DVD and steel security pins included



| Item Number | Length | Width | Weight  | Capacity* |
|-------------|--------|-------|---------|-----------|
| SL330       | 3'     | 30"   | 15 lbs. | 800 lbs.  |
| SL336       | 3'     | 36"   | 18 lbs. | 800 lbs.  |
| SL430       | 4'     | 30"   | 20 lbs. | 800 lbs.  |
| SL436       | 4′     | 36"   | 23 lbs. | 800 lbs.  |
| SL530       | 5′     | 30"   | 25 lbs. | 800 lbs.  |
| SL536       | 5′     | 36"   | 30 lbs. | 800 lbs.  |

Not sure which ramp length you need? Refer to pages 13-14 to help you.

\* Maximum ramp capacity: 400 lbs. (one axle) – 800 lbs. (two axles).