


DEALER			COLOURS INFORMATION For Office Use Only		
Dealer:			Order Date:		S/N:
Account:		P.O.#:	RTS		
Ordered By:			RTS:		Phone:
Address:			Email:		
City:		State:	SHIP TO		
Phone:		Fax:	Name:		Ship Via:
Email:			Address:		
Customer Name:			City:		State: Zip:
Customer Weight:			Phone:		Fax:
CHAIR INFORMATION			PRICE		
RAZOR BLADE K0005 Fixed FSH-RSH Max. Capacity 250Lbs			Base Wheel Chair Price K0005 \$ 2,400		
			Total Section A \$		
			Total Section B \$		
NOTE <u>Main Frame : ZK55 & ZK60 High Grade A7000 Series Aluminum</u> <u>Main Frame : Lifetime Warranty</u> *See warranty policy			Total Section C \$		
			Total Section D \$		
			Total Section E \$		
			Total Section F (Custom Parts Color) \$		
			ORDER GRAND TOTAL \$		

A FRAME COLORS & PARTS COLORS

FRAME COLORS						
Powder-Co	☐ Gloss Black	☐ Twilight Flake	☐ Yellow	☐ Magenta	☐ Fire Engine Red	☐ White
	☐ Candy Blue	☐ Candy Purple	☐ Candy Green	☐ Candy Red	☐ Candy Raspberry	☐ Cobalt Blue
	☐ Matte Gray	☐ Matte Purple	☐ Matte Blue	☐ Matte Green	☐ Matte Black	☐ Matte Red
	☐ Neon Red	☐ Neon Green	☐ Orange Sparkle	☐ Mellon Yellow	☐ Gun Smoke Silver	☐ Pigeon Blue
	☐ Granny Sparkle Green	☐ Sparkle Sky Blue	☐ Vermilion	☐ Reed Green	☐ Blue Lilac	STD
Special	☐ Polished					\$ 800
	☐ 24 karat Gold Plated					\$ CFP
PARTS COLORS (Custom Parts Color Section, Please see Section F)						
STD Color	Black Color (*Suspension Springs-Silver Color)					\$ STD
Custom Color Main Bolt on Parts One Color Package						\$ 500
Including:	Backrest, Foot Plate, Footrest Rings, Backrest Brackets, Saddle Stops,					\$ 750
Forks & Knurl Nuts, Caps, Front Suspension Springs, Axle Tube, Axle Clamps, Axle Plates,	☐ Powder-Coated ()					\$ CFP
Axle Bases, Camber Discs ☆marked parts	☐ Polished					
	☐ 24 karat Gold Plated					
The diagram shows a side view of a bicycle frame assembly. Labeled components include: ☆Saddle Stops (pointing to the rear rack), ☆Backrest (pointing to the rear seat post), ☆Backrest Brackets (pointing to the brackets holding the backrest), ☆Axle Plate & Axle Base (pointing to the rear axle assembly), ☆Camber Plate (pointing to the plate between the frame and the rear axle), ☆Axle Tube & Axle Clamps (pointing to the rear axle tube), ☆Caps (pointing to the caps on the front fork), ☆Front Springs *Std: Silver color (Powder-Coated Only) (pointing to the front suspension springs), ☆Forks & Knurl Nuts (pointing to the front fork assembly), ☆Footrest Rings (pointing to the rings on the footplate), ☆Foot Plate (pointing to the main platform), and ☆Axle Bases (pointing to the base of the front wheel).						
Total Section A						\$

B FRAME MEASUREMENTS 1/2



1. SEAT WIDTH	☐ 14" ☐ 15" ☐ 16" ☐ 17" ☐ 18" 	\$	STD								
2. SEAT DEPTH * Only Choose One Option	☐ Net Seat Depth ☐ 13" ☐ 14" ☐ 15" ☐ 16" ☐ 17" ☐ 18" ☐ 19" 	\$	STD								
	☐ 20" E2203	\$	300								
	☐ Depth Adj. Frame ☐ Short Frame (16" Depth) Set at ☐ 13" ☐ 14" ☐ 15" ☐ 16"	\$	STD								
	☐ Long Frame (19" Depth) Set at ☐ 16" ☐ 17" ☐ 18" ☐ 19"	\$	STD								
	☐ Positioning Frame Add (")	\$	150								
	☐ Ergo Seat Tubing Depth (") *Measurement must be 4" longer than your Center of Gravity	\$	200								
3. FOOTREST WIDTH	☐ STD (Seat Width Minus 5") ☐ Taper on Top ☐ Taper on Bottom	\$	STD								
	☐ Custom (") ☐ Taper on Top ☐ Taper on Bottom	\$	150								
	☐ No Taper (Seat Width Minus 2")	\$	STD								
FOOT PLATE STYLE	☐ Standard single loop with adjust. height & angle. STD Foot Plate K0040	\$	STD								
	☐ Single Flip-Back-single loop (no riser avail) with adjust. height & angle. STD Foot Plate K0040	\$	150								
	STD Foot Plate Size Chart <table border="1"> <thead> <tr> <th>Your Footrest Width</th> <th>Foot Plate Size</th> </tr> </thead> <tbody> <tr> <td>8.5"-10.5"</td> <td>W6.5"xD5"</td> </tr> <tr> <td>11.0"-13.0"</td> <td>W8.5"xD6.5"</td> </tr> <tr> <td>13.5"-17.0"</td> <td>W11"xD6.5"</td> </tr> </tbody> </table>	Your Footrest Width	Foot Plate Size	8.5"-10.5"	W6.5"xD5"	11.0"-13.0"	W8.5"xD6.5"	13.5"-17.0"	W11"xD6.5"		
Your Footrest Width	Foot Plate Size										
8.5"-10.5"	W6.5"xD5"										
11.0"-13.0"	W8.5"xD6.5"										
13.5"-17.0"	W11"xD6.5"										
CUSTOM FOOT PLATE	☐ Custom Foot Plate. with adjust. height & angle K0040 *Available Solid Plate Only	\$	75								
	☐ Drill 2Sets Extra Holes in Footplate *Measurement must be 2.5" shorter than your Footrest Width										
	☐ Custom Tubular Footrest	\$	250								
	☐ With ABS Foot Plate (\$65 up charge)										
4. FRAME ANGLE	☐ 75° ☐ 80° ☐ 85° ☐ 90° to the ground	\$	STD								
	☐ Custom (°) to the ground	\$	NCO								

B FRAME MEASUREMENTS 2/2

5. BACKREST HEIGHT * Rigidizer bar 3.25" deep std.		☐ 8" ☐ 9" ☐ 10" ☐ 11" ☐ 12" ☐ 13" ☐ 14" ☐ 15" ☐ 16" ☐ 17" ☐ 18"	\$ STD														
TYPE	☐ Locking mechan. with adjust. height & angle		\$ STD														
	☐ Double Locking mechan. with adjust. height & angle		\$ 25														
	☐ Folding with adjust. height & angle		\$ NCO														
	☐ Welded Backrest Angle Set at (") to the ground *Non Folding Backrest only ..		\$ 150														
CUSTOM	☐ Welded Side Guards *Non Folding Backrest only		\$ 450														
	☐ Custom Backrest (ie: double rigidizer bars, custom depth, custom height, etc)		\$ 150														
6. BACKREST ANGLE		Set at (°) to the ground *STD 90°	\$ STD														
7. FRONT SEAT HEIGHT		☐ 16" ☐ 16.5" ☐ 17" ☐ 17.5" ☐ 18" ☐ 18.5" ☐ 19" ☐ 19.5" ☐ 20" ☐ 20.5" ☐ 21"	\$ STD														
STD Front Seat Height Ranges <table border="1"> <thead> <tr> <th>Caster Wheel Size</th> <th>Front Seat Height</th> </tr> </thead> <tbody> <tr> <td>3" Caster</td> <td>16.0" - 19.5"</td> </tr> <tr> <td>4" Caster</td> <td>16.0" - 20.0"</td> </tr> <tr> <td>5" Caster</td> <td>16.0" - 21.0"</td> </tr> </tbody> </table>				Caster Wheel Size	Front Seat Height	3" Caster	16.0" - 19.5"	4" Caster	16.0" - 20.0"	5" Caster	16.0" - 21.0"						
Caster Wheel Size	Front Seat Height																
3" Caster	16.0" - 19.5"																
4" Caster	16.0" - 20.0"																
5" Caster	16.0" - 21.0"																
☐ CUSTOM FRONT SEAT HEIGHT (")		\$ 250															
8. REAR SEAT HEIGHT		☐ 14.5" ☐ 15.0" ☐ 15.5" ☐ 16.0" ☐ 16.5" ☐ 17.0" ☐ 17.5" ☐ 18.0" ☐ 18.5" ☐ 19.0" 	\$ STD														
STD Rear Seat Height Ranges *No Custom Modifications <table border="1"> <thead> <tr> <th>Rear Tire Size</th> <th>Rear Seat Height</th> </tr> </thead> <tbody> <tr> <td>24"x1"(540)</td> <td>14.5" - 18.0"</td> </tr> <tr> <td>24"x1 3/8"(540)</td> <td>15.0" - 18.0"</td> </tr> <tr> <td>25"x1"(559)</td> <td>15.5" - 18.0"</td> </tr> <tr> <td>24"x2"(540)</td> <td>15.5" - 18.5"</td> </tr> <tr> <td>26"x1"(590)</td> <td>15.5" - 19.0"</td> </tr> <tr> <td>25"x1.95"(559)</td> <td>16.0" - 19.0"</td> </tr> </tbody> </table>				Rear Tire Size	Rear Seat Height	24"x1"(540)	14.5" - 18.0"	24"x1 3/8"(540)	15.0" - 18.0"	25"x1"(559)	15.5" - 18.0"	24"x2"(540)	15.5" - 18.5"	26"x1"(590)	15.5" - 19.0"	25"x1.95"(559)	16.0" - 19.0"
Rear Tire Size	Rear Seat Height																
24"x1"(540)	14.5" - 18.0"																
24"x1 3/8"(540)	15.0" - 18.0"																
25"x1"(559)	15.5" - 18.0"																
24"x2"(540)	15.5" - 18.5"																
26"x1"(590)	15.5" - 19.0"																
25"x1.95"(559)	16.0" - 19.0"																
9. CENTER OF GRAVITY		☐ 1.0" ☐ 1.5" ☐ 2.0" ☐ 2.5" ☐ 3.0" ☐ 3.5" ☐ 4.0" ☐ 4.5" ☐ 5.5" ☐ 6.0" ☐ 6.5" ☐ 7.0"	\$ STD														
10. CAMBER ANGLE		☐ 2° ☐ 4° ☐ 6° 	\$ STD														
11. FRONT OF SEAT SLING TO FRONT OF FOOTREST		★ ☐ 15.0" ★ ☐ 15.5" ★ ☐ 16.0" ☐ 16.5" ☐ 17.0" ☐ 17.5" ☐ 18.0" ☐ 18.5" ☐ 19.0" ☐ 19.5" ★ Not Available w/ 19.0", 19.5", 20.0", 20.5", 21.0" Front Seat Height *Measurement must be 1.5" shorter than your Front Seat Height	\$ STD														
HIGH MOUNT ☐ 10.0" ☐ 10.5" ☐ 11.0" ☐ 11.5" ☐ 12.0" ☐ 12.5" ☐ 13.0" ☐ 13.5" ☐ 14.0" ☐ 14.5" ☐ 15.0" ☐ 15.5" ☐ 16.0" With Risers or Clamps		\$ 65															
12. INSIDE WHEEL CLEARANCE		☐ 3/4 ☐ 1" ☐ 1 1/4" ☐ 1 1/2" ☐ 1 3/4" ☐ 2" *3/4" Inside wheel clearance not available with 2 degree camber. (with the exception of Y-360 wheels) *3/4" Inside wheel clearance not available with Armrest *3/4" Inside wheel clearance not available with D's Locks	\$ STD														
NOTE Razor blade is a custom made wheelchair. If you would like to customize your Razor Blade beyond what is available on this order form please submit your request here. Colours will determine as to whether this request can be accomplished while maintaining the safety and integrity of the wheelchair. If it can be accomplished the quote will reflect any and all custom charges, thank you.																	
Sub Total Section B 1/2			\$														
Sub Total Section B 2/2			\$														
Total Section B			\$														




C WHEELS, TIRES & CASTERS, FORKS 1/2

[illegible]

C WHEELS, TIRES & CASTERS, FORKS 2/2

TIRES					
☐ High Pressure Tire	☐ 24"x1", ☐ 25"x1" or ☐ 26"x1"	Gray Tire with black wall Only		E2211/E2212	\$ STD
☐ High Pressure Tire	☐ 24"x1", ☐ 25"x1" or ☐ 26"x1"	☐ Blue ☐ Red		E2211/E2212	\$ 65
☐ Low Pressure Tire	☐ 24"x1-3/8" or ☐ 26"x1-3/8"	Gray Only		E2211/E2212	\$ STD
☐ E2211Schwalbe Marathon Plus High Pressure Tire	☐ 24"x1", ☐ 25"x1" or ☐ 26"x1"	Black Only		E2211/E2212	\$ 99
☐ Solid Tire	☐ 24"x1", ☐ 25"x1" or ☐ 26"x1"	☐ Dark Gray ☐ Black(\$30up charge)		E2218	\$ 135
☐ E2211 Moderate V-Track Knobby Tire	24"x1-3/8" (37-540)	Gray Only		E2211/E2212	\$ 117
☐ E2211 Aggressive Knobby Tire	24"x 2" (50-540)	Gray Only		E2211/E2212	\$ 117
☐ E2211 Aggressive Knobby Tire	26"x 1.95" (50-559)	Black Only		E2211/E2212	\$ 117
☐ E2213 Flat-Free Inserts For 24"x1-3/8" Tires Only				E2211/E2213	\$ 117
CASTERS					
SIZE	TIRE WIDTH	HUB COLOR	TIRE COLOR	DESCRIPTION	
☐ 3"	1"	Black Hub	Black Tire Only	Polyeurethane	\$ STD
☐ 3"	1"	Plastic Black Hub	Clear Tire Only	Lighted, Polyeurethane	\$ 65
☐ 3"	1"	Red Hub	Red Tire Only	Polyeurethane	\$ STD
☐ 3"	1"	Black Hub	Clear Tire Only	Polyeurethane	\$ STD
☐ 3"	1"	Aluminum Silver Hub	Black Tire Only	Made By Flog Legs	\$ 165
☐ 3"	1"	Aluminum Black Hub	Black Tire Only	Made By Flog Legs	\$ 165
☐ 3"	1.5"	Aluminum Silver Hub	Black Tire Only	Soft Roll Made By Flog Legs	\$ 165
☐ 3"	1.5"	Aluminum Black Hub	Black Tire Only	Soft Roll Made By Flog Legs	\$ 165
☐ 4"	1"	Plastic Black Hub	Black Tire Only		\$ STD
☐ 4"	1"	Plastic Black Hub	Clear Only	Lighted, Polyeurethane	\$ 65
☐ 4"	1"	Plastic Black Hub	Black Tire Only	Made By Flog Legs	\$ 125
☐ 4"	1"	Aluminum Silver Hub	Black Tire Only	Made By Flog Legs	\$ 165
☐ 4"	1"	Aluminum Black Hub	Black Tire Only	Made By Flog Legs	\$ 165
☐ 4"	1.5"	Composite Black Hub	Tire ☐ Black ☐ Blue ☐ Red	Soft Roll Made By Flog Legs	\$ 135
☐ 4"	1.5"	Aluminum Silver Hub	Tire ☐ Black ☐ Blue ☐ Red	Soft Roll Made By Flog Legs	\$ 165
☐ 4"	1.5"	Aluminum Black Hub	Black Tire Only	Soft Roll Made By Flog Legs	\$ 165
☐ 5"	1"	Composite Black Hub	Black Tire Only		\$ STD
☐ 5"	1"	Plastic Black Hub	Clear Tire Only	Lighted, Polyeurethane	\$ 65
☐ 5"	1"	Composite Black Hub	Black Tire Only	Made By Flog Legs	\$ 125
☐ 5"	1"	Aluminum Silver Hub	Black Tire Only	Made By Flog Legs	\$ 165
☐ 5"	1"	Aluminum Black Hub	Black Tire Only	Made By Flog Legs	\$ 165
☐ 5"	1.5"	Composite Black Hub	Tire ☐ Black ☐ Blue ☐ Red	Soft Roll Made By Flog Legs	\$ 135
☐ 5"	1.5"	Aluminum Silver Hub	Tire ☐ Black ☐ Blue ☐ Red	Soft Roll Made By Flog Legs	\$ 165
☐ 5"	1.5"	Aluminum Black Hub	Black Tire Only	Soft Roll Made By Flog Legs	\$ 165
FORKS					
☐ STD Aluminum Forks	Black Only				\$ STD
Extra Options for Aluminum Forks	☐ With Suspension Springs		E1015	\$ 55	
	☐ With Quick Release Axles		K0073	\$ 107	
☐ Frog Leg Suspension Forks Std 1" Wide (For 3", 4" x 1" width Caster and 5" STD Caster, 5" lighted Caster)	Color ☐ Anodized Black ☐ Anodized Silver		E1015	\$ 575	(Not available with Quick Release Axles)
☐ Frog Leg Suspension Forks 1.5" Wide (For 3", 4" x 1.5" width Caster and 5" x 1" & 5" x 1.5" Casters)	Color ☐ Anodized Black ☐ Anodized Silver		E1015	\$ 575	(Not available with Quick Release Axles)
Sub Total Section C 1/2					\$
Sub Total Section C 2/2					\$
Total Section C					\$


D COMPONENTS

WHEEL LOCKS		
☐ Push to Lock / Anodized Black Only		\$ STD
☐ NEW Compact Push to Lock / Anodized Black Only ☐ w/ 2" Extension Handle (\$65up)		\$ 110
☐ NEW Compact Scissor Lock / Anodized Black Only		\$ 110
☐ D's Locks (*Not available with X-Core, Topolino WX2.5, Pro Wheels, Magic wheels) (*Does not include installation fee)		\$ CFP
PUSH HANDLES *If you are ordering push handles, You must select the backrest locking mechanism.		
☐ Integrated w/Fixed Backrest Height / Anodized Black Only		\$ STD
☐ Bolt-on / Anodized Black Only Specify Height from floor -()		\$ 142
☐ Bolt-on w/Quick Release / Anodized Black Only Specify Height from floor -()		\$ 190
☐ NEW Height Adj. / Anodized Black Only Specify Maximum height from floor -()		\$ 295
NEW ARM RESTS Anodized Black Only,		
☐ Fixed Tubular Armrest	Adj. Angle and Height Height from seat sling Set at : ☐ 8" ☐ 9" ☐ 10" ☐ 11" ☐ 12" E0973	\$ 250
☐ Swing Away Tubular Armrest	Adj. Angle and Height Height from seat sling Set at : ☐ 8" ☐ 9" ☐ 10" ☐ 11" ☐ 12" E0973	\$ 250
☐ Flip-Up Armrest	☐ w/Fixed Tubular Armrest Height from seat sling Set at : ☐ 8" ☐ 9" ☐ 10" ☐ 11" ☐ 12" E0973	\$ 250
Adj. Angle ☐ w/ Swing Away Tubular Armrest	Height from seat sling Set at : ☐ 8" ☐ 9" ☐ 10" ☐ 11" ☐ 12" E0973	\$ 250
☐ Custom Length & Height (\$150 up charge) *Specify Height & Length → 		
ANTI-TIPPERS		
☐ Anti-Tippers (pair) w/ Quick Release / Anodized Black Only		E0971 \$ 112
SIDE GUARDS		
☐ Rigid Side Guards - Black ABS	☐ Non-Removable ☐ Removable ☐ Removable w/Hole (\$10 up charge) Size: ☐ STD ☐ Medium ☐ Cut to Wheel (\$98 up charge) Measurement	\$ 142
☐ Rigid Side Guards - Aluminum	☐ Non-Removable ☐ Removable Size: ☐ Medium Only ☐ Cut to Wheel (\$98 up charge) Measurement Color ☐ Black (STD) ☐ Pol. Aluminum (STD) ☐ Powder Coated () (\$50 up charge) ☐ Gold Plate (+CFP)	\$ 197
		
☐ PRO-TEX Fender Guards / (Black ABS -Only) ☐ Cut to Fit Specs ☐ Do Not Install		\$ 175
☐ Fabric Clothing Guards Black Nylon w/Velcro		\$ 78
SEAT CUSHION *Please see section E for other cushion choices.-		
☐ Foam Cushion w/Nylon Black Cover Thickness ☐ 1" ☐ 2" ☐ 3" ☐ 4"		\$ 78
☐ Custom Cushion Cover (Std. Foam Cushion Only) -Select Fabric () (\$90 up charge)		
MISCELLANEOUS		
☐ Frame Protectors / Black Leather/Zipper Type		\$ 78
☐ K0065 Spoke Protectors / Black Only ☐ 24" ☐ 25" ☐ 26"		K0065 \$ 107
☐ E0978 Seat Belt w/Auto Buckle (BP) 1.5" Black Belt Only		E0978 \$ 69
☐ Tie Down Brackets 4pcs.		\$ 210
SPECIALTY PARTS ☐ Wheelchair Tray Size: ☐ W14"x D12" ☐ W16"x D14" ☐ W18"x D14" Color: ☐ Raw Alum. ☐ Black Alum. ☐ Polish Alum. (\$250 up charge)		
\$ 250		
UPHOLSTERY CHOICE & COLOR *Positioning Back Available. Please see Section E		
☐ Standard Adjustable Ht Backrest Black Only		\$ STD
☐ Adjustable Tension Backrest Black Only (*Not available with 8" Backrest Height)		\$ 135
Custom Options	☐ Fabrics for Fun ☐ 001 ☐ 002 ☐ 003 ☐ 004 ☐ 006 ☐ 010 ☐ 011 . . . (\$95 up charge) For Std. or Adj. Style ☐ 014 ☐ 015 ☐ 016 ☐ 017 ☐ 018 ☐ 020 ☐ 021 ☐ 022 ☐ 023 ☐ 024 ☐ 025 ☐ 026 ☐ 027 ☐ 028 ☐ Backrest Leather Trim (outside panels) (\$77 up charge)	
☐ Calf Strap ☐ Black ☐ Matching Fabric (\$55 up charge) -Select Fabric()		K0038 \$ 24
☐ Matching Fabric on Padded Side Guards -Select Fabric ()		\$ 210
Total Section D \$		

CHAIR OPTIONS 1/2

SEAT CUSHION									
SUPRACOR® WHEELCHAIR CUSHIONS									
Stimulite® Contoured Cushion		☐ W14"xD16" (CD1416)						E2603	\$ 525
*Black Upholstery only		☐ W16"xD16" (CD1616)						E2603	\$ 525
		☐ W16"xD18" (CD1618)						E2603	\$ 525
		☐ W18"xD16" (CD1816)						E2603	\$ 525
		☐ W18"xD18" (CD1818)						E2603	\$ 525
		☐ W18"xD20" (CD1820)						E2603	\$ 540
		☐ W20"xD16" (CD2016)						E2603	\$ 540
		☐ W20"xD18" (CD2018)						E2603	\$ 540
		☐ W20"xD20" (CD2020)						E2603	\$ 555
		☐ Extra Sensitive (\$30 up charge)							
		☐ Sling Style Bottom (\$45 up charge)							
		☐ Custom Size W(") x D(") (\$- CFP)							
Stimulite® Slimline Cushion		☐ W14"xD16" (SM1416)						E2603	\$ 495
*Black Upholstery only		☐ W16"xD16" (SM1616)						E2603	\$ 495
		☐ W16"xD18" (SM1618)						E2603	\$ 495
		☐ W18"xD18" (SM1818)						E2603	\$ 495
		☐ W18"xD20" (SM1820)						E2603	\$ 495
		☐ W20"xD16" (SM2016)						E2603	\$ 510
		☐ W20"xD18" (SM2018)						E2603	\$ 510
		☐ W20"xD20" (SM2020)						E2603	\$ 510
		☐ Extra Sensitive (\$30 up charge)							
		☐ Sling Style Bottom (\$45 up charge)							
		☐ Custom Size W(") x D(") (\$- CFP)							
Stimulite® Corbee Pediatric		☐ W10"xD10" (PD1010)						E2603	\$ 450
*Black Upholstery only		☐ W12"xD12" (PD1212)						E2603	\$ 450
		☐ W14"xD14" (PD1414)						E2603	\$ 450
		☐ Extra Sensitive (\$30 up charge)							
		☐ Sling Style Bottom (\$45 up charge)							
		☐ Custom Size W(") x D(") (\$- CFP)							
ROHO WHEELCHAIR CUSHIONS									
☐ 4" HIGH PROFILE® QUADTRO SELECT® CUSHION		*Black Upholstery						E2622-2623	\$ 525
Choose Size	☐ W12"xD12"	☐ W12"xD13"	☐ W13"xD13"	☐ W13"xD15"	☐ W13"xD16"	☐ W15"xD13"	☐ W15"xD15"		
	☐ W15"xD16"	☐ W15"xD18"	☐ W16"xD15"	☐ W16"xD16"	☐ W16"xD18"	☐ W16"xD20"	☐ W18"xD15"		
	☐ W18"xD16"	☐ W18"xD18"	☐ W18"xD20"	☐ W18"xD22"	☐ W20"xD16"	☐ W20"xD18"	☐ W20"xD20"		
	☐ W20"xD22"	☐ W22"xD18"	☐ W22"xD20"	☐ W22"xD22"	☐ Custom Size W(") x D(") (+\$85)				
☐ 3" MID PROFILE® QUADTRO SELECT® CUSHION		*Black Upholstery only						E2622-2623	\$ 525
Choose Size	☐ W15"xD15"	☐ W15"xD16"	☐ W16"xD16"	☐ W16"xD18"	☐ W18"xD16"	☐ W18"xD18"			
	☐ Custom Size W(") x D(") (+\$85)								
☐ 2" LOW PROFILE® QUADTRO SELECT® CUSHION		*Black Upholstery only						E2622-2623	\$ 525
Choose Size	☐ W12"xD12"	☐ W12"xD13"	☐ W13"xD13"	☐ W13"xD15"	☐ W13"xD16"	☐ W15"xD13"	☐ W15"xD15"		
	☐ W15"xD16"	☐ W15"xD18"	☐ W16"xD15"	☐ W16"xD16"	☐ W16"xD18"	☐ W16"xD20"	☐ W18"xD15"		
	☐ W18"xD16"	☐ W18"xD18"	☐ W18"xD20"	☐ W20"xD16"	☐ W20"xD18"	☐ W20"xD20"			
	☐ Custom Size W(") x D(") (+\$85)								
☐ 2" NEXUS SPIRIT® CUSHION		*Black Upholstery only						E2622	\$ 525
Choose Size	☐ W14"xD14"	☐ W14"xD16"	☐ W15"xD15"	☐ W16"xD16"	☐ W16"xD18"	☐ W17"xD17"	☐ W18"xD16"		
	☐ W18"xD18"						* Custom Size Not Available		
STAR CUSHIONS									
Galaxy Cushion *Black Upholstery only						Choose Zoning Options		E2624	\$ 498
☐ W14"xD14" (GXY 1414-2)		☐ W14"xD16" (GXY 1416-2)		⇒		☐ Standard - Dual Valve			
☐ W16"xD14" (GXY 1614-2)		☐ W16"xD16" (GXY 1616-2)				☐ Dual Valve - Front to Back			
☐ W16"xD18" (GXY 1618-2)		☐ W18"xD16" (GXY 1816-2)				☐ Four Valve			
☐ W18"xD18" (GXY 1818-2)		☐ W18"xD20" (GXY 1820-2)				☐ Simple - Single Valve			
☐ W20"xD16" (GXY 2016-2)		☐ W20"xD18" (GXY 2018-2)				☐ Four Valve - Enhanced			
☐ W20"xD20" (GXY 2020-2)						☐ Four Valve - Front to Back			
THE CURVE CUSHION		☐ W10"xD10"	☐ W10"xD11"	☐ W10"xD11"	☐ W10"xD12"	☐ W11"xD12"	☐ W12"xD12"	\$ 133	
*Black Upholstery only		☐ W12"xD13"	☐ W12"xD14"	☐ W13"xD11"	☐ W13"xD12"	☐ W13"xD13"	☐ W13"xD14"		
		☐ W14"xD14"	☐ W14"xD15"	☐ W14"xD16"	☐ W15"xD14"	☐ W15"xD15"	☐ W15"xD16"		
		☐ W16"xD14"	☐ W16"xD15"	☐ W16"xD16"	☐ W16"xD17"	☐ W16"xD18"	☐ W17"xD16"		
		☐ W17"xD17"	☐ W17"xD18"	☐ W18"xD16"	☐ W18"xD17"	☐ W18"xD18"			
NOTE									

E CHAIR OPTIONS 2/2

POSITIONING BACK		
*Note : Actual height of the backrest will increase 2" or more. 		
ADI CARBON BACKREST (Support Contour-ACTIVE 2.5") *Not Available with Welded Side Guards		
Mounting Hardware: ☐ 2 Point Pro ☐ Quick Release		
Height	☐ LOW (10" Tall) *Backrest Height requires 12" or more *Available 13"-19" width Width (") E2613 \$ 745 ☐ Medium (13" Tall) *Backrest Height requires 15" or more *Available 13"-19" width Width (") E2613 \$ 815 ☐ Standard (16" Tall) *Backrest Height requires 18" or more *Available 13"-19" width Width (") E2613 \$ 880 ☐ Custom Size (" Tall) x (" width) E2613 \$ CFP	
ADI ALUMINUM BACKREST (Support Contour-ACTIVE 2.5") *Not Available with Welded Side Guards		
Mounting Hardware: ☐ 2 Point Pro ☐ Quick Release		
Height	☐ LOW (10" Tall) *Backrest Height requires 12" or more *Available 10"-20" width Width (") E2613 \$ 610 ☐ Medium (13" Tall) *Backrest Height requires 15" or more *Available 12"-20" width Width (") E2613 \$ 625 ☐ Standard (16" Tall) *Backrest Height requires 18" or more *Available 14"-20" width Width (") E2613 \$ 640 ☐ Tall (18" Tall) *Backrest Height requires 20" or more *Available 16"-20" width Width (") E2613 \$ 655 ☐ X-Tall (20" Tall) *Backrest Height requires 22" or more *Available 16"-20" width Width (") E2613 \$ 670 ☐ Custom Size (" Tall) x (" width) E2613 \$ CFP	
ADI ALUMINUM BACKREST (Support Contour-DEEP 4") *Not Available with Welded Side Guards *Not Available with Lateral Supports		
Mounting Hardware: ☐ 2 Point Pro ☐ Quick Release		
Height	☐ LOW (10" Tall) *Backrest Height requires 12" or more *Available 10"-20" width Width (") E2615 \$ 710 ☐ Medium (13" Tall) *Backrest Height requires 15" or more *Available 12"-20" width Width (") E2615 \$ 725 ☐ Standard (16" Tall) *Backrest Height requires 18" or more *Available 14"-20" width Width (") E2615 \$ 740 ☐ Tall (18" Tall) *Backrest Height requires 20" or more *Available 16"-20" width Width (") E2615 \$ 755 ☐ X-Tall (20" Tall) *Backrest Height requires 22" or more *Available 16"-20" width Width (") E2615 \$ 770 ☐ Custom Size (" Tall) x (" width) E2615 \$ CFP	
ADI LATERAL SUPPORTS (Lateral Supports are compatible with only active 2.5" ADI positioning back.)		
☐ Stationary Pad Pad Size ☐ Small (3" H x 3" D) ☐ Medium(3" H x 4" D) ☐ Large(4.5" H x 4.5" D) ☐ Large(4" H x 6" D)		
	☐ Pair \$ 320 ☐ Left side only or ☐ Right side only \$ 160	
☐ Swing-Out Pad Size ☐ Small (3" H x 3" D) ☐ Medium(3" H x 4" D) ☐ Large(4.5" H x 4.5" D) ☐ Large(4" H x 6" D)		
	☐ Pair \$ 520 ☐ Left side only or ☐ Right side only \$ 260	
OTHER		
☐ ADI Transfer Board ☐ 21" (ADV21) ☐ 24" (ADV24) ☐ 29"		\$ 62
☐ ADI Push and Transfer Gloves Leather (one size fits)		\$ 60
☐ ADI Carbon Solid Seat (Wide: 18" Max, Depth: 20") ※Colours Side Guards Not Available with ADI Carbon Solid Seat		\$ 625
With / ☐ ADI Folding Carbon Side Guard		\$ 400
With / ☐ ADI Fixed Carbon Side Guard		\$ 375
VARILITE® Seating Systems Available		
Please visit for all options. www.varilite.com		CFP
NOTE		
Sub Total Section E 1/2 \$		
Sub Total Section E 2/2 \$		
Total Section E \$		

F SPECIAL PARTS COLORS (Please allow an additional 1-2 weeks delivery)

PARTS COLORS

Powder-Coated	Gloss Black	Twilight Flake	Yellow	Magenta	Fire Engine Red	White
	Candy Blue	Candy Purple	Candy Green	Candy Red	Candy Raspberry	Cobalt Blue
	Matte Gray	Matte Purple	Matte Blue	Matte Green	Matte Black	Matte Red
	Neon Red	Neon Green	Orange Sparkle	Mellon Yellow	Gun Smoke Silver	Pigeon Blue
	Granny Sparkle Green	Sparkle Sky Blue	Vermilion	Reed Green	Blue Lilac	


Special Colors Polished Aluminum
24 karat Gold Plated

Please Choose the Color for each Part

*Parts not selected will be black color. (Suspension Spring: Silver Color)

*Minimum 3 parts selected per order.

*PC=Powder-Coated



Saddle Stops

- ☐ PC () \$40
- ☐ Polished..... \$115
- ☐ Gold Plated..... \$CFP

Backrest

- ☐ PC () \$125
- ☐ Polished..... \$225
- ☐ Gold Plated..... CFP

Backrest Brackets

- ☐ PC () \$90
- ☐ Polished..... \$115
- ☐ Gold Plated..... \$CFP

Foot Plate

- ☐ PC () \$45
- ☐ Polished..... \$90
- ☐ Gold Plated..... \$CFP

Footrest Rings

- ☐ PC () \$20
- ☐ Polished..... \$40
- ☐ Gold Plated..... \$CFP

Foot Plate Clamps

- ☐ PC () \$50
- ☐ Polished..... \$75
- ☐ Gold Plated..... \$CFP

Footrest Tubes

- ☐ PC () \$65
- ☐ Polished..... \$85
- ☐ Gold Plated..... \$CFP

Caps

- ☐ PC () \$20
- ☐ Polished..... N/A
- ☐ Gold Plated..... N/A

Front Springs

- ☐ PC () \$20
- ☐ Polished..... N/A
- ☐ Gold Plated..... N/A

Forks & Knurl Nuts

- ☐ PC () \$65
- ☐ Polished..... \$90
- ☐ Gold Plated..... \$CFP

Axle Plates

- ☐ PC () \$65
- ☐ Polished..... \$85
- ☐ Gold Plated..... \$CFP

Camber Plates

- ☐ PC () \$50
- ☐ Polished..... \$75
- ☐ Gold Plated..... \$CFP

Axle Bases

- ☐ PC () \$50
- ☐ Polished..... \$75
- ☐ Gold Plated..... \$CFP

Axle Clamp

- ☐ PC () \$50
- ☐ Polished..... \$75
- ☐ Gold Plated..... \$CFP


Axle Tube

- ☐ PC () \$50
- ☐ Polished..... \$75
- ☐ Gold Plated..... \$CFP


Total Section F \$

1. SEAT WIDTH


Measure from outside of seat tube to the outside of the opposite seat tube.

**2. SEAT DEPTH**


Measure from at the front bottom of backrest tube to the front edge of seat sling.

**3. FOOTREST WIDTH**


Measure from inside of front frame tube to inside of opposite front frame tube.

**4. FRAME ANGLE**


Measure angle frame front to floor.

**5. BACKREST HEIGHT**


Measure from top of back post to seat sling at front bottom of backrest tube.

**6. BACKREST ANGLE**


Measure angle from front of back post to floor.

**7. FRONT SEAT HEIGHT**

Measure from floor to top of seat sling at front edge.


**8. REAR SEAT HEIGHT**

Measure from floor to top of seat sling at front bottom of backrest tube.


**9. CENTER OF GRAVITY**

Measure from center of back tube at the folding pivot to the center of the rear axle.


(If welding backrest is selected, Measure from center of back tube at the seat tube to the center of the rear axle.)

**10. CAMBER ANGLE**

Measure angle from wheel to vertical at floor.

**11. FRONT OF SEAT SLING TO FOOTREST****ADULT CHAIR**

Measure from top of seat sling at front edge to top of foot plate at front edge.

**PEDIATRIC CHAIR**

Measure vertical from top of seat sling at front edge to top of foot plate.

**12. INSIDE WHEEL CLEARANCE**

Measure from outside of seat back tube to inside of rear tire.

(If side guards are selected, measure from outside of side guards to inside of rear tire.)

