

# PORTABLE VAN RAMP SYSTEM



## FEATURES & BENEFITS

- Folds For Easy Storage
- 10" to 12" Occupied Interior Space
- 1,000 lb Load Capacity
- All-Weather Traction
- Mounts to Floor of Many Standard Vans (Check for compatibility)



### Manual Folding Van Ramps

Installs into many standard factory vans. A more cost-effective alternative to a fully converted accessible van. Operated by raising and lowering the ramp from a folded position to the ground.

- Ramp is raised to an upright position inside the van door.
- Secured with included security strap.
- Ships as a complete unit ready for installation.
- Occupies approximately 10" (25.4 cm) depth (Removable if door space is needed).

Available ramp widths:  
26" (66 cm)  
30" (76 cm)  
36" (91.5 cm)

### REAR INSTALLATION



### Power Folding Van Ramps

An alternative to a fully-converted accessible van. Operated with the push of a button. Ramp deploys from a folded position inside the van down to the ground.

- Built-in safety feature: Ramp only moves when button is being pressed.
- Reliable direct-drive actuator motor.
- Ships as a complete unit with hardware.
- Pendant cord or wireless remote control options.
- Occupies approximately 12" (30 cm) depth in folding position.

Available ramp widths:  
26" (66cm)  
30" (76cm)

### SIDE INSTALLATION






CHECK COMPATIBILITY WITH VAN PRIOR TO ORDERING

1) Measure Door Width and Height 2) Measure Wheelbase of Van 3) Refer to Chart Below, or Call for Details

MINI VAN MANUAL RAMPS	
MF26	6' x 26" Manual Folding Ramp 40" door height
MF26.7	7' x 26" Manual Folding Ramp 46" door height
Requires 28" minimum side door opening	
MF30.6	6' x 30" Manual Folding Ramp 40" door height
MF30.7	7' x 30" Manual Folding Ramp 46" door height
Requires 32" minimum side door opening	
FULL SIZE VAN MANUAL RAMPS	
MF30.7	7' x 30" Manual Folding Ramp 46" door height
MF-30	8' x 30" Manual Folding Ramp 52" door height
MF30.9	9' x 30" Manual Folding Ramp 58" door height
Requires 32" minimum side door opening	
MF36.7	7' x 36" Manual Folding Ramp 46" door height
MF-36	8' x 36" Manual Folding Ramp 52" door height
MF36.9	9' x 36" Manual Folding Ramp 58" door height
Requires 38" minimum side door opening	
POWERED MINI VAN RAMPS	
MF1	6' x 26" Power Mini Fold Pendant Cord
MF1-R	6' x 26" Power Mini Fold Wireless Remote
Door opening: minimum 29" x 42"	
POWERED FULL SIZE VAN RAMPS	
AF1	7' x 30" Power Auto Fold Pendant Cord
AF1-R	7' x 30" Power Auto Fold Wireless Remote
Door opening: minimum 33" x 48"	



## VAN RAMPS ARE SHIPPED FREIGHT:

Commercial no liftgate: \$250  
Commercial w/ liftgate: \$335  
Residential w/ liftgate: \$375

## Roll-A-Ramp® Return Policy

\*20% Restocking Fee  
\*Outbound shipping is not refunded  
\*Return within 90 days from purchase  
\*Customer responsible for ramp return  
Refund issued upon return of items in acceptable condition.