

# rollz<sup>®</sup> Motion Performance

All terrain rollator and transport chair in one



## Smooth and safe ride on uneven terrain.

Thanks to the four pneumatic tyres, the Rollz Motion Performance provides a comfortable ride, even on bumpy roads. The extra thick padding of the armrests and the ergonomically shaped handles add to the comfort. Moreover, the footrests are adjustable in height.



### **Ergonomic grips**

Ergonomically shaped soft handles for optimal grip.

### **Soft padding**

Armrests with extra thick foam padding to increase comfort.

### **Firm seat**

Wedge-shaped cushion for a better sitting position.

### **Stable**

Walk tall and steadily, while safely supported by a sturdy design.

### **Visible**

Visible in the dark thanks to multiple reflective parts.

### **Eye-catching colour**

Modern matt coating in the new colour jungle green

### **Stylish 2 in 1**

Dutch design for going out in style, which transforms from rollator to transport chair within seconds.

### **Air filled tyres**

Sturdy pneumatic tyres enable more comfortable rides on bumpy roads.

## Accessories (Sold separately)



**Back Support**  
510-1020RM0012



**Cup Holder**  
510-1020RM0002



**Shopping Bag**  
510-1020RM0004



**Travel Cover**  
510-1020RM0010



**3-in-1 Cane, Chair & Bag Holder**  
510-2030RM0049

## Specifications



Scan to explore more!

Colour	Model Numbers
Jungle Green	2010RM0014
Description	
Seat Height	21.5" (55 cm)
Handle height range	34 - 39" (86 - 99cm)
Distance between rear wheels	22.5" (57cm)
Seat size	
Rollator	16.5" x 8" (42 x 20.3cm)
Transport chair	16.5 x 15" (42 x 38cm)
Overall width	25.5" (65 cm)
Overall length	24.8" (63 cm)
Depth folded	11.5" (29.2 cm)
Wheel diameter	
Front	8" (20.3 cm)
Rear	12" (30.5 cm)
Weight Capacity	275lb (125 kg)
Product weight	
Rollator	25.08 lb (11.4 kg)
Transport chair package	8.16 lb (3.7 kg)