

rollz[®] Motion²

Walker & Transport-chair in one



Now you can keep going...

You have already seen and done a lot. But your journey of discovery isn't over yet. Effortless staying power. Live life to the full. That's what Rollz is all about. Which is why we have designed the modern Rollz Motion2, a walker that converts into a transport-chair at a flick of the wrist. And of course, just as easily back again. Rollz takes you further.

Be it for an outdoor activity, a long holiday, a treasured exposition or a shopping afternoon, this 2-in-1 walker inspires you to keep on going without effort. In addition, the Rollz Motion2 folds in a compact package and fits into any car.

The Rollz Motion2 is an ergonomically designed walker with a strong but lightweight frame which enables you to walk steadily, upright and comfortably. And if you want to rest along the way, the Rollz Motion2 converts into a transport chair at a flick of the wrist. And of course, just as easily back again.

Rollz Motion² is available in four colours and two seat heights: Regular (21.5") and Small (19.5").



Black

Island Blue

Pebble White

Dark Purple

Compact

The Rollz Motion² is easily foldable into a compact package.

Ease

Cross doorsteps and curbs easily with the standard curb climbers.

Stable

Walking tall, steadily and supported between the ergonomic hand grips that are adjustable in height.

Comfortable

Comfortable ride and natural suspension thanks to the foam-filled rear tires.

Safe

Innovative drum brakes for safe and balanced braking.

Convenient

Resting on a comfortable cushioned seat.

Maneuverable

Light steering and easy motion with the maneuverable front wheels.

rollz[®] Motion²

Accessories (Sold Separately)




Cup Holder
510-1020RM0002



**3-in-1 Cane, Chair
& Bag Holder**
510-2030RM0049



Travel Cover
510-1020RM0010



Back Support
510-1020RM0012



Shopping Bag
510-1020RM0004

Rollz Motion² Specifications



product
design award
2012



reddot design award
honourable mention 2012

Colour	Model Numbers	
● Pebble White	2010RM0010	
● Island Blue	2010RM0011	2011RM0011
● Black	2010RM0012	2011RM0012
● Dark Purple	2010RM0013	
Description	Regular	Small
Seat height	21.5" (55 cm)	19.5" (50 cm)
Handle height range	34 - 39" (86 - 99 cm)	32 - 37" (81 - 94 cm)
Distance between rear wheels	22.5" (57 cm)	22.5" (57 cm)
Seat size		
Walker	16.5 x 8" (42 x 20.3 cm)	16.5 x 8" (42 x 20.3 cm)
Transport-chair	16.5 x 15" (42 x 38 cm)	16.5 x 15" (42 x 38 cm)
Overall width	25.5" (65 cm)	25.5" (65 cm)
Overall length	24.8" (63 cm)	24.8" (63 cm)
Depth folded	11.5 (29.2 cm)	11.5 (29.2 cm)
Wheel diameter		
Front	8" (20.3 cm)	8" (20.3 cm)
Rear	12" (30.5 cm)	12" (30.5 cm)
Weight capacity	275 lb (125 kg)	275 lb (125 kg)
Product weight		
Walker	23.7 lb (10.6 kg)	23.7 lb (10.6 kg)
Transport-chair Package	8.16 lb (3.7 kg)	8.16 lb (3.7 kg)