## Deck Profile Sheet POOL LIFT

\*Deck profile sheet must accompany your pool lift order

1 - PREFERRED LIFT PAL PAL2 Splash! Splash! 300 aXs2 multiLift multiLift2 ML300		
2 - GUTTER CONFIGURATION		
Image: Check One	Partially recessed gutter with or without bullnose coping	$A \rightarrow E$ $Bullnose coping$
A       Image: Constraint of the images above match your gutter configured	Flush gutter and deck with or without bullnose	Florida Rollout with curb
3 - FACILITY TYPE Health Club Hospitality Public Pool Residential Commercial		
CHECK ONE     Serial #      Pre-exisiting deck Replacing pre-existing lift- serial #		
<ul> <li>5 - DECK MATERIAL Concrete Pavers</li> <li>6 - DISTANCE FROM POOL DECK TO WATER LINE (A</li> <li>7 - HEIGHT OF COPING (B)</li> <li>8 - WIDTH OF COPING (C)</li> <li>9 - WIDTH OF GUTTER (D)</li> <li>10 - BULL NOSE RADIUS (E)</li> <li>11 - DECORATIVE STONE SETBACK</li> </ul>	"On a cle with mini mini	Check this box to confirm that you have verified the location of the lift will meet the ADA Clear Deck Space Requirement (1009.2.3): the side of the seat opposite the water, ear deck space shall be provided parallel the seat. The space shall be 36" wide mum and shall extend forward 48" mum from a line located 12" behind the edge of the seat."
12 - LOCATION OF POOL GROUND LEVEL	STORY LEVEL	edge of the seat.
NOTE: Pool lifts are application specic. Please provide accurate measurements for your pool in the space provided. S.R.Smith will conrm that the lift selected will meet the location and installation requirements based on the ADA Design Standards (2010) or suggest an alternative lift that will meet the requirements. S.R.Smith bears no responsibility due to misapplication of a lift without a completed Deck Prole Sheet on record.		
Name of distributor lift will be pruchased from	City	State
Your name	Email	Phone
Project name or property	City, State	PO Number
RETURN COMPLETED AND SIGNED FORM TO S.R.SMITH f: 503.266.4334 or e: customercare@srsmith.com		
S.R.SMITH USE ONLY         AXS2       MULTILIFT2       PAL HI/LO       SPLASH       SPLASH ER HI/LO       SPLASH 300         AXS2 ROUND POST       ML300       PAL SPA       SPLASH HI/LO       SPLASH ROUND POST       SPLASH 300 HI/LO         MULTILIFT       PAL       PAL2       SPLASH ER       SPLASH SPA       APPPOVAL #		
		APPROVAL #
SRSmith. A Fluidra Brand   srsmith.com   800.824.4		