Splash! Round Post



This fixed anchored lifting system can be installed in an existing round anchor.

- 3rd party tested & verified ADA compliant
- Integrated armrests
- · State of California compliant
- · 300 lb/159kg lifting capacity
- · Easily removable from stainless steel deck anchor
- Secure-It Kit[™] keeps the lift from rotating in the anchor socket
- LiftOperator® Intelligent Controller
- True operational redundancy lift can be operated from the control box and the handset
- 344° rotation ensures a safe transfer area
- Made from powder-coated stainless steel and aluminum
- Adapters available.
 For additional anchor sizes see table to the lower left



Lift Color



GRAY MIST

Due to printing technology actual color may differ.

Secure-It Kit™





Splash! Round Post Anchor Adapter



Adapters available for:

Model No.	Anchor Size
300-6210	23/8"
300-6220	21/2"
300-6230	29/16"





Note:

- Includes battery, charger, battery console cover, water resistant hand control, footrest, seat belt assembly, armrests and Secure-it-Kit
- * Package includes mast cover and caddy.

Parts & Accessories:

- 1001495 Battery
- 500-1000 Spineboard Attachment
- 920-5000T Cover
- 970-5000T Seat Saver Cover
- 900-1000 Seat Belt

- 400-0000 Caddy
- 900-6000 Stability Strap
- 300-1090A Splash! Round Post Secure-It Kit (standard)
- 170-3000A Armrest Replacement (pair)

