

THE SCOUT EXCEL



Item# F-SCTXL

The New Scout Excel pool lift is ADA compliant and UL Certified. It was designed to provide greater accessibility solutions with versatility, greater water drafts, and a setback up to 44 inches. Additional features of the new Scout Excel include an adjustable height seat pole, direct drive rotation, and it will accommodate most existing anchors by utilizing adapter sleeves. The optional pull-out leg rest may be used for greater clearance on raised pool and spa applications.

The Scout Excel is ideal for easy and affordable ADA compliance. Like all Aqua Creek lifts, the Scout Excel is constructed from all stainless steel with a durable powder coat finish in your choice of colors; standard white, gray, or tan. Seats and footrests are also available with color options in standard blue, white, gray, and tan.

Purchase the optional solar charging station to help extend the lift of your battery. The Scout Excel comes standard with the best warranty in the industry, a 5 year structural and 5 year prorated electronics warranty.

U.S. Patent Nos.
D507,769S



ACCENT COLORS

Available in a variety of accent and powder coat colors



*Blue seat standard (not a powder coat option)

OPTIONS/ACCESSORIES:

F-808SA	Anchor kit, standard, 1.9" x 6" deep
F-808SA-10	Anchor kit, paver applications, 1.9" ID x 10" deep
F-842SWDA	Anchor kit, wood deck applications
F-77SSA	Anchor kit, semi-recessed, shallow concrete deck
F-450SSC	Cover, standard lockable, blue
F-450SSC-G	Cover, standard lockable, gray
F-450SSC-T	Cover, standard lockable, tan
F-422PLH-2	Headrest
F-LRB	Pull-out leg rest assembly (blue, tan, gray, or white)
F-806SLS	Sling assembly
F-045SCH	Solar charging station
F-814SCTC	Transport cart


FEATURES:

- UL Certified
- 375 lb. (170 kg) weight capacity
- Setback range up to 44"
- 10" water draft
- Adjustable seat pole
- Direct drive rotation
- Flip-up armrests
- All stainless steel construction
- Durable powder coat finish
- Rechargeable battery & wall mount charger
- Five year structural warranty
- Five year pro-rated electronics warranty

PROUDLY SOLD BY:


THE SCOUT EXCEL

PRODUCT SPECIFICATIONS




APPLICABLE POOL PROFILES


ROLL OUT




PARTIALLY RECESSED




FULLY RECESSED



RAISED POOL / SPA




BULL NOSE




PARTIALLY RECESSED W/CURB



FLUSH W/ POOL DECK



FLORIDA ROLL OUT



*provide measurements indicated to confirm application or use online pool lift selector

Become a member of our new Dealer and Architect page. Access to CAD blocks, CSI specs and sales materials.