
AQUATIC LIFTS





The Traveler Long Reach BP 350 is a wireless battery powered lift that is self operable from the deck and water with a robust operating system. The Traveler can accommodate most decks and gutter setbacks for both in-ground and above ground pools.



Includes

- 5 year structural warranty
- 2 year electronics warranty
- Stabilizer bar for transfer
- Spectra Shield[®] coating for maximum corrosion protection
- 24 volt battery and charger

*Anchor sold separately

- 2 watertight remote controls
- Adjustable swing out footrest
- Flip-up armrests with non-slip grip
- Padded headrest
- Seat belt

Features

- 350 pound lifting capacity
- 90 lifts per battery charge
- Spacious seat with excellent lumbar support and ease of transfer
- Rotationally and vertically powered with 360° continuous power rotation
- Stainless steel construction (electropolished 304L grade)
- Accommodates setback of 6" to 40" and water draft up to 12"
- Easy to install, even into existing concrete decks and anchors
- Powder coating availability to match any pool setting

Accessories

- Wheelchair attachment
- Spine board
- Protective cover
- Bronze 1.90" ID anchor with threaded cap and key
- Extra 24 volt battery and charger
- Extra wireless remote
- Spine board attachment and stability strap
- Transport cart

800-791-8056 www.spectrumproducts.com