



THE PRIVESS™ BASIC

3-PANEL PRIVACY SCREEN

The Privess™ Basic provides portable privacy. Each screen includes casters for mobility and crutch tips for stationary use.

SPECIFICATIONS

Model Number	3130 / 3139
Overall Width (fully extended)	79" (200.7 cm)
Individual Panel Width	27" (68.6 cm)
Overall Height (with casters)	70" (177.8 cm)
Base Width	27" (68.9 cm)
Base Depth	20" (50.8 cm)
Individual Panel Length	49" (124.5 cm)



THE PRIVESS™ ELITE

3-PANEL DESIGNER PRIVACY SCREEN

The Privess™ Elite is the common choice for style-minded users who need a larger patient privacy solution matching their decor. Each screen includes casters for mobility and crutch tips for stationary use.

SPECIFICATIONS

Model Number	3170
Overall Width (fully extended)	84" (213.4 cm)
Individual Panel Width	28" (71.1 cm)
Overall Height (with casters)	70" (177.8 cm)
Base Width	31" (78.7 cm)
Base Depth	23" (58.4 cm)
Individual Panel Length	55" (139.7 cm)