Treatment Tables

champion.

FEATURES



WEIGHT CAPACITY: 400 lb. (181.4 kg)



PATIENT SURFACE (LxW): 72" (183 cm) X 28" (71 cm)



DURABILITY: Powder-coated. steel frames



STABILITY: Adjustable leveling feet



FUNCTIONALITY: All models include a paper holder and are available with or without a shelf (SH)



The 8500

Flat top

The 8500 flat-top treatment table features a durable dove-gray powder-coated steel frame, paper holder with cutter, adjustable leveling feet, and comfortable patient surface. Standard height is 30". Specialty heights (19", 25", or 36") are also available.



The 8550

Recovery Couch

The 8550 is engineered to provide years of service. It is commonly used in extended care, long-term care, and clinics. The comfortable contoured headrest combined with upholstered 2.5" medium-density foam provides exceptional comfort. Includes paper holder and cutter.



The 8570

Adjustable Back

The 8570 adjustable treatment table includes features such as a powdercoated steel frame, paper holder, adjustable leveling feet, and comfortable 2" medium-density foam. The headrest is gas spring operated and uses a 'one touch' adjustment control. Standard height is 30". Specialty heights (25" or 36") are also available.