

COLORS



Due to printing technology actual color may differ.



The multiLift2 with its field reversible seating design brings increased performance and flexibility to our best selling flange mount series.

- Field reversible seat with integrated armrests
- Rotomolded seat and footrest
- Powder-coated stainless steel and aluminum construction
- 350 lb/159kg lifting capacity
- Retrofit anchor jig is standard
- LiftOperator[®] Intelligent Controller
- California CEC Compliant
- Third-party tested and verified ADA compliant
- Made in the USA

Part #	Description	Shipping
580-0000	multiLift2, Gray Mist	
580-0000-ST	multiLift2, Slate	Box 1 100 lbs 60.5"L x 24.25"W x 18.5"H
580-0000-TP	multiLift2, Taupe	
580-0000N	multiLift2 w/o anchor, Gray Mist	Box 2
580-0000N-ST	multiLift2 w/o anchor, Slate	40 lbs 24" x 24" x 30"
580-0000N-TP	multiLift2 w/o anchor, Taupe	

Includes battery, charger, battery console cover, water-resistant hand control, footrest, seat belt assembly, armrests and retrofit anchor jig.

Parts and Accessories

Part #	Description
1001495	Battery
900-1000	Seat Belt
300-6700A	Anchors, set of 4
300-6900	Retrofit Anchor Jig
300-6800A	Anchor Bolts, set of 4
970-5000T	Seat Saver Cover
900-6000	Stability Strap
170-3000A	Armrest Replacement (pair) - Gray Mist
170-3000A-TP	Armrest Replacement (pair) - Taupe
70-3000A-ST	Armrest Replacement (pair) - Slate
1001601	Hand Control

New Construction Guidelines

Scenario	Pool Lift	Anchor Jig to Order
Order a Pool Lift and new construction jig at the same time	580-0000- X *	500-5000 (no anchors)
Order a new construction jig ahead of pool lift	580-0000N- X *	500-5000A (anchors included)

*Replace ${f X}$ listed above with the color

A completed Deck Profile form is required with your pool lift order. See srsmith.com/configurators



New Construction Jig with Anchors 580-5000A

800.824.4387 srsmith.com

No part of this may be copied or reproduced without the expressed written consent of S.R.Smith, LLC ©Copyright 2022 S.R.Smith, LLC. All rights reserved. 0122

