

aXs2 Pool Lift



The aXs2 is a low profile pool lift that is an ideal access solution for community and hospitality swimming pools.

- Compact design requires minimal deck space
- 300 lb/136kg lifting capacity
- 360° rotation for transfer safety and flexibility
- Integrated armrests
- Sturdy and generous rotomolded seat
- LiftOperator® Intelligent Controller with 24V battery
- Full range of operation from handset and control box
- Conveniently ships in two boxes



Lift Color



GRAY MIST

Due to printing technology actual color may differ

Caddy

Only available with pool lift kits – models:
310-0000C
310-0000CN



Model No.	Description	Shipping (2 – 3* boxes)			
		Weight	Length	Width	Height
310-0000	aXs2	102 lbs	55"	19"	14"
310-0000N	aXs2 without anchor	46kg	140cm	48cm	36cm
310-0000R	aXs2 Round Post without anchor	36 lbs	24"	24"	27"
310-0000C	aXs2 with Caddy	16kg	61cm	61cm	69cm
310-0000CN	aXs2 with Caddy, without anchor	50 lbs	30"	20"	57"
		23kg	76cm	51cm	145cm

Notes:

- Model 310-0000 includes battery, charger, battery console cover, water resistant hand control, LiftLock stainless steel anchor socket with cover, footrest, armrests and seat belt assembly
- Models 310-0000N and 310-0000CN come without a LiftLock Anchor, and the weight of Box 2 is 34 lbs/15kg
- Model 310-0000R is designed to fit in an existing round anchor with dimensions 2.38" ID x 6" deep
- aXs2 Caddy only available with pool lift purchase and ships in Box 3*

Parts & Accessories:

- 1001495 LiftOperator Battery
- 300-6200L LiftLock Anchor & Cap
- 970-5000T Seat Saver Cover
- 170-3000A Armrest Replacement (pair)
- AX9006 Cover
- 900-6000 Stability Strap
- 900-1000 Seat Belt
- 300-6200-2L LiftLock Anchor Replacement Cap

LiftLock

Locking anchor secures lift to deck

